



**QUALITY CONTROL RANGES
PYY (TOTAL) AND PYY (3-36) RIA KIT
CATALOG #: PYYT-66HK AND PYY-67HK**

**STANDARD, LYOPHILIZED
CATALOG #: 8066-K**

STANDARD	PYY VALUE
Lot Number. 14T66A	1,303 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED
CATALOG #: 6066-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

Kit	LOT NUMBER	EXPECTED RANGE	UNITS
PYY Total	HPYY-114A	47 – 98	pg/mL
	HPYY-214A	164 – 340	pg/mL
PYY 3-36	HPYY-114A	43 – 89	pg/mL
	HPYY-214A	165 – 343	pg/mL

6066-K: HPYY-114A.214A (02 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.