

QC RANGE SHEET
MILLIPLEX® MAP Human CVD Magnetic Bead Panel 6
Kit Catalog# HCVD6MAG-67K
QC Catalog# HCVMG6-6067-6
Lot# HCVM6-104 and HCVM6-204



<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CD14	Control 1	0.15 – 0.31	ng/mL
	Control 2	1.3 – 2.6	ng/mL
DPP4	Control 1	0.74 – 1.5	ng/mL
	Control 2	6.0 – 13	ng/mL
Endostatin	Control 1	0.04 – 0.08	ng/mL
	Control 2	0.31 – 0.64	ng/mL
Progranulin	Control 1	0.37 – 0.77	ng/mL
	Control 2	3.1 – 6.3	ng/mL
LRG1	Control 1	0.35 – 0.73	ng/mL
	Control 2	3.1 – 6.4	ng/mL
MCAM	Control 1	0.39 – 0.80	ng/mL
	Control 2	3.1 – 6.4	ng/mL
ucMGP	Control 1	0.79 – 1.6	ng/mL
	Control 2	6.2 – 13	ng/mL

HCVMG6-6067-6: HCVM6-104.204 (04 October 2022)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA