



## QUALITY CONTROL RANGES

### SARS-COV2 ANTI-RBD IgM ELISA CATALOG # EZRBDM-120K

#### STANDARD, LYOPHILIZED CATALOG # E8120-K

STANDARD	SARS-COV2 ANTI-RBD IgM VALUE
LOT NUMBER: 8120-K	50 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6120-K

Reconstitute each vial with **1.0mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	6120-101	3.1 – 6.5	ng/mL
QC 2	6120-201	9.7 - 20	ng/mL

E6120-K: 6120-101.201 (31 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)