

Calbiochem®

Antibody SourceBook™

Sixth Edition



Antibody SourceBook™

Sixth Edition


Primary Antibodies	5
Secondary Antibodies & Conjugates	57
Anti-Human Antibodies & Conjugates	56
Anti-Mouse Antibodies & Conjugates	56
Anti-Rabbit Antibodies & Conjugates	58
Anti-Rat Antibodies & Conjugates	59
Other Secondary Antibodies & Conjugates	59
Antibody Support Products	61
Adjuvants & Immunization Tools	61
Avidin, Streptavidin, Biotin & Conjugates	61
Extracts, Lysates, Cell Pellets & Blots	62
Markers	62
Immunodetection Reagents	62
Protein A & Protein G Reagents	63
Substrates for AP & HRP	63

Find the antibody you need quickly

Type the name of your target of interest into the search box -- and click.



findyourantibody.eu



Antibody guarantee

If our antibody fails to perform as described on the product data sheet we will offer you a replacement or a refund. Visit www.antibodyguarantee.eu

Antibody Sourcing Service

If we don't supply the antibody you need we will help you find an alternative supplier.

Visit www.FindYourAntibody.eu or contact our technical service (e-mail: techservice@merckbio.eu)

Primary Antibodies

Product	Applications	Species Reactivity	Cat No.	Size	Price
Primary Antibodies					
Anti-14-3-3 (221-242) Rabbit pAb	IB; NOT FS, NOT IF, NOT PS	human	PC70	100 µg	€319
Anti-14-3-3σ (69-81) Rabbit pAb	PS; NOT FS, NOT IB	human	PC482	100 µl	€294
Anti-19S Regulator ATPase subunit Rpt4 Mouse mAb (p42-23)	IB; NOT IP	human, mouse	ST1059	100 µl	€307
Anti-2,2,7-Trimethylguanosine Mouse mAb (K121)	IC, IP, PS	human, mouse, rat	NA02	100 µg	€338
Anti-2,2,7-Trimethylguanosine Mouse mAb (K121) Agarose Conjugate	AFF PUR, IB	human, mouse, rat	NA02A	0.5 ml	€458
Anti-5-Methylcytosine Mouse mAb (162 33 D3)	ELISA, FC, FS, IB, IF, PS, RIA	a broad range of species	NA81	50 µg	€313
Anti-A1 Adenosine Receptor (309-326) Rabbit pAb	IB, IH, IP	bovine (weakly), human, rat	119117	100 µl	€362
Anti-A1 Adenosine Receptor (309-326) Rabbit pAb	FS, IB, IP	bovine (weakly), human, rat	PC158	100 µl	€344
Anti-A20 Mouse mAb (59A426)	IB, IP	human	AM63	100 µg	€346
Anti-A2a Adenosine Receptor (Ab-1) (373-391) Rabbit pAb	FS, IB, IP	canine, human; not rat	PC159	100 µl	€344
Anti-c-Abl (Ab-1) Rabbit pAb	IB, IF, IP; NOT NS	human, mouse; not rat	PC01	100 µg	€319
Anti-c-Abl (Ab-2) Mouse mAb (19-110)	IB, IP, NT	human, mouse	OP19	100 µg	€338
Anti-c-Abl (Ab-3) Mouse mAb (24-21)	IB, IF, IP; NOT NS	human, mouse	OP20	100 µg	€338
Anti-c-Abl (Ab-3) Mouse mAb (24-21)	IB, IF, IP; NOT NS	human, mouse	OP20T	10 µg	€77
NEW Anti-c-Myb/v-Myb Mouse mAb (4E3)	IB, PS	human	CA1028	200 µl	€213
Anti-Acetylated Lysine Rabbit pAb	ELISA, IB, IP, PS	a broad range of species	ST1027	50 µl	€288
Anti-Acetylcholinesterase (Ab-2) Mouse mAb (ZR3)	FS, IP; NOT IB	feline, guinea pig, rabbit, rat; not chicken, not human, not mouse	NB04	100 µl	€319
Anti-ACINUS (Ab-2) (1329-1341) Rabbit pAb	IB	human	PC552	100 µg	€338
Anti-Acrp30 (18-32; 187-200) Rabbit pAb	IB, IP	mouse	107920	200 µg	€384
Anti-Actin (Ab-1) Mouse mAb (JLA20)	IB, IF	a broad range of species	CP01	1 each	€319
Anti-Actin (Ab-3) Mouse mAb (HHF-35)	FS, IB, PS	canine, human, monkey, mouse, rat, shrew	CP64	100 µl	€253
Anti-α-actin Mouse mAb (1A4)	FS, PS	baboon, human, rabbit	113200	500 µl	€378
Anti-Actin, Smooth Muscle Specific (Ab-2) Mouse mAb (1A4)	FS, IB, IF, PS	baboon, bovine, chicken, human, mouse, rabbit, rat	CP47	100 µg	€338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



← Apaf-1

NEW

Anti-Apaf-1 Mouse mAb (94408), Cat No. AP1038

Detection of Apaf-1 by immunoblotting. Sample: Whole cell lysate [2 x 10⁵ cells (~25 µg)] from human Jurkat cells. Primary antibody for immunoblotting: Anti-Apaf-1 (1 µg/ml). (Cat. No. AP1038) Detection: Chemiluminescence

NEW
NEW

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-ADAM10 (735-749) Rabbit pAb	IB	hamster, human, mouse	422751	100 µl	€309
Anti-ADAM10 (Ab-1) (732-748) Rabbit pAb	IB	human	PC528	100 µg	€338
Anti-ADAM9 Goat pAb	ELISA, IB, IP	human	PC669	100 µg	€508
Anti-ADAMTS4 Rabbit pAb	IB	human, mouse, rat	CA1027	100 µg	€305
Anti-Apaf-1 Mouse mAb (94408)	IB	human	AP1038	100 µg	€315
Anti-α-Adaptin (AP2) Mouse mAb (AP.6)	IB, IF, IP	bovine, human, monkey, mouse, rat	CP46	100 µl	€338
Anti-Adenine Nucleotide Translocase-1 (Ab-1) (2-11) Rabbit pAb	IB, IF, PS	human	PC522	100 µg	€294
Anti-Adenine Nucleotide Translocase-1 (Ab-2) (161-171) Rabbit pAb	IB, IC; NOT PS	human	PC523	100 µg	€294
Anti-Adenovirus 2 E1A Mouse mAb (M73)	IB, IC, IF, IP	adenovirus infected cells	DP11	100 µg	€307
Anti-Adenovirus 2 E1B 19 kDa (Ab-1) Mouse mAb (IF8)	IB, IF, IP	adenovirus infected cells	DP16	100 µg	€324
Anti-Adenovirus 2 E1B 19 kDa (Ab-2) Mouse mAb (3D11)	IB, IF, IP	adenovirus infected cells	DP17	100 µg	€324
Anti-Adenovirus 5 E1B 55 kDa Rat mAb (9C10)	IF, IP	adenovirus infected cells	DP08	100 µg	€324
Anti-Adrenergic α ₁ Receptor (339-349) Rabbit pAb	IB, IF	mouse, rat	PC160	100 µl	€319
Anti-Adrenergic α _{2A} Receptor (218-235) Rabbit pAb	IB, IH	mouse, rat	PC161	100 µl	€319
Anti-Adrenergic β ₁ Receptor (394-408) Rabbit pAb	IB, IF	mouse, rat	PC162	100 µl	€319
Anti-Agouti-Related Protein Rabbit pAb	ELISA, IB, IF, IH	mouse	122402	50 µg	€319
Anti-AIF (Ab-2) (589-600) Rabbit pAb	IB, PS; NOT IF	human, mouse, rat	PC536	100 µg	€294
Anti-Akt, PH Domain Mouse mAb (SKB1)	FC, IB, IP	human, rat	ST1088	1 Set	€384
Anti-Akt1 (88-100) Rabbit pAb	ELISA, IB, IP	human, mouse, rat, <i>Xenopus</i>	530311	100 µg	€343
Anti-Akt1 (Ab-1) (135-145) Rabbit pAb	IC; NOT IB, NOT PS	human, mouse	PC510	50 µl	€294
PhosphoDetect™ Anti-Akt1 (pSer ⁴⁷³) Mouse mAb (11E6)	ELISA, IB	human, mouse	124003	1 Set	€350
PhosphoDetect™ Anti-Akt1 (pThr ³⁰⁸) Rabbit pAb	IB	human, mouse	124001	10 T	€373
Anti-Akt2 Rabbit pAb	ELISA, IB	human, mouse, rat	124002	100 µl	€378
Anti-Albumin Goat pAb	IP	human	507313	1 ml	€62
Anti-Albumin Rabbit pAb	ELISA, IB	human	126584	1 ml	€176
Anti-Albumin Rabbit pAb	IEP, IP	bovine	126620	2 ml	€109
Anti-Albumin Rabbit pAb	IP, RID	chicken	126705	1 ml	€116
Anti-Alcohol Dehydrogenase Rabbit pAb	ELISA, IB	mouse, rat, yeast	126745	1 ml	€284
Anti-Aldehyde Dehydrogenase 1 (488-501) Goat pAb	IB	human, mouse	PC713	100 µg	€397
Anti-AML1 (Ab-1) Rabbit pAb	GS, IB, IP; NOT FS, NOT IF, NOT PS	human, mouse	PC284	100 µg	€319
Anti-AML1/RHD Domain (Ab-2) (50-177) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human, mouse	PC285	100 µg	€319
Anti-AML2 (38-54) Rabbit pAb	GS, IB; NOT FS, NOT IF, NOT IP, NOT PS	human, mouse, rat	PC286L	100 µg	€319
Anti-AML3 Rabbit pAb	GS, IB; NOT FS, NOT IF, NOT IP, NOT PS	human, mouse, rat	PC287L	100 µg	€319
Anti-AMPKα-2 Rabbit pAb	IB, IP	bovine, human, rat	ST1089	50 µg	€189
Anti-Amylase Rabbit pAb	ELISA, IP, RID	human	171534	1 ml	€192
Anti-β-Amyloid (1-40) Rabbit pAb	DB, ELISA, RIA	human	PC149	25 µg	€457
Anti-β-Amyloid (1-42) Rabbit pAb	DB, ELISA, RIA	human	PC150	25 µg	€457
Anti-β-Amyloid (1-43) Rabbit pAb	DB, ELISA, RIA	human	PC151	25 µg	€457
Anti-β-Amyloid (Asp-1) Rabbit pAb	IB, IP, PS	a broad range of species	PC729	25 µl	€278
Anti-β-Amyloid Mouse mAb (4G8)	ELISA, IB, IH, IP	human, mouse (weakly)	NE1002	100 µl	€435
Anti-β-Amyloid Mouse mAb (6E10)	ELISA, IB, IP, PS	human	NE1003	100 µl	€397
Anti-β-Amyloid Rabbit pAb	FS, IB	mouse, rat	NE1012	50 µl	€206
Anti-β-Amyloid, C-Terminal Mouse mAb (5C3)	ELISA, IB, IC	human	171604	100 µg	€409

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti- β -Amyloid, C-Terminal Mouse mAb (6G12)	ELISA, IB, IC	human	171607	100 μ g	€409
Anti- β -Amyloid, C-Terminal Mouse mAb (8G7)	ELISA, IB, IC	human	171606	100 μ g	€409
Anti- β -Amyloid, C-Terminal Mouse mAb (9F1)	ELISA, IB, IC	human	171605	100 μ g	€409
Anti- β -Amyloid, N-Terminal Mouse mAb (19H11)	ELISA, IB, IC	human	171603	100 μ g	€409
Anti-Amyloid Precursor Protein Mouse mAb (3E9)	IB, PS	human	171537	100 μ g	€378
Anti-Amyloid Precursor Protein, C-Terminal (681-695) Goat pAb	ELISA, IB	human	171599	100 μ l	€312
Anti-Amyloid Precursor Protein, C-Terminal (751-770) Rabbit pAb	IB, IF, IH, IP	human	171610	50 μ l	€304
Anti-Amyloid Precursor Protein, N-Terminal (44-63) Goat pAb	ELISA, FS, PS	human	171598	100 μ l	€281
Anti- β -Amyloid Precursor-Like Protein 1, C-Terminal (643-653) Rabbit pAb	IB, IF, IP	human, mouse, rat	171615	100 μ l	€345
Anti- β -Amyloid Precursor-Like Protein 2, C-Terminal (752-763) Rabbit pAb	IB, IF, IP	human, mouse, rat	171616	100 μ l	€345
Anti- β -Amyloid Precursor-Like Protein 2, Full-length Rabbit pAb	FS, IB, IF, IP, PS	human, mouse, rat	171617	100 μ l	€345
Anti- β -Amyloid ₄₀ (FCA3340) Rabbit pAb	ELISA, IB, IF, IP, PS	canine, human, Microcebus murinus, mouse	171608	50 μ l	€409
Anti- β -Amyloid ₄₂ (FCA3542) Rabbit pAb	ELISA, IB, IF, IP, PS	canine, human, Microcebus murinus, mouse	171609	50 μ l	€409
Anti-Androgen Receptor Rat mAb (AN1-15)	FS, IB, IC, IP, PS	human, mouse, primate, rat	CA1022	50 μ g	€378
Anti-Angiogenin Goat pAb	ELISA, IB, NT	human	PC317L	100 μ g	€338
Anti-Angiopoietin-2 Rabbit pAb	ELISA, IB, IF, IH	mouse	176002	50 μ g	€316
Anti-Angiostatin (Ab-1) Rabbit pAb	IB, IF	human	PC371	100 μ g	€338
Anti-Angiostatin (Ab-1) Rabbit pAb	IB, IF	human	PC371T	20 μ g	€86
Anti-Angiostatin (Ab-2) Mouse mAb (Angio53 1-188)	IB, IP; NOT PS	human; not mouse, not rat	GF47	100 μ g	€338
Anti-Ankyrin ₉ Mouse mAb (2.20)	FS, IB, IF	human, mouse, rat	NB19	100 μ g	€338
Anti-Ankyrin ₉ Mouse mAb (463)	FS, IB, IF	canine, human, mouse, rat	NB20	100 μ g	€338
Anti-ANT Mouse mAb (5F51BB5AG7)	IB, IC, IP	bovine, human, rat	AP1034	50 μ g	€294
Anti- α_1 -Antitrypsin Mouse mAb (TMF1#4B5)	EIA, RIA	human	178260	1 mg	€206
Anti-AP-2 (120-134) Rabbit pAb	ELISA, IB	human	PC692	100 μ l	€409
NEW Anti-Apaf-1 Rabbit pAb	IB	human	AP1052	50 μ l	€135
Anti-APC (Ab-1) Mouse mAb (FE9)	IB	human, mouse, rat	OP44	100 μ g	€363
Anti-APC (Ab-2) Mouse mAb (IE1)	IB	human	OP47L	100 μ g	€346

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherteriony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-ANT Mouse mAb (5F51BB5AG7) Cat. No. AP1034

Detection of human ANT by immunocytochemistry. Sample: Cultured fibroblasts. Primary antibody: Anti-ANT Mouse mAb (5F51BB5AG7) (Cat. No. AP1034) (5 μ g/ml). Detection: Fluorescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-APC (Ab-5) Mouse mAb (CF11)	IP	human, mouse	OP62	100 µg	€346
Anti-APC (Ab-7) Mouse mAb (CC-1)	FS, IF, PS	human, mouse, rat	OP80	100 µg	€346
Anti-mAPH-1a ¹ (245-265) Rabbit pAb	IB, IP	human, mouse, rat	PC728	100 µl	€384
Anti-Apolipoprotein A-I Goat pAb	IEP, IP	human	178463	1 ml	€87
Anti-Apolipoprotein A-I Mouse mAb (412)	EIA, ELISA, RIA	human	178474	1 mg	€218
Anti-Apolipoprotein A-I Mouse mAb (513)	EIA, ELISA, RIA	human	178472	1 mg	€218
Anti-Apolipoprotein A-I Mouse mAb (6001)	EIA, ELISA, RIA	human	178470	1 mg	€218
Anti-Apolipoprotein A-I Rabbit pAb	ELISA, IB	human	178422	1 ml	€87
Anti-Apolipoprotein A-II Goat pAb	IEP, IP	human	178464	1 ml	€102
Anti-Apolipoprotein B Goat pAb	IEP, IP	human	178467	1 ml	€87
Anti-Apolipoprotein B Mouse mAb	EIA, ELISA, RIA	human	178444	1 mg	€285
Anti-Apolipoprotein B Mouse mAb (F ₂ C ₉)	EIA, ELISA, RIA	human	178446	1 mg	€285
Anti-Apolipoprotein C-II Rabbit pAb	ELISA, IB	human	178428	1 ml	€114
Anti-Apolipoprotein E Goat pAb	ELISA, IB, IP	human	178479	500 µl	€285
Anti-Apolipoprotein E Mouse mAb (E6D7)	ELISA, FS, IB, IP, PS	human; not mouse, not rabbit	NE1004	100 µl	€397
Anti-Apollon (4775-4829) Rabbit pAb	IB, IP	human	AP1031	50 µg	€183
Apoptosis Antibody Sampler Kit		human	ASK06	1 each	€270
Anti-Aquaporin 0 Rabbit pAb	ELISA, IB, IF, IH	bovine, human, mouse, rat	178610	50 µg	€319
Anti-Aquaporin 1 (243-261) Rabbit pAb	FS, IB	human, rat	178611	1 Set	€206
Anti-Aquaporin 2 (254-271) Rabbit pAb	IB, IH, IP	human, mouse, rat	178612	1 Set	€206
Anti-Aquaporin 3 (275-292) Rabbit pAb	IB, IH	rat	178613	1 Set	€206
Anti-Aquaporin 5 Rabbit pAb	ELISA, IB, IF, IH	human, mouse, rat	178615	50 µg	€319
Anti-ARC (2-18) Rabbit pAb	IB	human, mouse, rat	PC486	100 µg	€338
Anti-ARHI Mouse mAb (15E11)	ELISA, IB, PS	human	CA1011	50 µg	€319
Anti-Arrestin (347-363) Rabbit pAb	IB	bovine	178605	200 µg	€384
Anti-Arrestin (384-397) Rabbit pAb	IB	bovine, human, mouse, rat	178600	200 µg	€384
NEW Anti-ASC Rabbit pAb	IB, IP	human	ST1121	50 µl	€135
Anti-ASK1 (1360-1372) Rabbit pAb	IB	human	PC471	100 µg	€294
PhosphoDetect™ Anti-ATF2 (pSer ^{490/498}) Rabbit pAb	DB, IB, IH	human	PS1001	100 µl	€385
PhosphoDetect™ Anti-ATF-2 (pThr ⁷¹) Rabbit pAb	FFS, IB, IC, IP, PS	human, mouse, rat	ST1075	50 µl	€199
Anti-ATF-4 (339-351) Goat pAb	IB	human	PC715	100 µg	€397
Anti-ATM (Ab-1) (13-24) Rabbit pAb	IB	human, mouse	PC85	100 µg	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-ASC Rabbit pAb, Cat No. ST1121

Detection of human ASC by immunoblotting (Lane 1): Sample: Whole cell lysate (25 ug) from HT-29 cells. Primary antibody for immunoblotting: Anti-ASC (1:30000). (Cat. No. ST1121) Detection: Chemiluminescence

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-ATM (Ab-2) Mouse mAb (6B7)	IP; NOT IB	human	OP90	200 µg	€338
Anti-ATM (Ab-3) (819-844) Rabbit pAb	IB, IP	human	PC116	100 µg	€319
Anti-ATM (Ab-3) (819-844) Rabbit pAb	IB, IP	human	PC116T	10 µg	€77
PhosphoDetect™ Anti-ATM (pSer ¹⁹⁸¹) Mouse mAb (10H11.E12)	IB, IC	human, mouse	DR1002	50 µg	€170
Anti-ATR (Ab-2) (2122-2644) Rabbit pAb	IB	hamster, human, mouse	PC538	50 µl	€338
Anti-ATR (Ab-2) Rabbit pAb	IB	hamster, human, mouse	PC538T	5 µl	€77
Anti-Aurora-A Rabbit pAb	IP, PS	human	PC742	100 µl	€384
Anti-Avidin Goat pAb	ELISA, IH	chicken	189731	1 ml	€112
Anti-BACE (26-45) Rabbit pAb	IB, IP	human	195100	100 µl	€301
Anti-BACE (46-65) Rabbit pAb	IB, IP	human	195101	100 µl	€301
Anti-BACE (487-501) Rabbit pAb	IB, IP	hamster, human	195102	100 µl	€301
Anti-BACE (Ab-1) (44-59) Rabbit pAb	PS	human	PC478	50 µl	€338
Anti-BACE (Ab-2) (485-501) Rabbit pAb	IB	human, mouse	PC529	100 µg	€338
Anti-BACE1, C-Terminal (485-501) Rabbit pAb	IB, IF, IH, IP, RIA	human, mouse, rat	195111	100 µl	€413
Anti-BACE2 (Ab-1) (496-511) Rabbit pAb	IB	human	PC555	100 µg	€338
Anti-BACE2 (Ab-2) (44-59) Rabbit pAb	IB	human, mouse, rat	PC556	100 µg	€338
Anti-Bad (19-35) Rabbit pAb	IB, PS	bovine, canine, hamster, human, monkey, mouse, ovine, porcine, rabbit, rat	195872	100 µg	€378
PhosphoDetect™ Anti-Bad (pSer ¹¹²) (Ab-1) Rabbit pAb	IB, IP	human, mouse, rat	PC518	50 µl	€282
PhosphoDetect™ Anti-Bad (pSer ¹¹²) (Ab-3) Rabbit pAb	IB	human, mouse, rat	PC636	50 µg	€384
PhosphoDetect™ Anti-Bad (pSer ¹³⁶) (Ab-2) Rabbit pAb	IB	human, mouse, rat	PC519	50 µl	€282
PhosphoDetect™ Anti-Bad (pSer ¹³⁶) (Ab-4) Rabbit pAb	IB	human, mouse, rat	PC637	10 T	€373
Anti-BAFF (Ab-2) Rat mAb (Buffy-1)	FC, IB, IC	human	AM69	25 µg	€248
Anti-Bak (2-14) Rabbit pAb	ELISA, IB, IP	Chinese hamster, human, porcine	196150	100 µg	€378
Anti-Bak (Ab-1) Mouse mAb (TC-100)	IB	human; not mouse, not rat	AM03	20 µg 100 µg	€113 €338
Anti-Bak (Ab-1) Mouse mAb (TC-100)	IB	human; not mouse, not rat	AM03T	10 µg	€70
Anti-Bak (Ab-2) Mouse mAb (TC-102)	IB, IF	human, mouse	AM04	20 µg 100 µg	€86 €338
Anti-Bax (44-62) Rabbit pAb	IB	bovine, human, mouse, rabbit (weakly), rat	196821	100 µl	€378

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-BMI1 (1-326) Rabbit pAb, Cat No. AP1050

Detection of human BMI1 by immunoblotting. Sample: Whole cell lysate (25 µg) from U2OS cells. Primary antibody: Anti-BMI1 (1-326) Rabbit pAb (Cat. no. AP1050) (1:1000). Detection: Chemiluminescence

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Bax (Ab-1) (150-165) Rabbit pAb	FS, IB, PS	human, mouse, opossum, rat	PC66	100 µg	€319
Anti-Bax (Ab-1) Rabbit pAb	FS, IB, PS	human, mouse, opossum, rat	PC66T	10 µg	€77
Anti-Bax (Ab-2) Mouse mAb (ID3)	IB, IF, IP; NOT FS, NOT PS	human	AM13	100 µg	€338
Anti-Bax (Ab-3) Mouse mAb (2D2)	IB	human	AM32	100 µg	€338
Anti-Bax (Ab-4) (98-117) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC103	100 µg	€319
Anti-Bax (Ab-6) Mouse mAb (6A7)	IP	human, mouse, rat	AM44	100 µg	€338
Bax Peptide	CT, DB		PP51	100 µg	€84
Anti-Bcl-2 (Ab-1) Mouse mAb (100)	FS, IB, PS	human; not mouse	OP60	100 µg	€294
Anti-Bcl-2 (Ab-1) Mouse mAb (100)	FS, IB, PS	human; not mouse	OP60T	10 µg	€77
Anti-Bcl-2 (Ab-2) (20-34) Rabbit pAb	FS, IB, PS	human, mouse, opossum	PC68	100 µg	€319
Anti-Bcl-2 (Ab-2) Mouse mAb (8C8)	FS, IB, IP, PS	human, monkey; not mouse, not rat	AM59	100 µg	€338
Anti-Bcl-2 (Ab-2) Rabbit pAb	FS, IB, PS	human, mouse, opossum	PC68T	10 µg	€77
Anti-Bcl-2 (Ab-3) Mouse mAb (4D7)	FC, IB, IF, IP, PS	human; not mouse, not rat	OP91	100 µg	€294
Anti-Bcl-2 (Ab-4) Mouse mAb (10C4)	IB	mouse, rat	AM43	100 µg	€338
Anti-Bcl-2 (Ab-4) Mouse mAb (10C4)	IB	mouse, rat	AM43T	10 µg	€86
PhosphoDetect™ Anti-Bcl-2 (pSer ⁸⁷) Rabbit pAb	PS	human	PC502	25 µg	€282
Bcl-2 Family Antibody Sampler Kit		human	ASK12	1 each	€461
Bcl-2 Peptide	CT, DB		PP52	100 µg	€84
Anti-Bcl-2α (Ab-1) Mouse mAb (100/D5)	FS, IB, IP, PS	human	AM58	100 µg	€338
Anti-Bcl-w (16-29) Rabbit pAb	IB, IH, IP	bovine, human, monkey, mouse, rat	197209	100 µg	€378
Anti-Bcl-w Rat mAb (16H12)	IB, IF, IP	human, mouse	AM54	100 µg	€338
Anti-Bcl-x (Ab-1) (201-216) Rabbit pAb	FS, IB, PS	human, mouse, opossum, rat	PC67	100 µg	€319
Anti-Bcl-x (Ab-3) (150-169) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC104	100 µl	€319
Anti-Bcl-x _L (Ab-2) Mouse mAb (2H12)	IB	human	AM05	100 µg	€338
Anti-Bcl-x _S (Ab-1) (119-131) Rabbit pAb	IB, PS	human; not mouse, not rat	PC89	100 µg	€319
Anti-Bcr (Ab-2) Mouse mAb (7C6)	IB, IF, IP	human; not mouse, not rat	OP26	100 µg	€324
Anti-BCRP Mouse mAb (BXP-21)	FS, IB, IC, PS	human	OP191	200 µl	€248
Anti-BDNF Mouse mAb (35928.11)	ELISA, IB, IH, NT	human	GF35L	100 µg	€270
Anti-Bid (Ab-1) (76-85) Rabbit pAb	PS; NOT IB	human; not mouse	PC562	100 µg	€294
Anti-Cleaved Bid (Ab-1) Rabbit pAb	IB	mouse; not human	PC645	10 T	€457
Anti-Bim (22-40) Rabbit pAb	IB, IC	human, mouse, rat	202000	100 µg	€350
Anti-Bim Rat mAb (14A8)	IB, IF, IP	human, mouse	AM53	50 µg	€384
Anti-BLM Rabbit pAb	IB, IP	human	DR1034	50 µg	€189
Anti-Blood Group B Antigen Mouse mAb (CLCP-19B)	ELISA, IB, IH	human	202769	1 ml	€314
Anti-Bmf (2-14) Rabbit pAb	IB	human, mouse	PC685	100 µg	€319
NEW Anti-BMI1 (1-326) Rabbit pAb	IB	human	AP1050	50 µl	€135
Anti-BMP-1 (Ab-1) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC571	50 µg	€253
Anti-BMP-2/4, C-Terminal Rabbit pAb	ELISA, IB, IP	human	203640	100 µg	€362
Anti-BNF1 (150-250) Rabbit pAb	IB	human	CA1021	50 µg	€176
Anti-Bnip3L (77-92) Rabbit pAb	IB	human	PC525	100 µg	€338
Anti-Bone Sialoprotein Mouse mAb (ID1.2)	ELISA, IB, PS	human	203637	100 µg	€350
Anti-BRCA1 (Ab-1) Mouse mAb (MS110)	FS, IB, IF, IP, PS	human	OP92	100 µg	€338
Anti-BRCA1 (Ab-1) Mouse mAb (MS110)	FS, IB, IF, IP, PS	human	OP92T	10 µg	€77
Anti-BRCA1 (Ab-2) Mouse mAb (MS13)	IB, IF, IP, PS	human	OP93	100 µg	€338


10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-BRCA1 (Ab-2) Mouse mAb (MS13)	IB, IF, IP, PS	human	OP93T	10 µg	€77
Anti-BRCA1 (Ab-4) Mouse mAb (SD118)	FS, IB, IF, IP; NOT PS	human	OP107	100 µg	€338
PhosphoDetect™ Anti-BRCA1 (pSer ¹³⁸⁷) Rabbit pAb	ELISA, IB, IP	human	PC631	50 µl	€233
PhosphoDetect™ Anti-BRCA1 (pSer ¹⁴²³) Rabbit pAb	ELISA, IB, IP	human	PC632	50 µl	€221
PhosphoDetect™ Anti-BRCA1 (pSer ¹⁵²⁴) Rabbit pAb	IB	human	PC732	100 µl	€385
Anti-BRCA2 (Ab-1) Mouse mAb (2B)	IB, IP	human	OP95	100 µg	€338
Anti-BRCA2 (Ab-2) (3245-3418) Rabbit pAb	IB, PS	human	PC146	100 µg	€314
Anti-BRCA2 (Ab-2) Rabbit pAb	IB, PS	human	PC146T	10 µg	€77
Anti-BrdU (Ab-2) Mouse mAb (ZBU30)	FS, IF, PS	a broad range of species	NA20	100 µg	€275
Anti-BrdU (Ab-3) Mouse mAb (Mobu-1)	FS, IC, PS; NOT IB	a broad range of species	NA61	100 µg	€338
Anti-BRST-2 Mouse mAb (D6)	PS	human	CA1001	200 µl	€338
Anti-Bruton's Tyrosine Kinase Mouse mAb (10D11)	IB, IP; NOT IF	human	203665	50 µl	€231
Anti-CA 19-9 Mouse mAb (116-NS-19-9)	PS	human	CA1003	200 µl	€338
Anti-CA 125 Mouse mAb (OC 125)	PS	human	CA1004	200 µl	€338
Anti-Ca _v 1.2 (818-835) Rabbit pAb	IB, IC, IP	mouse, rabbit, rat	681507	1 Set	€212
Anti-Ca _v 1.3 (809-825) Rabbit pAb	IB	mouse, rat	681509	1 Set	€212
Anti-Ca _v 2.2 (851-867) Rabbit pAb	IB, IC, IH, IP	mouse, rat	681505	1 Set	€212
Anti-Ca ²⁺ /Calmodulin Kinase II α-Subunit (Ab-1) Mouse mAb (6G9)	IB, IF	rat	NB12	100 µl	€319
PhosphoDetect™ Anti-Ca ²⁺ /Calmodulin Kinase II α-Subunit (pThr ²⁸⁶) (Ab-2) Mouse mAb (22B1)	ELISA, IB, IF	rat	NB13	100 µl	€319
PhosphoDetect™ Anti-Ca ²⁺ /Calmodulin Kinase II (pThr ²⁸⁶) Rabbit pAb	DB, IB	human	208741	100 µl	€412
Anti-Ca ²⁺ /Calmodulin Kinase II, N-Terminal (7-20) Rabbit pAb	IB	mouse, rat	208703	100 µg	€397
Anti-Ca ²⁺ /Calmodulin Kinase IV, C-Terminal (465-474) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	208709	100 µg	€378
Anti-Ca ²⁺ /Calmodulin Kinase Kinase, C-Terminal (490-503) Rabbit pAb	ELISA, IB	human, mouse, rat	208717	100 µg	€384
NEW Anti-Caldesmon Mouse mAb (TD107)	FS, IB, IC	human	ST1109	50 µg	€171
Anti-CAD (205-221) Rabbit pAb	IB	mouse	205543	100 µg	€319
Anti-E-Cadherin (ECCD-1) Rat mAb (ECCD-1)	CIA	mouse; not bovine, not canine, not chicken, not human	205603	100 µg	€390
Anti-E-Cadherin (ECCD-2) Rat mAb (ECCD-2)	FS, IB; NOT PS	human, mouse	205604	100 µg	€426
Anti-E-Cadherin Mouse mAb (HECD-1)	CIA, CT, ELISA, FC, IB, IF, IH	human	205601	100 µg	€426
Anti-E-Cadherin Mouse mAb (SHE78-7)	CIA, ELISA, FC, FS, IB, IF, PS	human; not rat	205602	100 µg	€390
Anti-M-Cadherin Mouse mAb (12G4)	ELISA, IB, IF	mouse; not human	205610	100 µg	€397
Anti-N-Cadherin (808 - 827) Rabbit pAb	IB, PS	bovine, chicken, human, mouse, rat	205606	400 µg	€422
Anti-N-Cadherin Rat mAb (NCD-2)	CIA, FS, IB, IH; NOT PS	chicken, <i>Xenopus</i>	205605	100 µg	€390
Anti-P, E, N-Cadherin (853-866) Rabbit pAb	IB, PS	chicken, human, mouse	205609	100 µg	€390
Anti-CAIN (1881-2173) Rabbit pAb	ELISA, IB, IP	human, rat	207002	100 µl	€346
Anti-Calbindin D-28K (Ab-1) Rabbit pAb	FS, VS	baboon, human, monkey, rat	PC253L	100 µl	€409
Anti-Calceineurin B (Ab-1) Rabbit pAb	IB, IP	bovine, human, mouse	PC359	100 µl	€319
Anti-Calcitonin Gene Related Peptide (CGRP) Rabbit pAb	FFS, FS, IF, PS	mouse, rat	PC205L	100 µl	€338
Anti-Calmodulin (Ab-1) Mouse mAb (2D1)	ELISA, IB, IC	bovine, chicken, rat	208693	200 µl	€384
Anti-Calmodulin (Ab-2) Mouse mAb (6D4)	IB, IC	bovine, chicken, rat	208696	200 µl	€384
Anti-Calnexin, C-Terminal (575-593) Rabbit pAb	IB, IC, IP	bovine, canine, guinea pig, hamster, human, monkey, mouse, rabbit, rat, sheep	208880	50 µl	€340

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Calnexin, N-Terminal (50-68) Rabbit pAb	IB, IC, IP	bovine, canine, guinea pig, hamster, human, monkey, mouse, porcine, rabbit, rat, sheep	208882	50 µl	€340
Anti-Calpain-1, C-Terminal Rabbit pAb	IB	human, mouse, rat	208753	100 µg	€378
Anti-Calpain-1, Domain I Rabbit pAb	IB	human, mouse, rat	208751	100 µg	€378
Anti-Calpain-1, Domain II Mouse mAb (2H2A7C2)	IB, IC	bovine, human, rat	208727	100 µl	€384
Anti-Calpain-1, Domain III Mouse mAb (9A4H8D3)	IB, IC	bovine, human, porcine, rat	208728	100 µl	€384
Anti-Calpain-1, Large Subunit Rabbit pAb	ELISA, IB	human	208732	100 µl	€301
Anti-Calpain-1/2, Small Subunit Mouse mAb (28F3)	ELISA, IB, IP	bovine, human; not mouse, not rat	208730	100 µg	€340
Anti-Calpain-2, Domain I Rabbit pAb	IB	human, mouse, rat	208754	100 µg	€378
Anti-Calpain-2 Domain III Rabbit pAb	ELISA, IB, IH, IP	rat	208737	100 µg	€397
Anti-Calpain-2, Domain III Rabbit pAb	IB	human, mouse, rat	208755	100 µg	€378
Anti-Calpain-2, Domains III/IV Mouse mAb (107-82)	IB, IC, IF	bovine, human, porcine, rat	208729	100 µl	€384
Anti-Calpain-2, Domain IV Rabbit pAb	IB	human, mouse, rat	208756	100 µg	€378
Anti-Calpain-3, Domain I Rabbit pAb	IB	human	208757	100 µg	€378
Anti-Calpain-3, Domain III (Insert #2) Rabbit pAb	IB	human	208759	100 µg	€378
Anti-Calpastatin, Domain IV Mouse mAb (1F7E3D10)	IB, IF	bovine, human, rat	208903	100 µl	€384
Anti-Calpastatin Mouse mAb (CSL1-5)	ELISA, FS, IB	human	208905	100 µg	€413
Anti-Calreticulin (405-417) Rabbit pAb	IB, IC, IP	bovine, canine, hamster, human, monkey, mouse, porcine, rat	208910	50 µl	€290
Anti-Calreticulin Mouse mAb (FMC75)	EIA, FC, IB, IC, IH	human, monkey	208912	50 µg	€340
Anti-Calretinin Rabbit pAb	FS, VS	rat	PC254L	100 µl	€409
Anti-cAMP Rabbit pAb	RIA		116820	1 Set	€385
Anti-CAND1 (1211-1230) Rabbit pAb	IB, IF, IP	human, mouse, rat	PC745	100 µl	€384
Anti-Cannabinoid Receptor CB ₁ (1-77) Rabbit pAb	IB, IC, PS	chicken, fish, human, monkey, mouse, rat, salamander	209550	100 µl	€362
Anti-Cannabinoid Receptor CB ₂ (1-32) Rabbit pAb	IC	human, rat	209554	100 µl	€362
Anti-Cannabinoid Receptor CB ₂ (1-33) Rabbit pAb	IB, IC	human; not mouse, not rat	209552	100 µl	€362
Anti-Cannabinoid Receptor CB ₂ (20-33) Rabbit pAb	IB, IH	human	216407	200 µg	€384
Anti-Capsaicin Receptor (Ab-1) (824-838) Rabbit pAb	FS, IB, IF, PS	mouse, rat	PC420	5 µg 25 µg	€86 €338
Anti-Capsaicin Receptor (Ab-2) (1-21) Rabbit pAb	FS	rat	PC547	50 µg	€253
Anti-Casein Kinase 1δ (310-327) Goat pAb	PS	human	ST1006	100 µl	€319
Anti-Casein Kinase II, α-Subunit (70-89) Rabbit pAb	IB	human, mouse, rat	231573	100 µg	€384
Anti-Casein Kinase II, α-Subunit Mouse mAb (1AD9)	IB, IC	bovine, hamster, human, mouse	218703	100 µg	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW


Anti-α-catenin Mouse Antibody (1G5), Cat No. CA1030

Lane 1: Detection of human α-catenin by immunoblotting. Sample: Whole cell lysate (25µg) from HeLa cells. Primary Abs: anti-α-catenin mouse mAb (1G5) (Cat. No. 1030) (1:1000). Detection: Chemiluminescence

Lane 2: Detection of human α-catenin by immunoprecipitation followed by immunoblotting: Antibody for immunoprecipitation: anti-α-catenin mouse mAb (1G5) (Cat. No. CA1030) (1:100/500 µg total protein). Immunoblotting conditions: same as lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-Casein Kinase II Substrate Mouse mAb (9F4)	ELISA	human	231576	100 µg	€363
Anti-Casein Kinase II, β-Subunit Mouse mAb (6D5)	IB, IC	bovine, hamster, human, mouse	218712	50 µg	€319
Cleaved Caspase Antibody Sampler Kit	IB		AP1026	1 kit	€653
Anti-Caspase-1 Mouse mAb (14F468)	IB, PS	human, mouse	AP1043	50 µg	€198
Anti-Caspase-1 (31-45) Rabbit pAb	IB, PS	human	AP1044	50 µg	€176
Anti-Pan Caspase-1 (390-404) Rabbit pAb	IB, PS	human; not mouse	PC84	100 µg	€319
Anti-Caspase-2 (1-14) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC107	100 µg	€319
Anti-Caspase-3 (Ab-2) Mouse mAb (10C1.C9)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	AM34	100 µg	€338
Anti-Caspase-3 (Ab-3) Mouse mAb (AM1.4.1-1B)	IB; NOT FS, NOT IP, NOT PS	human	AM39	100 µg	€248
Anti-Caspase-3 (Ab-4) Mouse mAb (AM1.31-11)	IB, PS; NOT FS, NOT IP	human, mouse	AM65	100 µg	€338
Anti-Cleaved Caspase-3 (Asp ¹⁷⁵) Rabbit pAb	FC, IB, IC, PS	human, mouse, rat	AP1027	50 µl	€231
Anti-Caspase-3 Rabbit pAb	IB	bovine, canine, hamster, human, monkey, mouse, porcine, rabbit, rat, yeast	235412	100 µl	€404
Anti-Caspase-3, Cleaved (Ab-2) Rabbit pAb	FS, IB, IC, PS	human, mouse, rat	PC679	50 µg	€319
Anti-Caspase-4 (Ab-1) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC109	100 µl	€319
Anti-Caspase-6 (250-264) Rabbit pAb	IB	bovine, canine, hamster, human, monkey, mouse, porcine, rabbit, rat, sheep	218774	100 µg	€384
Anti-Cleaved Caspase-6 (Asp ¹⁶²) Rabbit pAb	IB	human, mouse, rat	AP1012	50 µl	€278
Anti-Caspase-7 (Ab-1) (186-198) Rabbit pAb	IB; NOT IF, NOT IP	human	PC334	100 µg	€338
Anti-Active Caspase-7 Rabbit pAb	IB, IF, IP	human, mouse, rat	PC520	50 µl	€338
Anti-Caspase-8 (Ab-1) (236-247) Rabbit pAb	IB; NOT IF, NOT IP	human	PC335	100 µg	€338
Anti-Caspase-8 (Ab-3) Mouse mAb (1-3)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	AM46	100 µg	€338
Anti-Caspase-8 (Ab-3) Mouse mAb (1-3)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	AM46T	10 µg	€77
Anti-Cleaved Caspase-8 (Asp ³⁸⁴) Mouse mAb (11G10)	IB	human	AP1013	100 µl	€413
Anti-Caspase-9 (1-134) Rabbit pAb	IB	bovine, chicken, guinea pig, human, monkey, mouse, rat, sheep	218794	100 µg	€413
Anti-Caspase-9 (Ab-2) Mouse mAb (1-2)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human; not mouse, not rat	AM47	100 µg	€338
Anti-Caspase-9 (Ab-2) Mouse mAb (1-2)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human; not mouse, not rat	AM47T	10 µg	€77
Anti-Cleaved Caspase-9 (Asp ³¹⁵) Rabbit pAb	IB, IP	human	AP1014	50 µl	€278
Anti-Cleaved Caspase-9 (Asp ³⁵³) Rabbit pAb	FFS, IB, IC	rat	AP1015	50 µl	€278
Anti-Cleaved Caspase-9 (Asp ³⁵³) Rabbit pAb	IB, IC	mouse	AP1028	50 µl	€231

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-Caspase-1 Mouse mAb (14F468) Cat. No. AP1043

Detection of human caspase-1 by staining paraffin sections. Sample: Colon carcinoma, formalin-fixed (panels A and B). Primary antibody: Anti-Caspase-1 Mouse mAb (14F468) (Cat. No. AP1043). Detection: DAB with hematoxylin counterstain. B is higher magnification than A.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Caspase-9, Cleaved (Ab-2) Rabbit pAb	IB, IC	human	PC680	50 µg	€319
Anti-Caspase-10 (Ab-1) (221-234) Rabbit pAb	IB	human	PC332	100 µg	€324
Anti-Caspase-12 (Ab-1) (100-116) Rabbit pAb	IB	human, mouse, rat	PC557	100 µg	€338
Anti-Caspase-12 (Ab-1) Rabbit pAb	IB	human, mouse, rat	PC557T	10 µg	€77
Anti-Caspase-12 (Ab-2) (2-17) Rabbit pAb	IB	human, mouse, rat	PC558	100 µg	€338
Anti-Caspase-14 (Ab-1) Mouse mAb (70A1426)	IB; NOT FS, NOT IP, NOT PS	mouse	AM64	100 µg	€346
Anti-Catalase Rabbit pAb	ELISA, IB, IEP	human, mouse	219010	1 ml	€124
Anti-α-Catenin Mouse mAb (1G5)	IB, IP	human	CA1030	50 µl	€135
Anti-β-Catenin C-Terminal Mouse mAb (10H8)	ELISA, IB, IC, IP	human, mouse, rat	219356	1 Set	€390
Anti-β-Catenin, Dephospho Mouse mAb (7A7)	ELISA, IB, IP	human, mouse, rat	219349	100 µg	€413
Anti-β-Catenin, Dephospho Mouse mAb (8E4)	ELISA, IB, IP	human, mouse, rat	219350	1 Set	€409
Anti-β-Catenin Mouse mAb (9G10)	ELISA, IB, IC, IP	human, mouse, rat	219353	1 Set	€390
Anti-β-Catenin N-Terminal Mouse mAb (7D11)	ELISA, IB, IP	human, mouse, rat	219357	1 Set	€390
Anti-Cathepsin B (Ab-1) Mouse mAb (CA10)	IB, IF, IP	bovine, human, porcine, not rat	IM27L	20 µg 100 µg	€86 €346
Anti-Cathepsin B (Ab-3) Rabbit pAb	IB, PS	human	PC41	100 µg	€319
Anti-Cathepsin B Rabbit pAb	ELISA, IB	human	219408	1 ml	€330
Anti-Cathepsin D (Ab-1) Mouse mAb (BC011)	IB, PS	human	IM03	100 µg	€338
Anti-Cathepsin D (Ab-2) Rabbit pAb	FS, IB, PS	human	IM16	100 µg	€338
Anti-Pro-Cathepsin D (Ab-1) Rabbit pAb	IB, PS	human; not bovine, not porcine, not rat	IM04L	100 µg	€319
Anti-Cathepsin D Rabbit pAb	ELISA, IEP	human	219361	1 ml	€340
Anti-Cathepsin F Rabbit pAb	IB	human, mouse	219359	100 µl	€319
Anti-Cathepsin G Rabbit pAb	ELISA, IB, IEP	human	219358	1 ml	€122
Anti-Cathepsin K Mouse mAb (182-12G5)	IB	human	IM55L	100 µg	€338
Anti-Cathepsin K Rabbit pAb	IB	human	219386	100 µl	€244
Anti-Cathepsin L Rabbit pAb	IB	human	219387	100 µl	€244
Anti-Cathepsin S Mouse mAb (2B6)	DB, ELISA, IB	human	IM1003	100 µg	€338
Anti-Cathepsin S Rabbit pAb	IB	human	219384	100 µl	€231
Anti-Cathepsin Z Rabbit pAb	IB	human	219378	100 µl	€319
Anti-Caveolin-2 (1-19) Rabbit pAb	IB, IP	human, mouse, rat	219383	100 µg	€362
Anti-Caveolin-2 (18-31) Rabbit pAb	IB	human	CB1004	50 µg	€390
PhosphoDetect™ Anti-Caveolin-2 (pSer ²³) Rabbit pAb	IB	human	CB1005	50 µg	€353
PhosphoDetect™ Anti-Caveolin-2 (pSer ³⁶) Rabbit pAb	IB	human	CB1006	50 µg	€353
Anti-CBFβ Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human, mouse	PC288	100 µg	€319
Anti-c-Cbl Mouse mAb (5C4)	ELISA, IB	human, mouse	219439	100 µg	€409

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-Caspase-1 (31-45) Rabbit pAb Cat. No. AP1044

Detection of human caspase-1 by staining paraffin sections. Samples: Lymph node, formalin-fixed. Primary antibody: Anti-Caspase-1 (31-45) Rabbit pAb (Cat. No. AP1044). Detection: DAB with hematoxylin counterstain. B is higher magnification than A.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-CCR1 (7-24) Rabbit pAb	ELISA, IB, IC	human	227020	100 µl	€413
Anti-CCR2 (Ab-1) Goat pAb	FC, IB, IC, PS	human	227001	100 µg	€409
Anti-CCR4 (17-39) Goat pAb	ELISA, IB, IC, PS	human	227003	100 µg	€409
Anti-CCR5 (6-20) Rabbit pAb	IB	human	PC559	100 µg	€338
Anti-CCR7 (2-22) Goat pAb	ELISA, IB, IC	mouse	227006	100 µg	€409
Anti-CCR10 (4-31) Goat pAb	ELISA, IB, IC	mouse	227011	100 µg	€445
Anti-CCR10 (9-31) Goat pAb	ELISA, IB, IC	human	227010	100 µg	€409
Anti-CD3 Mouse mAb (UCHT1)	FC, FS, IP, TCELL	human	217570	100 µg	€158
Anti-CD4 Mouse mAb (QS4120)	2D PAGE, ELISA, FC, FS	human	217575	100 µg	€198
Anti-CD4 Mouse mAb (QS4120) Fluorescein Conjugate	ELISA, FC	human	217577	120 T	€231
Anti-CD8 Mouse mAb (UCHT4)	FC, FS, IP	human	217580	100 µg	€198
Anti-CD10 Mouse mAb (SN5c)	FC, FS, IP; NOT PS	human	OP161	100 µg	€338
Anti-CD11a Mouse mAb (38)	CIA, FC, FS, IP	human	217640	100 µg	€273
Anti-CD11b Mouse mAb (ICRF44)	CIA, FC, FS, IP	human	217642	100 µg	€260
Anti-CD18 Mouse mAb (IB4)	CIA, FC, FS, IP	human	217660	100 µg	€229
Anti-CD28 Mouse mAb (ANC28.1/5D10)	FC, FS, IP, TCELL	human	217669	100 µg	€314
Anti-CD29 Mouse mAb (4B7R)	FC, FS, IP, PS	human	217648	100 µg	€314
Anti-CD34 Mouse mAb (QBEnd/10)	FC, PS	human; not rat	OP164	100 µg	€338
Anti-CD40 Mouse mAb (EA-5)	FC, FS, IP, PS	human	217590	100 µg	€252
Anti-CD40L Mouse mAb (24-31)	FC, FS, IB, IP, NT	human	217595	100 µg	€351
Anti-CD40L Mouse mAb (24-31) Fluorescein Conjugate	FC, FS, NT	human	217597	120 T	€450
Anti-CD44 Rat mAb (A020)	ELISA, FC, IB, IF, IH, IP, RIA	bovine, canine, human, mouse, porcine, rabbit	217594	100 µl	€377
Anti-CD44v10 Rabbit pAb	ELISA, FC, FS, IB, IH, IP, PS, RIA	bovine, canine, human, mouse, porcine, rabbit	217685	100 µg	€381
Anti-CD44v3 Rabbit pAb	ELISA, FC, IB, IF, IH, IP, RIA	bovine, canine, human, mouse, porcine, rabbit	217599	100 µl	€390
Anti-CD44v6 Rabbit pAb	ELISA, FC, FS, IB, IP, PS, RIA	bovine, canine, human, mouse, porcine, rabbit	217604	100 µl	€390
Anti-CD45 Mouse mAb (Bra55)	ELISA, FC, FS, IP, NT, PS	human; not Old World monkeys	OP166	100 µg	€338
Anti-CD47 Mouse mAb (2D3)	APOP, FC, IF, IP	human	OP167	100 µg	€338
Anti-CD49d Mouse mAb (BU49)	CIA, FC, FS, IP	human	217600	100 µg	€240
Anti-CD49d Mouse mAb (BU49) Fluorescein Conjugate	FC, FS	human	217602	120 T	€283
Anti-CD49f Mouse mAb (BQ16)	FC, FS, IB, IP	human	217657	100 µg	€314
Anti-CD54 Mouse mAb (15.2)	FC, FS, IB, IP, NT	human	217605	100 µg	€278
Anti-CD54 Mouse mAb (8.4A6)	CIA, ELISA, FC, FS, IP	human	217677	100 µg	€350
Anti-CD55 Mouse mAb (55C02)	FC, IF	human	OP168	100 µg	€338
Anti-CD62L Mouse mAb (LAM1-116)	CIA, FC, IH, IP	human	217683	100 µg	€314
Anti-CD63 Mouse mAb (NK1/C-3)	IB, PS	human	OP171	100 µg	€338
Anti-CD64 Mouse mAb (10.1)	FC, FS, IB, IP, NT	human	217620	100 µg	€244
Anti-CD69 Mouse mAb (HP-4B3)	FC, FS, IP	human	217850	100 µg	€299
Anti-CD86 Mouse mAb (BU63) Fluorescein Conjugate	FC, FS, NT	human	217632	120 T	€365
Anti-CD95L Mouse mAb (Alf-2.1a) Fluorescein Conjugate	FC, NT	human	217952	120 T	€332
Anti-CD106 Mouse mAb (1.G11B1)	FC, FS, IP	human	217635	100 µg	€301
Anti-CD106 Mouse mAb (1.G11B1) Fluorescein Conjugate	FC, FS, IP	human	217637	120 T	€377
Anti-CD117 Mouse mAb (57A5)	FC, FS, IP	human	218150	100 µg	€299
Anti-CD162 Mouse mAb (PL1) R-Phycoerythrin Conjugate	CIA, FC, FS	human	218131	120 T	€469


10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-CD68 Mouse mAb (KP1)	IH	feline, human, monkey, rat	CB1014	500 µl	€253
PhosphoDetect™ Anti-Cdc2 (pTyr ¹⁵) Rabbit pAb	DB, IB	rat	219437	100 µl	€363
Anti-Cdc2/Cdk1 (Ab-1) Rabbit pAb	IB, IF, IP	human, mouse, rat	PC25	100 µg	€207
Anti-Cdc2/Cdk1 (Ab-2) Mouse mAb (cdc2.1)	IB, IP	human, mouse, rat, yeast	CC01	100 µg	€206
Anti-Cdc2/Cdk1 (Ab-3) Mouse mAb (A17.1.1)	FC, FS, IB, IF, IP	chicken, guinea pig, human, mouse, rat, woodchuck, <i>Xenopus</i>	CC16	100 µg	€206
Anti-Cdc6 (Ab-3) Mouse mAb (DCS 180.3)	IB, PS; NOT FS, NOT IF	human, mouse, rat (weakly)	CC30	100 µg	€206
Anti-Cdc25A Rabbit pAb	IB, IP	human	PC733	100 µl	€384
Anti-Cdc25B (Ab-1) Mouse mAb (DCS.162.2)	IB; NOT FS, NOT PS	human, mouse (weakly)	CC25	100 µg	€206
Anti-Cdc25C (Ab-5) Mouse mAb (25C07)	IP, PS; NOT IB	human; not mouse, not rat	CC42	100 µg	€338
Anti-Cdh1 (Ab-2) Mouse mAb (DH01)	IB, IP	human	CC43	100 µg	€338
PhosphoDetect™ Anti-Cdk1 (pThr ¹⁴ /pTyr ¹⁵) Rabbit pAb	IB, IP	human, mouse, rat	219441	10 T	€373
PhosphoDetect™ Anti-Cdk1 (pTyr ¹⁵) Rabbit pAb	IB, IC, IP	human, rat	219440	100 µl	€353
PhosphoDetect™ Anti-Cdk1 Substrate (Ab-1) Rabbit pAb	IB, PS	human	PC456	25 µl	€282
Anti-Cdk2, C-Terminal (283-298) Rabbit pAb	IB, IP	human, mouse, rat	219443	100 µg	€384
Anti-Cdk4 (Ab-2) Mouse mAb (DCS32.1)	IB; NOT IF, NOT IP	human, mouse	CC19L	100 µg	€206
Anti-Cdk5, C-Terminal (268-283) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	219449	100 µg	€397
Anti-Centrosome Mouse mAb (6C6)	IB, IC, IF	bovine, human	CB1003	100 µg	€319
Anti-CETP Mouse mAb (ATM192)	ELISA, IB, IP	human, rabbit	219656	100 µl	€384
Anti-CHIP Rabbit pAb	IB, IF, IP	<i>Drosophila</i> , hamster, human, monkey, mouse, rat	PC711	100 µl	€351
Anti-Chk1 (1-474) Sheep pAb	IB	human	PC423	250 µg	€338
PhosphoDetect™ Anti-Chk1 (pSer ³¹⁷) Rabbit pAb	IB, IC, PS	human, monkey, mouse, rat	DR1025	50 µl	€199
Anti-Chk2 (Ab-2) Mouse mAb (2CHK01)	IB, IF, IP	human	CC44	100 µg	€338
PhosphoDetect™ Anti-Chk2 (pThr ⁶⁸) Rabbit pAb	FC, IB, IC, IH, IP	human, monkey	DR1026	50 µl	€199
Anti-Cholecystokinin-8 (26-33) Rabbit pAb	FS, IF	rat	PC206L	100 µl	€319
Anti-Cholera Toxin, B-Subunit Goat pAb	IB, IH, RID		227040	100 µl	€64
Anti-Choline Acetyltransferase Rat mAb	FFS, FS, VS	canine, feline, ferret, guinea pig, monkey, mouse, rat, turtle	NB05L	100 µl	€543
Anti-β-Chorionic Gonadotropin Mouse mAb (P3G7)	EIA, RIA	human	230746	1 mg	€168
Anti-CIN85, C-Terminal (646-665) Rabbit pAb	ELISA, IB, IH, IP	human, mouse, porcine, rat	231006	100 µl	€397
Anti-Citrulline Rabbit pAb	ELISA	a broad range of species	231246	50 µl	€294
Anti-Caspin Rabbit pAb	IB, IP	human	DR1006	100 µl	€397
Anti-Clathrin Heavy Chain Mouse mAb (X22)	IB, IF, IP	bovine, human, mouse, rat	CP45	100 µl	€315
Anti-CLOCK (11-30) Rabbit pAb	IB, IH	hamster, human, mouse, rat	233170	200 µl	€362
Anti-Coagulation Factor XIII, S-Subunit Rabbit pAb	ELISA, IB	human	233503	1 ml	€397
Collagen Antibody Sampler Kit I			ASK19	1 each	€270
Anti-Collagen, Type I, Cyanogen Bromide Fragments Rabbit pAb	IB	bovine (weakly), human (weakly), mouse, rat	234168	500 µl	€319
Anti-Collagen Type I Mouse mAb (I-8H5)	FS, IB	human, rabbit	CP17	20 µg 100 µg	€86 €338
Anti-Collagen, Type I Rabbit pAb	FS, IB, IF, IP	human, mouse, rat	234167	500 µl	€319
Anti-Collagen, Type I/III, Cyanogen Bromide Fragments Rabbit pAb	FS, IB, IF, IP	human, mouse, rat	234169	500 µl	€319
Anti-Collagen, Type II (Ab-1) Mouse mAb (II-4C11)	FS, IB, PS	bovine, human, rabbit, rat	CP18L	100 µg	€338
Anti-Collagen, Type II Rabbit pAb	FS, IB, IP	human, mouse, rat	234187	500 µl	€319
Anti-Collagen, Type II, Cyanogen Bromide Fragments Rabbit pAb	FS, IB, IF, IP; NOT PS	human, mouse, rat	234188	500 µl	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Collagen, Type III Mouse mAb (III-53)	FS, IB	human, rabbit	CP19	20 µg 100 µg	€86 €338
Anti-Collagen, Type III, Cyanogen Bromide Fragments Rabbit pAb	FS, IB, IF, IP	human, marmoset, mouse, rat	234189	500 µl	€319
Anti-Collagen, Type IV (Ab-1) Mouse mAb (IV-4H12)	IB	human	CP20	100 µg	€338
Anti-Collagen, Type IV (Ab-2) Mouse mAb (CIV22)	FS, PS	bovine, human; not canine, not porcine, not rat	CP56	500 µl	€339
Anti-Collagen, Type IX Rabbit pAb	FS, IB, IP	bovine, human, mouse, rat	234194	500 µl	€384
Anti-Collagen, Type VI Rabbit pAb	FS, IB, IF	human, marmoset, rat	234191	500 µl	€319
Anti-Collagen, Type VII Rabbit pAb	FS, IB, IF, IP	human, marmoset	234192	500 µl	€319
Anti-Collagen, Type X Rabbit pAb	FS, IB, IP	human, mouse, rat	234196	500 µl	€384
Anti-Complement 1q Goat pAb	IH, OT	baboon, canine, guinea pig, hamster, human, mouse, rabbit, rat; not bovine, not chicken, not sheep	234390	1 ml	€85
Anti-Complement 2 Goat pAb	IEP, OT	human	234365	1 ml	€231
Anti-Complement 3 Goat pAb	FC, IB, IC, IEP	baboon, guinea pig, hamster, horse, human, mouse, rabbit, rat; not bovine, not canine, not chicken, not feline, not sheep	204869	1 ml	€92
Anti-Complement 3a Rabbit pAb	ELISA, IB, IEP	human	204859	1 ml	€150
Anti-Complement 4 Goat pAb	IEP, IH	human	204894	1 ml	€92
Anti-Complement 5a Rabbit pAb	ELISA, IB, OT	human; not rat	204889	1 ml	€239
Anti-Complement 5b-9 Rabbit pAb	FS, IF, RID	human, mouse	204903	1 mg	€278
Anti-Complement 9 Goat pAb	IEP, RID	human	204912	1 ml	€92
Anti-β-COP (496-513) Rabbit pAb	IB, IC, IF	mouse, rat	PC175	100 µl	€304
PhosphoDetect™ Anti-Cortactin (pTyr ⁴²¹) Rabbit pAb	IB, IC	mouse; not human	CB1000	10 T	€373
Anti-Corticotropin Releasing Factor Rabbit pAb	FS, IF	rat, sheep	PC208L	100 µl	€319
Anti-Cortisol Mouse mAb (F4P1A3)	EIA, ELISA, RIA	a broad range of species	235042	1 mg	€192
Anti-COX2, C-Terminal Rabbit pAb	ELISA, IB	bovine, human, mouse, rat, sheep	236002	50 µg	€319
Anti-COX2 Mouse mAb (5E10/D10)	ELISA, IB	human, mouse	236004	100 µg	€378
Anti-CREB (123-136) Rabbit pAb	IB	rat	238465	100 µg	€384
PhosphoDetect™ Anti-CREB (pSer ¹³³) Mouse mAb (10E9)	ELISA, IB	human, mouse	238467	1 set	€385
PhosphoDetect™ Anti-CREB (pSer ¹³³) Rabbit pAb	IB	human, rat	238466	100 µl	€350
Anti-CREB Rabbit pAb	IB, IC, IP	finch, human, mouse, rat	238461	100 µl	€340
Anti-CRM1 Rabbit pAb	IB, IP	human, <i>Xenopus</i>	ST1100	50 µg	€145

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-CtIP Rabbit pAb Cat. No. AP1046

Lane 1: Detection of CtIP by immunoblotting. Sample: Whole cell lysate (50 µg) from MCF7 cells. Primary antibody: Anti-CtIP Rabbit pAb (Cat. No. AP1046) (1 µg/ml). Detection: Chemiluminescence. Lane 2, 3: Detection of CtIP by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from MCF7 cells. Antibodies for IP: Anti-CtIP Rabbit pAb (Cat. No. AP1046) (3 µg/mg protein) (lane 2) and Rabbit IgG (3 µg/mg protein) (lane 3). Immunoblotting conditions: same as above.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti- α B-Crystallin Peptide (163-175) Rabbit pAb	IB, IC	bovine, human, mouse, rat	238702	50 μ l	€294
NEW Anti-CtIP Rabbit pAb	IB, IP	human	AP1046	50 μ g	€145
Anti-CSF-1 Receptor/c-fms (Ab-2) Rat mAb (3-4A4-E4)	FC, IP	human	GR12L	100 μ g	€338
Anti-CSN3 Rabbit pAb	IB	human	ST1043	100 μ g	€290
Anti-CSN4 Rabbit pAb	IB	human	ST1044	100 μ g	€290
Anti-Cu ²⁺ -Oxidized Low Density Lipoprotein Rabbit pAb	ELISA, IH	a broad range of species	428033	100 μ l	€319
Anti-Cullin-1 (742-752) Rabbit pAb	FS, IB, IP, PS	human, mouse	PC441	500 μ g	€338
Anti-CX ₃ CR1 (2-21) Rabbit pAb	IB	human, mouse, rat	PC443	100 μ g	€338
Anti-CXCR4 (Ab-1) (182-196) Rabbit pAb	IB	human	PC389	100 μ g	€319
Anti-CXCR4 (Ab-2) (1-14) Rabbit pAb	IB	human	PC390	100 μ g	€319
Anti-Cyclin A (Ab-2) Mouse mAb (BF683)	IB, IP, PS	human; not mouse, not rat	CC02	100 μ g	€225
Anti-Cyclin A (Ab-3) Mouse mAb (E23)	FC, IB, IF; NOT FS, NOT IP, NOT PS	bovine, human, mouse, rat	CC17	100 μ g	€225
Anti-Cyclin A (Ab-5) Rabbit pAb	ELISA, IB, PS	human, mouse, rat	PC132	100 μ l	€307
Anti-Cyclin B1 (Ab-2) Mouse mAb (2H1-H6)	IF, IP; NOT IB, NOT PS	human; not mouse, not rat	CC03	100 μ g	€225
Anti-Cyclin B1 (Ab-3) Rabbit pAb	ELISA, IB	human, mouse, rat	PC133	100 μ l	€307
Anti-Cyclin D1 (Ab-1) Rat mAb (17A6-4)	IB, IP	human, mouse; not rat	CC04	100 μ g	€225
Anti-Cyclin D1 (Ab-2) Mouse mAb (HD11)	IB, IP	human; not mouse, not rat	CC11	100 μ g 200 μ g	€225 €378
Anti-Cyclin D1 (Ab-3) Mouse mAb (DCS-6)	FC, FS, IB, IF, IP, NT, PS	human, mouse, rat	CC12	100 μ g	€225
Anti-Cyclin D1 (Ab-3) Mouse mAb (DCS-6) Fluorescein Conjugate	FC, IF	human, mouse, rat	CC12F	100 μ g	€351
Anti-Cyclin D3 (Ab-2) Mouse mAb (DCS-22)	IB, PS; NOT FS	human	CC13	100 μ g	€225
Anti-Cyclin E (Ab-1) Mouse mAb (HE12)	IB; NOT IP, NOT PS	human; not mouse, not rat	CC05	20 μ g 100 μ g	€86 €225
Anti-Cyclin E (Ab-2) Rabbit pAb	ELISA, IB	human, mouse, rat	PC135	100 μ l	€307
Anti-Cyclin E (Ab-3) (381-395) Rabbit pAb	IB, IP	human	PC438	500 μ g	€338
Anti-Cyclin L Rabbit pAb	IB, IP	human	PC705	50 μ l	€378
Anti-Cyclophilin 40 Rabbit pAb	IB	human, rat	239780	100 μ l	€381
Anti-Cyclophilin A, Human, HuCAL [®] Fab	ELISA, IB	bovine, human	ST1071	50 μ g	€189
Anti-Cyclophilin A Rabbit pAb	IB	human, mouse, rat	239778	100 μ l	€384
Anti-Cyclophilin D Mouse mAb (E11AE12BD4)	IB, IC	bovine, human, rat	AP1035	50 μ g	€294
Anti-Cyclophilin D (42-54) Rabbit pAb	IB	canine, hamster, human, mouse, rat	239784	100 μ g	€362
Anti-Cytochrome c (Ab-1) Sheep pAb	IB, IF, IP	canine, human, rabbit, rat	PC323	50 μ g 100 μ g	€231 €338
NEW Anti-Cytochrome c Mouse mAb (6H2.B4)	FC, IC, IP; NOT IB	human, mouse, rat	AP1030	50 μ g	€183

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TECELL:** t cell co-stimulation; **VS:** vibratome sections


**NEW****Anti-Cyclophilin D Mouse mAb (E11AE12BD4)****Cat. No. AP1035**

Detection of human cyclophilin D by immunocytochemistry. Sample: Cultured fibroblasts. Primary antibody: Anti-Cyclophilin D Mouse mAb (E11AE12BD4) (Cat. No. AP1035) (5 μ g/ml). Detection: fluorescence.

NEW

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Cytochrome c Mouse mAb (7H8.2C12)	FC, IB, IC, PS; NOT IP	horse, human, mouse, rat	AP1029	50 µg	€183
Anti-Cytochrome p450 1A1 Mouse mAb (1-7-1)	NT	human, rat	250302	50 µl	€378
Anti-Cytochrome p450 2E1 Rabbit pAb	IB	human	250306	100 µl	€355
Anti-Cytochrome p450 3A4 Rabbit pAb	IB	human	250308	100 µl	€355
Anti-Cytokeratin 18 (Ab-2) Mouse mAb (DC10)	FS, IB, PS	human	IF14	100 µg	€338
Anti-Cytokeratin 18 Rabbit pAb	IB	human	AP1021	50 µg	€176
Anti-Cytokeratin 19 Mouse mAb (A53-B/A2.26)	FS, IB, IF, PS	guinea pig (weakly), human, mouse (weakly), rat (weakly)	IF15	100 µg	€338
Anti-Cytokeratin 20 Mouse mAb (Ks20.8)	FS, PS	human	CP67	100 µl	€253
Anti-Pan-Cytokeratin Mouse mAb (B311.1)	ELISA, FS, IB, IP, PS	bovine, human, mouse, rat	250400	500 µl	€319
PhosphoDetect™ Anti-DARPP32 (pThr ⁷⁵) Rabbit pAb	DB, IB	human, mouse	PS1005	100 µl	€385
Anti-Daxx (722-740) Rabbit pAb	IB	human	PC394	100 µg	€319
Anti-DCC Mouse mAb (AF5)	IB, IF, IP	human	OP45	100 µg	€346
Anti-DC-SIGN (Ab-1) (384-404) Rabbit pAb	IB	human	PC650	100 µg	€338
Anti-DDAH1 (272-284) Goat pAb	IB	human	PC716	100 µg	€397
Anti-DDAH2 (273-285) Goat pAb	IB	human, mouse	PC717	100 µg	€397
Anti-DDB1 (1128-1140) Goat pAb	IB	human, mouse	PC718	100 µg	€397
Anti-Decorin Goat pAb	ELISA, IB, IH	human	PC673	100 µg	€408
Anti-Desmin (Ab-1) Mouse mAb (DE-B-5)	FS, IB, IC, PS	human	IF02L	100 µg	€324
Anti-Desmin (Ab-2) Mouse mAb (D33)	FS, PS	baboon, bovine, canine, chicken, feline, hamster, human, monkey, rat	IF11	500 µl	€339
Anti-DFF40 (329-338) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC380	100 µg	€319
Anti-DFF45 (21-33) Rabbit pAb	IB; NOT PS	human, mouse	PC366	100 µg	€338
Anti-Dipeptidylpeptidase IV Rabbit pAb	FS, IB, IP	human	IM1004	100 µl	€390
Anti-DNA Methyltransferase 3a (86-100) Rabbit pAb	IB	human, mouse	317282	100 µg	€346
Anti-DNA Methyltransferase 3b (1-14) Rabbit pAb	IB	human, mouse	317284	100 µg	€346
Anti-DnaK, <i>E. coli</i> Mouse mAb (8E2/2)	EIA, ELISA, IB, IP	<i>E. coli</i>	323110	50 µg	€304
Anti-DNA-PK (Ab-1) (2015-2134) Rabbit pAb	IB, IF, IP, PS	horse, human, mouse, rat, <i>Xenopus</i>	PC127	100 µg	€319
Anti-DNA-PK (Ab-2) Mouse mAb (18-2)	FS, IB, IF, IP, NT, PS	human	NA57	100 µg	€338
Anti-DNase γ (Ab-1) Mouse mAb (hg302)	IH, IP	human, mouse, rat	AM78	100 µg	€384
Anti-DNase γ (Ab-2) Mouse mAb (hg303)	IB	human	AM79	100 µg	€384
Anti-DOPA Decarboxylase Rabbit pAb	IB, IC	bovine, canine, human, rabbit, sheep; not rat	324381	30 µg	€307
Anti-Dopamine D ₁ Receptor Rabbit pAb	ELISA, FFS, FS, IB; NOT PS	human, mouse, rat	324390	100 µl	€369

10T: 10 tests by immunoblotting; **AB:** adhesion blockade; **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW


Anti-CRM1 Rabbit pAb Cat. No. ST1100

Lane 1: Detection of human CRM1 by immunoblotting. Sample: Whole cell lysate (50 µg) from HeLa cells. Primary antibody: Anti-CRM1 Rabbit pAb (Cat. No. ST1100) (1 µg/ml). Detection: Chemiluminescence.

Lane 2: Detection of human CRM1 by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from HeLa cells. Antibody for IP: Anti-CRM1 Rabbit pAb (ST1100) (5 µg/500 µg total protein in 500 µl). Immunoblotting conditions: same as Lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Dopamine D ₂ Receptor Rabbit pAb	ELISA, IB, IH	human, rat	324393	100 µl	€369
Anti-Dopamine D ₂₅ Receptor Rabbit pAb	ELISA, IB, IH	human, rat	324396	100 µl	€369
Anti-Dopamine D ₃ Receptor Rabbit pAb	ELISA, FC, IB, IH	human, rat	324402	100 µl	€369
Anti-Dopamine D ₄ Receptor Rabbit pAb	ELISA, IB, IH	human, rat	324405	100 µl	€369
Anti-Dopamine D ₅ Receptor Rabbit pAb	ELISA, FC, IB, IH	human, rat	324408	100 µl	€369
Downstream Effector Caspase Sampler Kit	IB		ASK23	1 each	€406
Anti-DPC4/Smad4 Mouse mAb (DCS.46)	IB, IP	human, mouse	OP111	100 µg	€338
Anti-DR4/TRAIL Receptor-1, N-Terminal (1-20) Rabbit pAb	IB	human	318212	100 µg	€319
PhosphoDetect™ Anti-Dynamin (pSer ⁷⁷⁴) Sheep pAb	DB, IB	human, mouse, rat	PS1000	100 µl	€385
Anti-Dynamin II Rabbit pAb	IB	human, mouse, rat	PC407	100 µl	€346
Anti-Dynamin Mouse mAb (D5)	FS, IB, IF, IP	human, mouse, rat	NB28	100 µg	€338
Anti-Dynorphin A (1-17) Rabbit pAb	FS, RIA	porcine, rat	PC210	100 µl	€291
Anti-DYRK Sheep pAb	IB, IF, IP	human, mouse, rat	324446	100 µg	€319
Anti-E2F1 (Ab-1) Mouse mAb (KH95)	IB; NOT IF, NOT IP	human	NA49	100 µg	€338
Anti-E2F4 (Ab-2) Mouse mAb (WUF11)	IB	human	NA78	100 µl	€409
Anti-E2F6 Mouse mAb (TFE61)	IB	human	NA79	100 µl	€461
Anti-E47 (517-531) Rabbit pAb	ELISA, IB	human	PC695	100 µl	€409
Anti-E6-AP Rabbit pAb	IB, IP	human	PC543	50 µl	€253
Anti-Early Endosomal Antigen 1 (1391-1410) Rabbit pAb	IB, IF	canine, hamster, human, mouse, rat	324610	100 µg	€362
Anti-EBNA Mouse mAb (E8.26)	IB, IF, IP	EBV-infected mammalian B lymphocytes	DP15L	100 µg	€338
Anti-4EBP1/PHAS-I (101-118) Rabbit pAb	GS, IB	human, mouse, rat	516676	100 µg	€319
Anti-EBP50 Rabbit pAb	IB, IP, PS	human	324620	100 µg	€384
Anti-EDG-1, C-Terminal Rabbit pAb	IB	human	PC620	50 µg	€420
Anti-EDG-2 (328-344) Rabbit pAb	IB	human, rat	324635	50 µg	€243
Anti-EDG-3 (Ab-1) Mouse mAb	IB	human, rat	GR41	50 µg	€338
Anti-EDG-3 (Ab-2) Mouse mAb (AS60)	IB	human, rat	GR42	50 µg	€338
Anti-EDG-5 (Ab-1) Mouse mAb	IB	human, rat	GR44	50 µg	€338
Anti-EDG-7, C-Terminal Rabbit pAb	IB	human	PC621	50 µg	€420
Anti-EDG-8, C-Terminal Rabbit pAb	IB	human	PC622	50 µg	€420
Anti-Elastin Rabbit pAb	ELISA, IB, IH, IP	bovine, canine, human, rabbit; not rat	324756	500 µl	€445
PhosphoDetect™ Anti-eIF-2α (Ab-1) (pSer ⁵²) Rabbit pAb	IB	human, mouse, yeast	PC638	10 T	€458

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-ESPL1/Separase Rabbit pAb, Cat No. ST1118

Lane 1: Detection of human separase by immunoblotting Samples: Whole cell lysate (20 µg) from HepG2 cells. Primary antibody: Anti-Separase Rabbit pAb (Cat. No. ST1118) (1 µg/ml) Detection: Chemiluminescence


Lane 2: Detection of human separase by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate (500 µg) from HepG2 cells. Antibody for immunoprecipitation: Anti-Separase Rabbit pAb (Cat. No. ST1118) (5 µg/500 µg total protein). Immunoblotting conditions: Sample half of the immunoprecipitated separase. Primary antibody: Anti-Separase Rabbit pAb (Cat. No. ST1118) (2 µg/ml). Detection: Chemiluminescence

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-eIF-4E (Ab-1) (pSer ²⁰⁹) Rabbit pAb	IB	human, mouse, rabbit, rat	PC639	10 T	€458
PhosphoDetect™ Anti-Elk1 (pSer ³⁸³) Rabbit pAb	DB, IB	human	324643	100 µl	€455
Anti-Elk1 Rabbit pAb	IB, PS	human, mouse, rat	324641	100 µl	€355
Anti-Endonuclease G (56-70) Rabbit pAb	IB	human, mouse, rat	PC684	100 µg	€338
Anti-Endoplasmic Reticulum Protein 72 (623-638) Rabbit pAb	IB, IC, IP	bovine, canine, guinea pig, hamster, human, monkey, mouse, rabbit, rat	324734	25 µg	€290
Anti-Endostatin (Ab-2) Rabbit pAb	IB; NOT PS	human	PC541	100 µg	€294
Anti-Endostatin (Ab-4) (1-14) Rabbit pAb	IB; NOT PS	human	PC599	100 µl	€294
Anti-Endostatin (Ab-5) Mouse mAb (NE1.479)	IB; NOT PS	human	GF57	100 µg	€338
Anti-Endothelin (Ab-1) Mouse mAb (3B8C9)	ELISA	a broad range of species	CP43	100 µl	€315
Anti-Endothelin A Receptor (413-426) Rabbit pAb	IB	rat	324758	1 Set	€231
Anti-Endothelin B Receptor, C-Terminal Sheep pAb	IB, IF, IP	human, mouse, rat, sheep	324755	250 µg	€390
Anti-Endothelin-1 (Ab-1) Mouse mAb (TR.ET.48.5)	IC, PS, RIA	a broad range of species	CP44	100 µl	€315
Anti-EpCAM (Ab-1) Mouse mAb (1Dg.F3)	PS; NOT IB, NOT IF, NOT IP	human; not mouse, not rat	CP63	100 µg	€338
Anti-EpCAM-1 (Ab-1) Mouse mAb (VU-1D9)	FC, FS, IB, IP, PS	human	OP187	100 µg	€338
Anti-EphA2/Eck Mouse mAb (B2D6)	PS	human	OP147	100 µg	€338
Anti-EphB1, Cytoplasmic Tyrosine Kinase Domain (586-984) Sheep pAb	IP	human	324828	250 µg	€319
Anti-Epidermal Growth Factor (Ab-1) Mouse mAb (144-8)	FS, IP; NOT NS	human	GF01	200 µg	€338
Anti-Epidermal Growth Factor (Ab-3) Rabbit pAb	IB, PS	human	PC08	100 µg	€315
Anti-Epidermal Growth Factor Receptor (Ab-1) Mouse mAb (528)	FC, FS, IF, IP, NT, PS	human; not mouse, not rat	GR01	100 µg	€294
Anti-Epidermal Growth Factor Receptor (Ab-1) Mouse mAb (528)	FC, FS, IF, IP, NT, PS	human; not mouse, not rat	GR01L	100 µg	€294
Anti-Epidermal Growth Factor Receptor (Ab-2) Mouse mAb (455)	IB; NOT NS	human; not mouse, not rat	GR02L	100 µg	€338
Anti-Epidermal Growth Factor Receptor (Ab-3) Mouse mAb (225)	IP	human; not mouse, not rat	GR13	100 µg	€338
Anti-Epidermal Growth Factor Receptor (Ab-3) Mouse mAb (225)	IP, NT	human; not mouse, not rat	GR13L	100 µg	€338
Anti-Epidermal Growth Factor Receptor (Ab-5) Mouse mAb (EGFR.1)	FC, IF, IP, NT	human; not mouse, not rat	GR15	100 µg	€338
Anti-Epidermal Growth Factor Receptor (Ab-5) Mouse mAb (EGFR.1)	FC, IF, IP, NT	human; not mouse, not rat	GR15L	100 µg	€338
Anti-Epidermal Growth Factor Receptor (Ab-6) (985-996) Rabbit pAb	IB	human, mouse	PC98	100 µg	€319
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pSer ¹⁰⁴⁷) Mouse mAb (1H9)	ELISA, IB	human	324863	1 Set	€350
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pTyr ¹⁰⁶⁸) Rabbit pAb	IB	chicken, human, mouse, rat	324867	10 T	€373
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pTyr ¹¹⁴⁹) Rabbit pAb	IB	chicken, human, mouse, rat	324869	10 T	€373
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pTyr ¹¹⁷³) Mouse mAb (9H2)	ELISA, IB, IP	human, mouse, rat	324864	1 Set	€374
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pTyr ⁸⁴⁵) Rabbit pAb	IB	human	324865	10 T	€373
PhosphoDetect™ Anti-Epidermal Growth Factor Receptor (pTyr ⁹⁹²) Rabbit pAb	IB	human	ST1007	10 T	€350
Anti-α-Epithelial Cell Sodium Channel (20-42) Rabbit pAb	IB, IP	human, rat	324870	100 µg	€362
Anti-α-Epithelial Cell Sodium Channel (627-643) Rabbit pAb	IB, IF, IP	human	324906	50 µl	€231
Anti-β-Epithelial Cell Sodium Channel (619-640) Rabbit pAb	IP; NOT IB	human, rat	324874	100 µg	€368
Anti-EPLIN C-Terminal Rabbit pAb	IB	human	CA1014	50 µg	€189

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-EPLIN, Internal Rabbit pAb	IP; NOT IB	human	CA1019	50 µg	€189
Anti-ERAB/ABAD/HSD-10 (101-119) Rabbit pAb	IB	human	324918	100 µl	€301
Anti-c-ErbB2/c-Neu (Ab-1) (866-880) Rabbit pAb	IB, IF, IP, PS	human, rat	PC04	100 µg	€319
Anti-c-ErbB2/c-Neu (Ab-2) Mouse mAb (9G6)	FS, IP, PS	human	OP14	100 µg	€294
Anti-c-ErbB2/c-Neu (Ab-2) Mouse mAb (9G6)	FS, IP, PS	human	OP14L	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-3) Mouse mAb (3B5)	FS, IB, IF, IP, PS	human, mouse	OP15	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-3) Mouse mAb (3B5)	FS, IB, IF, IP, PS	human, mouse	OP15L	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-3) Mouse mAb (3B5)	FS, IB, IF, IP, PS	human, mouse	OP15T	10 µg	€77
Anti-c-ErbB2/c-Neu (Ab-3) Mouse mAb (3B5) Fluorescein Conjugate	FS, IF, PS	human, mouse	OP15F	100 µg	€346
Anti-c-ErbB2/c-Neu (Ab-4) Mouse mAb (7.16.4)	FS, GIA, IB, IC, IP; NOT NS	human, rat	OP16	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-4) Mouse mAb (7.16.4)	FS, GIA, IB, IC, IP; NOT NS	human, rat	OP16L	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-5) Mouse mAb (TA-1)	FC, FS, IB, IF, IP; NOT NS	human; not mouse, not rat	OP39	100 µg	€338
Anti-c-ErbB2/c-Neu (Ab-6) Mouse mAb (OD-3)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human; not mouse, not rat	OP67	100 µg	€338
PhosphoDetect™ Anti-c-ErbB2/c-Neu (pSer ¹¹¹³) (Ab-1) Mouse mAb (a-PS1113)	ELISA, IB	human	OP195	100 µl	€458
Anti-c-ErbB3 (Ab-1) (860-873) Rabbit pAb	IB	human	PC27L	100 µg	€319
Anti-c-ErbB3 (Ab-4) Mouse mAb (H3.90.6 or Ab90)	FC, IF, IP; NOT IB, NOT NS	human	OP119	100 µg	€338
Anti-c-ErbB3 (Ab-5) Mouse mAb (H3.105.5)	FC, IF, IP; NOT FS, NOT IB, NOT PS	human, mouse	OP120	100 µg	€338
Anti-c-ErbB4 (Ab-1) Mouse mAb (H4.77.16)	FC, IP; NOT IB	human, mouse, rat	OP121	100 µg	€338
Anti-c-ErbB4 (Ab-2) Mouse mAb (HFR-1)	FS, IP, PS; NOT IB	human	OP122	100 µg	€338
PhosphoDetect™ Anti-ERK5 (pThr ²¹⁹ /Tyr ²²⁰) Rabbit pAb	IB	human	442688	10 T	€389
Anti-ERK5 Rabbit pAb	IB, IP	human, monkey, mouse, rat	ST1080	50 µl	€189
Anti-Erythropoietin (1-12) Rabbit pAb	IB	human	PC18L	100 µg	€319
NEW Anti-ESPL1/Separase Rabbit pAb	IB, IP	human	ST1118	50 µg	€135
Anti-Estrogen Receptor (287-300) Mouse mAb (H-151)	GS, IB, IH, IP	bovine, human, mouse, rat	330659	200 µl	€397
Anti-Estrogen Receptor (Ab-1) Mouse mAb (TE-111)	IB, IP, PS	human, mouse; not rat	GR17	100 µg	€291
Anti-Estrogen Receptor, C-Terminal (582-595) Mouse mAb (C-542)	GS, IB, IC, IH, IP	human, mouse, rat	330658	200 µl	€384
Anti-Estrogen Receptor β (Ab-1) (485-503) Rabbit pAb	FS, GS, IB, PS	mouse, rat	PC168	50 µg	€319
Anti-Estrogen Receptor β (Ab-2) Mouse mAb (9.88)	IB, IP; NOT PS	human, mouse, rat	GR39	100 µl	€270
Anti-Estrogen Receptor β (Ab-3) Mouse mAb (10.76)	DB, IP; NOT IB	human	GR40	100 µg	€338
Anti-ETO (Ab-1) Rabbit pAb	IB, IP; NOT FS, NOT IF, NOT PS	human, mouse	PC283	100 µg	€319
Anti-c-Ets-1 (201-214) Rabbit pAb	ELISA, IB	human	PC694	100 µl	€409
Anti-Excitatory Amino Acid Transporter (GLT-1) Rabbit pAb	IB, IP	rat	PC154	100 µl	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-F₁F₀-α Mouse mAb (7H10BD4F9) Cat. No. AP1036

Detection of human F₁F₀-α by immunocytochemistry. Sample: 143B osteosarcoma cells. Primary antibody: Anti-F₁F₀-α Mouse mAb (7H10BD4F9) (Cat. No. AP1036) (5 µg/µl). Detection: Fluorescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Excitatory Amino Acid Transporter 2 (EEAT2) (Ab-2) Rabbit pAb	IB; NOT FS	mouse, rat	PC339	100 µg	€338
Anti-F1F0-a Mouse mAb (7H10BD4F9)	CIA, IB, IC	bovine, human, mouse, rat	AP1036	50 µg	€294
Anti-Factor B Goat pAb	IEP, OT, RID	human	341272	2 ml	€153
Anti-Factor H Goat pAb	OT, RID	human	341276	1 ml	€92
Anti-Factor VIII Related Antigen/vWF (Ab-1) Rabbit pAb	FS, IF, PS	human	PC313	500 µl	€338
Anti-Factor VIII Related Antigen/vWF (Ab-2) Mouse mAb (GMA-022)	ELISA, IB	human	GF54	100 µg	€202
Anti-FADD (125-140) Rabbit pAb	IB	human, monkey, mouse (weakly), rat (weakly)	341282	100 µg	€384
Anti-FAK (853-1052) Rabbit pAb	IB, IF, IP	amphibian, canine, human, mouse, rat	PC314	500 µg	€319
PhosphoDetect™ Anti-FAK (pSer ⁷²²) Rabbit pAb	IB	chicken, human, mouse, rat	ST1008	10 T	€373
PhosphoDetect™ Anti-FAK (pTyr ³⁹⁷) Rabbit pAb	IB, IH	chicken, <i>Drosophila</i> , frog, human, mouse, rat	341292	10 T	€373
PhosphoDetect™ Anti-FAK (pTyr ⁴⁰⁷) Rabbit pAb	IB, IC	chicken, frog, human, mouse	ST1010	10 T	€373
PhosphoDetect™ Anti-FAK (pTyr ⁵⁷⁶) Rabbit pAb	IB	human	341294	10 T	€373
PhosphoDetect™ Anti-FAK (pTyr ⁵⁷⁷) Rabbit pAb	IB, IC	chicken, human, mouse	ST1011	10 T	€373
PhosphoDetect™ Anti-FAK (pTyr ⁸⁶¹) Rabbit pAb	IB, IH	chicken, human, mouse	PS1008	100 µl	€373
Anti-FAP Mouse mAb (F11-24)	FC, FS, IP	human	OP188	100 µg	€338
Anti-Farnesyl Rabbit pAb	ELISA, IF	a broad range of species	341286	100 µg	€384
Anti-Farnesyltransferase Rabbit pAb	IB	mammalian	341285	50 µl	€390
Anti-Fas (29-44) Rabbit pAb	IB	mouse, rat	341289	100 µg	€384
Anti-Fas (Ab-1) (321-335) Rabbit pAb	FS, IB, PS	human, mouse, opossum, rat	PC69	100 µg	€319
Anti-Fas (Ab-2) Mouse mAb (DX2)	APOP, IB, IF; NOT FS, NOT IP, NOT PS	human; not mouse, not rat	AM01	100 µg	€278
Anti-Fas (Ab-2) Mouse mAb (DX2)	APOP, IB, IF; NOT FS, NOT IP, NOT PS	human; not mouse, not rat	AM01L	100 µg	€338
Anti-Fas Ligand (Ab-1) (261-277) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human, mouse	PC78	100 µg	€319
Anti-Fas Ligand (Ab-3) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human; not mouse, not rat	PC105	100 µg	€319
Anti-Fas Ligand (Ab-5) Mouse mAb (8B8)	IF; NOT FS, NOT IB, NOT IP, NOT PS	human; not mouse, not rat	AM29	100 µg	€338
Anti-v-Fes Rat mAb (F-113)	IP, NT	human, mouse, rat	OP02	100 µg	€338
Anti-Fibrinogen Rabbit pAb	ELISA, IEP, RID	human	341552	1 ml	€100
Anti-Fibroblast Antigen Mouse mAb (AS02)	FC, FS, IB, IF, IP	human, monkey; not mouse, not rat	CP28	200 µg	€378
Anti-Fibroblast Growth Factor Receptor (Ab-1) Mouse mAb (VBS-7)	FS, IB, IP, PS	bovine, chicken, guinea pig, human, mouse, rat	GR21	100 µg	€346
Anti-Fibroblast Growth Factor Receptor (Ab-2) (351-366) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC99	100 µl	€319
Anti-Fibroblast Growth Factor, Basic (Ab-1) (147-153) Rabbit pAb	ELISA, PS	human	PC15	100 µg	€319
Anti-Fibroblast Growth Factor, Basic (Ab-2) (40-63) Rabbit pAb	IB, PS	human, rat	PC16	20 µg 100 µg	€86 €319

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections



NEW

Anti-GFAP Cocktail Mouse mAb (SMI-22) Cat. No. NE1015


Detection of rat glial fibrillary acidic protein by staining frozen sections. Sample: Rat brain.
Primary antibody: Anti-GFAP Cocktail Mouse mAb (SMI-22) (Cat. No. NE1015) (1:1000).
Detection: Fluorescence (red) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Fibroblast Growth Factor, Basic (Ab-3) Mouse mAb (3H3)	ELISA, IB, IP, NT, PS	human; not mouse, not rat	GF22L	200 µg	€338
Anti-Fibronectin (Ab-1) Mouse mAb (N-294)	CIA, FS, IB, IF, PS	human	CP13L	100 µg	€338
Anti-Fibronectin (Ab-3) Mouse mAb (FBN11)	FS, IB, IF, IP, PS	human, mouse, rat	CP70	100 µg	€338
Anti-FKBP59/HSP56 (441-458) Rabbit pAb	IB	human, mouse, rabbit	343250	100 µg	€381
Anti-Fli-1 (1-239) Sheep pAb	IB	human, mouse	PC425	250 µg	€338
Anti-FLIP (447-464) Rabbit pAb	IB	human	343002	100 µg	€319
Anti-FLIPγ/δ (191-209) Rabbit pAb	IB	human, mouse	343006	100 µg	€319
Anti-Cleaved α-Fodrin (Asp ¹¹⁸⁹) Rabbit pAb	IB, IC	human, mouse	AP1011	50 µl	€278
Anti-Fodrin Mouse mAb (Fod009)	ELISA, FC, IB, IH, IP	bovine, canine, human, mouse, porcine, rabbit, rat	344050	50 µg	€360
Anti-c-Fos (Ab-1) Mouse mAb (2G9C3)	IB, IC; NOT PS	human, mouse, rat	OP17	200 µg	€338
Anti-c-Fos (Ab-1) Mouse mAb (2G9C3)	IB, IC; NOT PS	human, mouse, rat	OP17L	200 µg	€338
Anti-c-Fos (Ab-2) (4-17) Rabbit pAb	FS, IB, IF, IP, PS; NOT FFS	human, mouse, rat	PC05	100 µg	€319
Anti-c-Fos (Ab-2) (4-17) Rabbit pAb	FS, GS, IB, IF, IP, PS; NOT FFS	human, mouse, rat	PC05L	100 µg	€319
Anti-c-Fos (Ab-2) Rabbit pAb	FS, IB, IF, IP, PS; NOT FFS	human, mouse, rat	PC05T	10 µg	€77
Anti-c-Fos (Ab-3) Mouse mAb (CS93.1.18)	IB, PS	human, mouse	OP53	100 µg	€338
Anti-c-Fos (Ab-5) (4-17) Rabbit pAb	FFS	human, mouse, rat	PC38	100 µl	€248
Anti-c-Fos (Ab-5) (4-17) Rabbit pAb	FFS	human, mouse, rat	PC38T	20 µl	€77
PhosphoDetect™ Anti-Fos (pSer ³⁷⁴) Mouse mAb (34E4)	ELISA, IB	human	ST1029	1 Set	€373
c-Fos Peptide-2	FFS, FS, IB, IP, PS		PP10	100 µg	€84
PhosphoDetect™ Anti-FOXO1/FKHR (pThr ²⁴) Sheep pAb	IB	human, mouse, rat	342575	50 µg	€319
Anti-Fra2 (289-305) Rabbit pAb	ELISA, IB	human, mouse	PC696	100 µl	€409
Anti-Gα-Subunit, Internal (40-54) Rabbit pAb	IB	mammalian	371737	50 µl	€409
Anti-Gβγ-Complex Rabbit pAb	ELISA, IB	bovine, human, rat	371821	50 µl	€390
Anti-Gβ-Subunit, Internal (127-139) Rabbit pAb	IB	human, mouse	371738	50 µl	€390
Anti-Gβ-1-Subunit, Internal (26-39) Rabbit pAb	ELISA, IB, IH	mammalian	371813	50 µl	€390
Anti-Gβ-3-Subunit, Internal (26-39) Rabbit pAb	ELISA, IB, IH	human, mouse, rat	371817	50 µl	€390
Anti-Gγ-2-Subunit Rabbit pAb	ELISA, IB	human	371847	50 µl	€390
Anti-G ₁₂ α-Subunit, C-Terminal (327-342) Rabbit pAb	ELISA, IB, IP	mammalian	371778	50 µl	€390
Anti-G ₁₆ α-Subunit, C-Terminal (360-373) Rabbit pAb	ELISA, IB	human, mouse, rat	371783	50 µl	€384
Anti-G _q α-1 and G _q α-2-Subunits (Ab-4) Rabbit pAb	IB	mammalian	PC61	100 µl	€162
Anti-G _q α-1 and G _q α-2-Subunits, C-Terminal (345-354 and 346-355) Rabbit pAb	IB	mammalian	371723	50 µl	€384
Anti-G _q α-1-Subunit, Internal (159-168) Rabbit pAb	IB	mouse	371720	50 µl	€384
Anti-G _q α-2-Subunit, Internal (85-100) Rabbit pAb	ELISA, IB	human, mouse, rat	371727	50 µl	€390
Anti-G _q α-3- and G _o α-Subunits, C-Terminal (345-354) Rabbit pAb	IB	mammalian	371726	50 µl	€397
Anti-G _q α-3-Subunit, C-Terminal (345-354) Rabbit pAb	IB, IP	mammalian	371729	50 µl	€384
Anti-G _o α-Subunit Rabbit pAb	ELISA, IB	mouse, rat	371728	50 µl	€390
Anti-G ₁₂ α- and G ₁₁ α-Subunits, Internal (283-300) Rabbit pAb	IB, NT	human, mouse, rat	371754	50 µl	€384
Anti-G _q α- and G ₁₁ α-Subunits, C-Terminal Rabbit pAb	ELISA, IB, IP	mammalian	371751	50 µl	€409
Anti-G _q α-Subunit Rabbit pAb	ELISA, IB	mammalian	371730	50 µl	€390
Anti-G _q α-Subunit, Internal (115-133) Rabbit pAb	IB		371752	50 µl	€413
Anti-G _q α-Subunit, C-Terminal (385-394) Rabbit pAb	IB	mammalian	371732	50 µl	€413
Anti-G _s α-Subunit Rabbit pAb	ELISA, IB		371731	50 µl	€390
Anti-G _i α-1-Subunit Mouse mAb (TF15)	ELISA, IB	mammalian	371740	100 µl	€390
Anti-G _i α-Subunit, N-Terminal (3-18) Rabbit pAb	IB	mammalian	371741	50 µl	€390

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchtterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
G-Protein Subunit Antibody Set	IB		371770	1 Set	€501
Anti-GABA Rabbit pAb	FS	rat	PC213L	100 µl	€397
Anti-GABA Transporter GAT-1 (588-299) Rabbit pAb	FS, VS	rat	PC250L	100 µl	€390
Anti-GABA _A Receptor γ2-Subunit Rabbit pAb	IB	human	345231	50 µl	€390
PhosphoDetect™ Anti-GABA _B Receptor (pSer ⁸⁹²) Rabbit pAb	IB, IF, IH	rat	345238	100 µl	€455
Anti-GADD34 (664-674) Goat pAb	IB	human	PC720	100 µg	€397
Anti-GAL4, DNA Activation Domain Mouse mAb (5C8)	IB	yeast	345765	100 µg	€362
Anti-GAL4, DNA Binding Domain Mouse mAb (8C-1)	IB	yeast	345766	100 µg	€362
Anti-β-Galactosidase Mouse mAb (200-193)	IB, IP	<i>E. coli</i>	OB02	100 µg	€231
Anti-Ganglioside GD ₃ Mouse mAb (R24)	ELISA, IB	a broad range of species	345756	50 µl	€196
Anti-Ganglioside GM ₁ Rabbit pAb	ELISA, IB	a broad range of species	345757	100 µl	€242
Anti-Ganglioside GM ₁ , Asialo Rabbit pAb	ELISA, IB	a broad range of species	345758	100 µl	€242
Anti-Ganglioside GM ₂ Rabbit pAb	ELISA, IB	a broad range of species	345759	100 µl	€242
Anti-GAP-43 Mouse mAb (5E7)	FS, IB, IF, PS	human, rat	CP09	100 µg	€338
Anti-GAPDH Mouse mAb (6C5)	ELISA, IB, IC	canine, chicken, fish, frog, human, mouse, porcine, rabbit, rat; not bovine, not yeast	CB1001	500 µg	€347
Anti-Giantin Mouse mAb	IB, IC, IH	rat	324450	100 µl	€260
Anti-Glial Fibrillary Acidic Protein (GFAP) Mouse mAb (G-A-5)	FS, IB, IF, IP, PS	chicken, human, porcine, rat	IF03L	100 µg	€319
NEW Anti-Glial Fibrillary Acidic Protein (GFAP) Cocktail Mouse mAb (2C13B9)	ELISA, FS, IB, IC, PS	bovine, canine, chicken, guinea pig, human, mouse, porcine, rat, sheep	NE1015	100 µl	€283
Anti-Glial Fibrillary Acidic Protein Rat mAb (2.2B10)	ELISA, IB, IH, IP	bovine, human, rodent	345860	100 µg	€334
Anti-Glucocorticoid Receptor β (Ab-1) Rabbit pAb	IB, IC, IP	human	PC171	50 µl	€304
Anti-Glucocorticoid Receptor (Ab-1) (150-175) Rabbit pAb	GS, IB, IC, IP	human	PC170	200 µl	€304
Anti-Glucocorticoid Receptor (Ab-2) Mouse mAb (BuGR2)	FS, IB, IC, IP, PS	guinea pig, mouse, rabbit, rat, sheep; not amphibian, not avian, not primate	GR32L	100 µg	€319
Anti-Glucose Transporter-1 Rabbit pAb	ELISA, FS, IB, IF, IP	human, mouse, rabbit, rat	400060	50 µg	€319
Anti-Glucose Transporter-2 Rabbit pAb	ELISA, IB	human, mouse, rat	400061	50 µg	€319
Anti-Glucose Transporter-3 Rabbit pAb	ELISA, FS, IB, IF, IP	human; not mouse, not rat	400062	50 µg	€319
Anti-Glucose Transporter-4 Rabbit pAb	ELISA, FS, IB, IP	human, mouse, rat	400064	50 µg	€319
Anti-Glucose Transporter-5 Rabbit pAb	ELISA, FS, IB, IF	rat; not human	400066	50 µg	€319
PhosphoDetect™ Anti-GluR1 (pSer ⁸³¹) Rabbit pAb	DB, IB, IH	rat	400068	150 µl	€455
PhosphoDetect™ Anti-GluR1 (pSer ⁸⁴⁵) Rabbit pAb	IB	mouse, rat	400059	50 µl	€233
PhosphoDetect™ Anti-GluR1 (pSer ⁸⁴⁵) Rabbit pAb	DB, IB	rat	400069	100 µl	€455

10T: 10 tests by immunoblotting; **AB:** adhesion blockade; **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** outherlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW


Anti-IDE/Insulysin Rabbit pAb, Cat No. ST1120

Lane 1: Detection of human IDE/insulysin by immunoblotting. Sample: Whole cell lysate (20 µg) from HT1080 cells. Primary antibody: Anti-IDE/Insulysin Rabbit pAb (Cat. No. ST1120) (0.5 µg/ml). Detection: Chemiluminescence.

Lane 2: Detection of human IDE/insulysin by immunoprecipitation follow by immunoblotting. Sample: Whole cell lysate (500 µg) from HT1080 cells (lane 2) and negative control sample (lane 3). Antibody for immunoprecipitation: Anti-IDE/Insulysin Rabbit pAb (Cat. No. ST1120) (5 µg). Immunoblotting conditions: same as lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-GluR1 Rabbit pAb	IB, IF	rat	PC246	100 µg	€350
Anti-Glutathione Peroxidase Mouse mAb (347)	IB, IC, PS	human	ST1000	100 µg	€397
Anti-Glutathione Rabbit pAb	ELISA	a broad range of species	354106	50 µl	€294
Anti-Glutathione S-Transferase (Ab-1) Rabbit pAb	AC, ELISA, IB, IP	<i>S. japonicum</i> ; not <i>E. coli</i> , not human	PC53	100 µg	€231
Anti-Glutathione S-Transferase (Ab-2) Mouse mAb (8-326)	IB	<i>S. japonicum</i>	OB03	20 µg 100 µg	€86 €338
Anti-Glutathione S-Transferase (Ab-3) Mouse mAb (TC112)	IB	<i>S. japonicum</i> glutathione-S-transferase	OB07	100 µg	€338
Anti-Glutathione S-Transferase A1-1 Rabbit pAb	IB	human, rat	354206	100 µl	€390
Anti-Glutathione S-Transferase M1-1 Rabbit pAb	IB	human, rat	354209	100 µl	€390
Anti-Glutathione S-Transferase P1-1 Rabbit pAb	IB	human, rat	354212	100 µl	€390
Anti-Glycine Receptor Rabbit pAb	IB, IP	human, mouse	357001	200 µg	€416
PhosphoDetect™ Anti-Glycogen Synthase Kinase 3β (pSer ⁹) Mouse mAb (2D3)	ELISA, IB	human, mouse	361527	1 Set	€350
Anti-Glycogen Synthase Kinase 3β, C-Terminal (334-348) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	361528	100 µg	€404
PhosphoDetect™ Anti-Glycogen Synthase Kinase-3α/β (pTyr ^{279/216}) Rabbit pAb	DB, ELISA, IB, IC	human, mouse, rat	ST1013	10 T	€350
Anti-Glycogen Synthase Kinase-3α/β Mouse mAb (1H8)	ELISA, IB	bovine, canine, hamster, human, mouse, ovine, porcine, rabbit, rat	368662	100 µg	€384
PhosphoDetect™ Anti-Glycogen Synthase Kinase-3β (pSer ⁹) (Ab-1) Rabbit pAb	DB, ELISA, IB	human, rat, <i>Xenopus</i>	PC242	50 µg	€490
PhosphoDetect™ Anti-Glycogen Synthase Kinase-3β, (pSer ⁹) Rabbit pAb	IB	human, rat	PS1018	10 T	€350
Anti-cGMP Rabbit pAb	RIA	a broad range of species	370659	1 Set	€385
Anti-Golgi Apparatus Mouse mAb (AE-6)	IB, IF; NOT IP, NOT PS	human; not mouse, not rat	NB37	10 µg 100 µg	€86 €338
Anti-trans-Golgi Network 38 Mouse mAb (2F7.1)	FC, IB, IF	rat	NB10	100 µl	€338
Anti-Granulocyte Colony Stimulating Factor (Ab-1) Mouse mAb (5.24)	IB, IF, PS	human	GF05	200 µg	€338
Anti-Granulocyte-Macrophage Colony Stimulating Factor (Ab-1) Mouse mAb (GM4.1.9)	IP, NT	human	GF13L	100 µg	€338
Anti-Granzyme B Mouse mAb (2C5/F5)	IB, IP; NOT FS, NOT PS	human, rat; not mouse	AM52	50 µg	€294
Anti-Grb7 (522-532) Goat pAb	IB	human	PC721	100 µg	€397
Anti-Grb10 Rabbit pAb	IB, IH, IP	human, mouse	PC687	100 µl	€319
Anti-Green Fluorescent Protein (26-39) Rabbit pAb	IB, IC		PC408	100 µl	€346
Anti-GRIM-19 Mouse mAb (6E1BH7)	IB, IC	bovine, human, mouse, rat	AP1033	50 µg	€294
Anti-GroEL Mouse mAb (9A1/2)	IB, IP	<i>E. coli</i>	CB1007	50 µl	€291
Anti-Grp78 (645-354) Rabbit pAb	IB, IP	hamster, mouse, rat	368673	50 µl	€307
Anti-Guanylate Binding Protein-1 Rat mAb (1B1)	IB, IF, PS	human	IM1011	50 µg	€220

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections




NEW

Anti-GRIM-19 Mouse mAb (6E1BH7) Cat. No. AP1033

Detection of human GRIM-19 by immunocytochemistry. Sample: MRC5 fibroblast cells. Primary antibody: Anti-GRIM-19 Mouse mAb (6E1BH7) (Cat. No. AP1033) (1 µg/ml). Detection: Fluorescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Guanylate Cyclase α, β, γ , Soluble Chicken pAb	IB	human, rat	371716	100 μ l	€384
Anti-H ⁺ ,K ⁺ -ATPase, α -Subunit (1019-1034) Rabbit pAb	ELISA, IB, IC, PS	canine, human, mouse, rat	119101	50 μ l	€371
Anti-H ⁺ ,K ⁺ -ATPase, α -Subunit (1-17) Rabbit pAb	ELISA, IB, IC, PS	canine, human, mouse	119102	50 μ l	€371
Anti-H ⁺ ,K ⁺ -ATPase, α -Subunit Mouse mAb (HK 12.18)	ELISA, FS, IB, IP, PS	canine, gerbil, guinea pig, human, porcine, rabbit, rat	119100	50 μ l	€371
PhosphoDetect™ Anti-H2AX (pSer ¹³⁹) Rabbit pAb	IB, IC	human, mouse	DR1017	100 μ g	€350
Anti-H2AX Rabbit pAb	IB, IP	human, mouse	DR1016	100 μ g	€360
Anti-HAUSP Rabbit pAb	IB	human	DR1014	50 μ g	€189
Anti-Heme Oxygenase-1 (1-261) Rabbit pAb	FC, IB, IC, IH, IP	canine, guinea pig, hamster, human, monkey, mouse, rabbit, rat, sheep	374089	200 μ l	€390
Anti-Heme Oxygenase-1 (1-30) Rabbit pAb	IB, IH, IP	canine, hamster, human, monkey, mouse, rat	374090	100 μ l	€390
Anti-Heme Oxygenase-1 (Ab-1) (2-23) Rabbit pAb	FFS, FS, IB	mouse, rat	PC340	100 μ g	€338
Anti-Heme Oxygenase-1 Mouse mAb (HO-1-1)	IB, IC, IP	bovine, canine, human, monkey, mouse, rat	374087	100 μ g	€390
Anti-Heme Oxygenase-2 (21-50) Rabbit pAb	ELISA, IB, IC, IP	canine, human, mouse, rat	374094	100 μ g	€390
Anti-Hemocyanin, Keyhole Limpet Rabbit pAb	ELISA, IB		374806	500 μ l	€131
Anti-Hemoglobin F Sheep pAb	IEP	human	374824	1 ml	€206
Anti-Heparin Binding Epidermal Growth Factor (Ab-1) Goat pAb	ELISA, IB, NT	human	PC319L	100 μ g	€338
Anti-Heregulin (Ab-1) Goat pAb	ELISA, IB, NT, PS	human	PC185L	100 μ g	€344
Anti-HERG (1106-1159) Rabbit pAb	IB	horse, human, mouse	373500	1 Set	€278
Anti-HGF (Scatter Factor) (Ab-1) Mouse mAb (4C12.1)	IB, IP	human	GF39	100 μ g	€338
Anti-Hip-2 (1-12) Rabbit pAb	FS, IB, IC, PS	human, mouse, rat	NE1011	50 μ l	€150
Anti-Histidine-Tagged Protein Mouse mAb (13/45/31/2)	ELISA, IB, IF, IP	all His•Tag® proteins	OB05	100 μ g	€338
Anti-Histone Deacetylase 1 (467-482) Rabbit pAb	IB, IF	canine, hamster, human, mouse, rat	PC544	50 μ g	€253
Anti-Histone Deacetylase 2 (475-488) Rabbit pAb	IB, IC	canine, hamster, human, mouse, rat	382153	100 μ g	€346
Anti-Histone Deacetylase 3 (415-428) Rabbit pAb	IB	canine, hamster, human, rat	382154	100 μ g	€346
Anti-Histone Deacetylase 4 (129-137) Rabbit pAb	IB	human, mouse	382161	50 μ g	€301
Anti-Histone Deacetylase 5 (194-206 and 536-545) Rabbit pAb	IB	human, mouse (weakly), rat (weakly)	382162	50 μ g	€301
Anti-Histone Deacetylase 6 (1-16) Rabbit pAb	IB	human, mouse	382163	100 μ g	€384
PhosphoDetect™ Anti-Histone H1 Rabbit pAb	IB, IC; NOT IP	bovine, human	382156	50 μ g	€221
Anti-Histone H3 (1-20) Rabbit pAb	IB, IC	human, mouse, rat	382157	50 μ g	€309

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-IKIP1 Goat pAb Cat. No. AP1047

Lane 1: MWM

Lane 2: Detection of human IKIP1 by immunoblotting. Sample: Whole cell lysate (50 μ g) from U2OS cells. Primary antibody: Anti-IKIP1 Goat pAb (Cat. No. AP1047) (2 μ g/ml). Detection: Chemiluminescence.


Lane 3: Detection of human IKIP1 by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from U2OS cells. Antibody for IP: Anti-IKIP1 Goat pAb (Cat. No. AP1047) (5 μ g/500 μ g total protein). Immunoblotting conditions: Sample: One-half of the total amount of immunoprecipitated IKIP1. Primary antibody: Anti-IKIP1 Goat pAb (Cat. No. AP1047) (2 μ g/ml). Detection: Chemiluminescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-Histone H3 (pSer ¹⁰) (7-20) Rabbit pAb	IB, IC, IH, IP	<i>C. elegans</i> , <i>Drosophila</i> , human, <i>Tetrahymena</i> , <i>Xenopus</i>	382159	50 µg	€221
Anti-Histone H3, Acetylated (1-20) Rabbit pAb	IB, IC, IP	human, <i>Tetrahymena</i>	382158	50 µg	€309
Anti-Histone H4, Acetylated (2-19) Rabbit pAb	IB, IC, IP	human, <i>Tetrahymena</i>	382160	50 µg	€301
Anti-HMG-1/HMG-2 Mouse mAb (854E10/D11)	GS, IB	bovine, canine, hamster, human, mouse, rabbit, rat, sheep	387800	100 µl	€384
Anti-HNE Fluorophore Rabbit pAb	ELISA, IB, IH	a broad range of species	393206	100 µl	€384
Anti-HNE-Michael Adducts, Reduced Rabbit pAb	ELISA, FS, IB, PS	a broad range of species	393207	100 µl	€409
Anti-HOCl-Oxidized Low Density Lipoprotein Rabbit pAb	ELISA, FS, PS	a broad range of species	428035	100 µl	€350
Anti-HPV 16-E7 Mouse mAb (TVG-701Y)	FC, IF		DP18	100 µg	€338
Anti-HPV 18/16-E6 Mouse mAb (C1P5)	IB, IF, IP, PS		DP12	100 µg	€338
Anti-HPV 18/16-E6 Mouse mAb (C1P5)	IB, IF, IP, PS		DP12L	100 µg	€338
Anti-HSF1 Rabbit pAb	ELISA, IB, IP	guinea pig, human, monkey, mouse, porcine, rat	385580	25 µl	€231
Anti-Hsp25 Rabbit pAb	IB, IC, IH, IP	bovine, canine, guinea pig, hamster, mouse, rat	385877	25 µl	€231
PhosphoDetect™ Anti-Hsp27 (pSer ¹⁵) (10-21) Rabbit pAb	IC, PS; NOT IB	human, rat	ST1005	100 µg	€350
PhosphoDetect™ Anti-Hsp27 (pSer ⁸²) Mouse mAb (5B9)	IB	human	ST1030	1 Set	€373
Anti-Hsp27 Mouse mAb (G3.1)	IC, IH, IP	human, monkey	386011	50 µg	€291
Anti-Hsp27 Rabbit pAb	ELISA, IB, IC, IP	human, monkey	386013	50 µl	€290
Anti-Hsp27 Rabbit pAb	IB	human	CA1015	50 µl	€189
Anti-Hsp40 Rabbit pAb	IB, IP	bovine, canine, chicken, fish, guinea pig, hamster, human, monkey, mouse, porcine, rat, sheep, <i>Xenopus</i>	386020	25 µl	€319
Anti-Hsp47 Mouse mAb (M16.10A1)	EIA, ELISA, IB, IC, IF, IP	bovine, canine, chicken, guinea pig, hamster, human, monkey, mouse, porcine, rat, sheep	386023	50 µg	€319
Anti-Hsp60 Mouse mAb (LK-1)	EIA, ELISA, FC, IB, IC, IH, IP	bovine, canine, <i>Drosophila</i> , fish, guinea pig, hamster, human, monkey, mouse, porcine, rabbit, rat, yeast	386028	50 µg	€290
Anti-Hsp60 Mouse mAb (LK-2)	ELISA, IB, IC	a broad range of species	386029	50 µg	€290
Anti-Hsp70 Mouse mAb (C92F3A-5)	ELISA, FC, FS, IB, IC, IP, PS	a broad range of species	386032	50 µg	€307
Anti-Hsp70 Rabbit pAb	ELISA, IB, IC, IP	a broad range of species	386035	25 µl	€290
Anti-Hsp70/Hsc70 Mouse mAb (N27F3-4)	FC, IB, IH, IP	a broad range of species	386037	50 µg	€291
Anti-Hsp72/73 (Ab-1) Mouse mAb (W27)	IB, IF, IP, PS	human, mouse, rat	HSP01	100 µg	€338
Anti-Hsp90 Mouse mAb (AC88)	IB, IC, IP	a broad range of species	386040	50 µg	€291
Anti-Hsp90 Rat mAb (16F1)	ELISA, IB, IC, IH, IP	a broad range of species	386041	50 µg	€291
Anti-Hsp90α Rabbit pAb	IB	human	CA1016	50 µl	€189
Anti-Hsp90α Rat mAb (9D2)	IB, IC, IP	human	386043	50 µg	€290
Anti-Hsp90β Rabbit pAb	IB	human	CA1017	50 µl	€189
Anti-5-HT _{1A} -Receptor (248-262) Guinea Pig pAb	FS	mouse, rat	387000	100 µl	€362
Anti-8-Hydroxyguanosine Goat pAb	PS	a broad range of species	DR1001	100 µl	€319
Anti-Hypoxia-Inducible Factor 1α Mouse mAb (H1α67)	IB, IP, PS	ferret, human, monkey, mouse, rat, sheep	400080	100 µl	€390
Anti-IκBα (Ab-1) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP	human, mouse	PC381	100 µg	€243
Anti-IκBα (Ab-2) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	PC142	100 µl	€390
PhosphoDetect™ Anti-IκBα (pSer ^{32/36}) Mouse mAb (39A1431)	IB, IP	human	OP142	100 µg	€335
Anti-IκBα Rabbit pAb	IB, IC, IP	human, mouse, rat	400001	100 µl	€344
Anti-IκBα, Cleaved Mouse mAb (5D1623)	IB	human	OP197	100 µg	€325
Anti-IκBβ (Ab-1) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	PC143	100 µl	€390

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-ICAD (2-21) Rabbit pAb	IB	mouse	400100	100 µg	€319
Anti-ICAM-1 (Ab-2) Mouse mAb (W-CAM-1)	CIA, FC, FS, IB, IF, IP; NOT PS	human	CP53	100 µg	€338
NEW Anti-ICMT (170-279) Rabbit pAb	IB	human, mouse, rat	CA1031	100 µl	€271
Anti-IDE, N-Terminal (97-273) Rabbit pAb	ELISA, IB, IC	hamster, human, mouse, rat	PC730	100 µl	€378
NEW Anti-IDE/Insulysin Rabbit pAb	IB, IP	human	ST1120	50 µl	€135
PhosphoDetect™ Anti-IFNAR1 (pSer ^{535,539}) Rabbit pAb	DB, IB	bovine, human, mouse, rat	PS1003	100 µl	€385
PhosphoDetect™ Anti-IFNAR1 (pSer ^{535,539}) Rabbit pAb	DB, IB	bovine, human, mouse, rat	PS1003	100 µl	€385
Anti-IGFBP-3 Mouse mAb (84728.111)	ELISA, IB, PS	human	GF60	500 µg	€420
Anti-IGF-I (Ab-2) Goat pAb	ELISA, IB, NT	human	PC195L	100 µg	€294
Anti-IGF-I Receptor (Ab-1) Mouse mAb (αIR3)	IF, IP; NOT IB, NOT PS	human	GR11	100 µg	€338
Anti-IGF-I Receptor (Ab-1) Mouse mAb (αIR3)	IF, IP, NT; NOT IB, NOT PS	human	GR11L	100 µg	€338
Anti-IGF-I Receptor (Ab-1) Mouse mAb (αIR3)	IF, IP; NOT IB, NOT PS	human	GR11T	10 µg	€86
Anti-IGF-I Receptor (Ab-2) Goat pAb	ELISA, IB, NT, PS	human	PC196L	100 µg	€390
Anti-IGF-I Receptor (Ab-3) Mouse mAb (33255.111)	ELISA, IB, NT	human	GR31L	100 µg	€260
Anti-IGF-I Receptor (Ab-4) Mouse mAb (24-31)	ELISA, IP; NOT IB	human, rabbit (weakly); not rat	GR38	100 µg	€338
Anti-IKIP1 Goat pAb	IB, IP	human	AP1047	100 µl	€145
Anti-IKKα Mouse mAb (14A231)	IB	human, monkey, mouse	OP133	100 µg	€315
PhosphoDetect™ Anti-IKKα/β (pSer ^{176/180}) Rabbit pAb	IB	human	ST1078	50 µl	€199
Anti-IKKβ Mouse mAb (10AG2)	IB	human	OP134	100 µg	€255
Anti-IKKγ (400-416) Rabbit pAb	IB	human, mouse, rat	400003	100 µg	€243
PhosphoDetect™ Anti-IKKγ/NEMO (pSer ³⁷⁶) Rabbit pAb	IB	human	ST1079	50 µl	€199
Anti-Inositol Hexakisphosphate Kinase (216-233) Rabbit pAb	IB, IP	rat	407265	100 µg	€362
Anti-Insulin Mouse mAb (E2E3C2)	PS; NOT IB, NOT IF, NOT IP	human; not mouse, not rat	GF44	100 µg	€338
Anti-Insulin Receptor (β-Subunit) Mouse mAb (CT-3)	FS, IB, PS; NOT IP	human, mouse, rat	GR36	100 µg	€338
Anti-Insulin Receptor (Ab-1) Mouse mAb (CII 25.3)	IP; NOT IB, NOT NS, NOT PS	human, Old World monkeys; not mouse, not New World monkeys, not rat	GR05	100 µg	€338
Anti-Insulin Receptor (Ab-3) Mouse mAb (29B4)	IP, NT; NOT IB, NOT PS	human, mouse, rat	GR07	100 µg	€338
PhosphoDetect™ Anti-Insulin Receptor (pTyr ^{1162/1163}) Rabbit pAb	ELISA, IB, IC	human	407707	10 T	€413
PhosphoDetect™ Anti-Insulin Receptor (pTyr ⁸⁷²) Rabbit pAb	IB	human, mouse	GF1000	10 T	€373
Anti-Integrin α _{IIb} Mouse mAb (283.16B7)	ELISA, FC, PS	human	407313	100 µl	€301
Anti-Integrin α _v (Ab-1) Mouse mAb (272-17E6)	CIA, ELISA, FC, IP, PS	human	407286	100 µg	€301
Anti-Integrin α _v (Ab-5) Mouse mAb (313.6F8)	CIA, ELISA, FC, IB, PS	human	407311	100 µl	€301

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCCELL:** t cell co-stimulation; **VS:** vibratome sections


**NEW****Anti MAVS Rabbit pAb, Cat No. ST1116**

Lane 1: Detection of human MAVS by immunoblotting. Sample: Whole cell lysate (20 µg) from Saos-2. Primary antibody: Anti-MAVS (Cat. No. ST1116) (0.1 µg/ml). Detection: Chemiluminescence

Lane 2: Detection of MAVS by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from Saos-2 cells. Antibody for immunoprecipitation: Anti-MAVS pAb (Cat. No. ST1116) (5 µg/500 µg total protein). Immunoblotting conditions: Sample: half of the immunoprecipitated MAVS. Primary antibody: Anti-MAVS Rabbit pAb (Cat. No. ST1116) (0.1 µg/ml). Detection: Chemiluminescence

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Integrin $\alpha 3$ Mouse mAb (P1B5)	CIA, FC, IF, IP; NOT IB, NOT PS	human	CP11L	100 μ g	€338
Anti-Integrin $\alpha 5$ /FnR Mouse mAb (P1D6)	CIA, FC, IF, IP	human	CP12L	100 μ g	€338
Anti-Integrin β_3 (Ab-1) Mouse mAb (290.5C10)	ELISA, FC, IB, IP, PS	human	407314	100 μ l	€301
Anti-Integrin β_5 Mouse mAb (343.11D1)	ELISA, FC, IB, IP, PS	human	407316	100 μ l	€301
Anti-Integrin β_6 Mouse mAb (442.5C4)	ELISA, FC, IB, IP, PS	human	407317	100 μ l	€301
Anti-Integrin $\beta 1$ Mouse mAb (4B7)	FS, IF, IP, PS; NOT IB	human; not mouse, not rat	CP26	100 μ g	€338
Anti- α -Interferon Mouse mAb (MMHA-2)	ELISA, IB, IH, IP, NT	human	407290	500 μ g	€330
Anti- α -Interferon Rabbit pAb	ELISA, IB, NT	mouse	407291	20 KU*	€319
Anti- β -Interferon Rabbit pAb	ELISA, IB, IH, IP, NT	human	407299	20 KU	€260
Anti- β -Interferon Rabbit pAb	ELISA, IB, IP, NT, PS	mouse	407307	20 KU	€319
Anti- γ -Interferon Mouse mAb (MMHG-1)	ELISA, FS, IB, IP, NT, PS	human	407309	500 μ g	€319
Anti- γ -Interferon Rat mAb (RMMG-1)	ELISA, IB, IH, IP, NT	mouse	407310	500 μ g	€319
Anti- α -Interferon Receptor Chain 2 Mouse mAb (MMHAR-2)	FC, FS, IP, NT, PS	human	407295	50 μ g	€319
Anti-Interleukin-1 Receptor-Associated Kinase-M (581-596) Rabbit pAb	IB	human, mouse, rat	407160	100 μ g	€243
Anti-Cleaved Interleukin-1 β (Asp ¹¹⁶) Rabbit pAb	IB	human	AP1016	50 μ l	€278
Anti-Interleukin-6 Rabbit pAb	ELISA, IB, NT, RIA	human	407670	1 ml	€371
Anti-IP ₃ Receptor (1829-1848) Rabbit pAb	IB, IH, IP	canine, human, mouse, rat	NR09	100 μ l	€304
Anti-IP ₃ Receptor Mouse mAb (IPR.1)	ELISA, FC, IF, IH, IP, RIA	bovine, canine, human, mouse, porcine, rabbit, rat	407140	50 μ g	€361
Anti-IP ₃ Receptor Rabbit pAb	ELISA, FC, IB, IF, IH, IP, RIA	bovine, canine, human, mouse, porcine, rat	407143	100 μ g	€390
Anti-IP ₃ Receptor, Type 1 (2732-2749) Rabbit pAb	IB, IH, IP	human, mouse, rat	407144	100 μ l	€384
PhosphoDetect™ Anti-IR/IGF1R (pTyr ^{1158/1162/1163}) Rabbit pAb	IB	human, mouse	PS1009	100 μ l	€350
PhosphoDetect™ Anti-IR/IGF1R (pTyr ¹¹⁵⁸) Rabbit pAb	IB	human	407706	10 T	€373
PhosphoDetect™ Anti-IRS (pTyr ⁸⁹⁶) Rabbit pAb	ELISA, IB	human	PC463	25 μ l	€338
Anti-IRS-1 (900-1235) Rabbit pAb	IB, IP	human, mouse, rat	420292	100 μ l	€351
PhosphoDetect™ Anti-IRS-1 (pSer ³⁰⁷) Rabbit pAb	IB	mouse	GF1005	100 μ g	€350
PhosphoDetect™ Anti-IRS-1 (pSer ⁶¹⁶) Rabbit pAb	IB	human	GF1001	10 T	€373
PhosphoDetect™ Anti-IRS-1 (pSer ⁶³⁶) Rabbit pAb	IB	mouse	GF1007	100 μ g	€350
PhosphoDetect™ Anti-IRS1 (pTyr ⁶¹²) Rabbit pAb	IB	human	PS1010	10 T	€350
PhosphoDetect™ Anti-IRS-1 (pTyr ⁸⁹⁶) Rabbit pAb	IB	human, rat	GF1003	10 T	€373
PhosphoDetect™ Anti-IRS-1 (pTyr ⁹⁴¹) Rabbit pAb	IB	human, rat	GF1004	10 T	€373
Anti-IRS-1 Mouse mAb (8-36)	IB, IP	human, mouse, rat	GR45	100 μ g	€338

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ousterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections




NEW

Anti-Nestin Mouse mAb (2C13B9), Cat.No. ST1111

Detection of human nestin by immunocytochemistry. Sample: U251 cells fixed in methanol. Primary antibody: Anti-Nestin Mouse mAb (2C13B9) (Cat. No. ST1111) (1:200). Detection: Fluorescence (green).

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-IRS-2 (1122-1321) Rabbit pAb	IB, IP	hamster, human, mouse, rat	420293	100 µl	€351
Anti-Itk Mouse mAb (2F12)	IB, IP	human, mouse	OP194	100 µg	€338
Anti-Jab1/CSN5 Rabbit pAb	IB	human	ST1045	100 µg	€290
PhosphoDetect™ Anti-JAK1 (pTyr ^{1022/1023}) Rabbit pAb	IB	human, mouse; not rat	420092	10 T	€363
* 1KU=1000 neutralization units					
PhosphoDetect™ Anti-JAK2 (pTyr ^{1007/1008}) Rabbit pAb	IB	human, mouse, rat	PS1014	100 µl	€350
PhosphoDetect™ Anti-JNK (pThr ¹⁸³ , Tyr ¹⁸⁵) Rabbit pAb	ELISA	human	PC452	25 µl	€282
PhosphoDetect™ Anti-c-Jun (pSer ⁷³) Rabbit pAb	IB, IC, IP	human, mouse	420114	100 µl	€350
Anti-c-Jun/AP-1 (Ab-2) (73-87) Rabbit pAb	FS, GS, IB, IF, IP, PS	human, mouse, rat	PC07	100 µg	€315
Anti-c-Jun/AP-1 (Ab-2) (73-87) Rabbit pAb	FS, GS, IB, IF, IP, PS	human, mouse, rat	PC07L	100 µg	€315
Anti-c-Jun/AP-1 (Ab-3) Mouse mAb (76)	IB, PS	human, mouse, rat	OP55	100 µg	€338
Anti-c-Jun/AP-1 (Ab-4) Mouse mAb (IB)	IF, PS	human, mouse, rat	OP56L	100 µg	€338
Anti-JunB (45-61) Rabbit pAb	IB, IP	human, mouse	PC28	100 µg	€319
Anti-KAIS0 (659-672) Goat pAb	IB	human, mouse	PC723	100 µg	€397
Anti-Kallikrein Rabbit pAb	IP, RIA	human; not guinea pig, not rat	420304	100 µl	€223
Anti-Kallikrein Rabbit pAb	IP, RIA	rat; not human, not rabbit	420308	100 µl	€244
Anti-KDEL Mouse mAb (10C3)	IB, IC, IP	canine, <i>Drosophila</i> , guinea pig, hamster, human, monkey, mouse, porcine, rabbit, rat, sheep, yeast	420400	50 µg	€291
Anti-KDEL Receptor Mouse mAb (KR-10)	IB, IC, IP	bovine, canine, chicken, <i>Drosophila</i> , hamster, human, monkey, mouse, porcine, rabbit, rat, sheep, <i>Xenopus</i>	421400	100 µl	€362
Anti-c-Kit (961-976) Rabbit pAb	IF, PS; NOT IB	human	PC34	100 µg	€315
Anti-Klotho Rabbit pAb	ELISA, IB	mouse	423500	50 µg	€384
NEW Anti-Kpna ₂ Rabbit pAb	IB, IP	human	ST1114	50 µl	€148
Anti-Ku (p80) (Ab-1) Mouse mAb (111)	FC, IB, IF, IP; NOT PS	human, monkey; not bovine, not mouse, not rabbit, not rat	NA52	100 µg	€338
Anti-Ku (p80) (Ab-3) Mouse mAb (S10B1)	IB, IF, IP, PS	human, monkey; not bovine, not mouse, not rabbit	NA54	100 µg	€338
Anti-Kv1.3 (471-523) Rabbit pAb	IB, IH	human, mouse, rat	421502	1 Set	€231
Anti-Cleaved Lamin A (Asp ²³⁰) Rabbit pAb	IB, IC	human, mouse, rat	AP1017	50 µl	€278
Anti-Lamin B Mouse mAb (101-B7)	FS, IB, IF; NOT PS	human; not mouse, not rat	NA12	100 µg	€390
Anti-Laminin Mouse mAb (HL-4H3)	ELISA, IB; NOT IP	human; not mouse, not rat	IF10L	100 µg	€338
Anti-Lamp1 Mouse mAb (LY1C6)	IB, IC, IF, IP	rat	428017	100 µg	€390
Anti-Lamp2 Rat mAb (GL2A7)	IB, IC, IP	mouse	428019	200 µl	€413
Anti-LDL Receptor Mouse mAb (15C8)	ELISA, FC, FS, IB, IC, IP	bovine, human, mouse; not rabbit	LP02	100 µg	€338
Anti-LEF-1 (Ab-1) Mouse mAb (REMB1)	IB, IF	human	NA64	100 µg	€338

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections




NEW

Anti-Microtubule Associated Protein 2 Mouse mAb (SMI-52)**Cat. No. ST1110**

Detection of rat microtubule associated protein 2 by staining frozen sections.
 Sample: Rat brain. Primary antibody: Anti-Microtubule Associated Protein 2 Mouse mAb (SMI-52) (Cat. No. ST1110) (1:1000). Detection: Fluorescence (green) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Leptin Rabbit pAb	ELISA, IB, RIA	mouse, rat; not human	431003	100 µl	€393
PhosphoDetect™ Anti-Leptin Receptor (pTyr ⁹⁸⁵) Rabbit pAb	IB, IP	human, mouse	431006	50 µg	€307
Anti-Leptin Receptor Rabbit pAb	ELISA, IB	mouse, rat	431005	50 µg	€384
Anti-Lewis a Mouse mAb (T174)	IH	human	434627	1 ml	€314
Anti-Lewis b Mouse mAb (T218)	ELISA, IB, IH	human	434629	1 ml	€314
Anti-Lewis x Mouse mAb (P12)	PS	human	434631	1 ml	€314
Anti-Lewis y Mouse mAb (F3)	IH	human	434636	1 ml	€314
Anti-Lipid Phosphate Phosphohydrolase 3, C-Terminal (296-310) Rabbit pAb	ELISA, IB, IF, IP, PS	human	IM1007	50 µg	€278
Anti-Lipoic Acid Rabbit pAb	ELISA, FS, IB, PS	a broad range of species	437695	100 µl	€384
Anti-5-Lipoxygenase (130-149) Rabbit pAb	IB, IC, IH	hamster, human, mouse, porcine, rat	438001	1 ea	€319
Anti-Livin (267-277) Goat pAb	IB	human, mouse (weakly)	PC724	100 µg	€397
Anti-LKB1 (120-160) Rabbit pAb	IB	human	ST1092	50 µg	€176
Anti-LMO1 (110-170) Rabbit pAb	IB, IP	human	CA1020	50 µg	€176
Anti-LRP, Heavy Chain Mouse mAb (8G1)	ELISA, FC, IB, IP, NT, PS	human	438190	100 µg	€294
Anti-LRP, Light Chain Mouse mAb (5A6)	FS, IB, PS	human, mouse, rabbit	438192	100 µg	€294
Anti-LRP-5/6 Mouse mAb (1A12)	IB	human, mouse	489475	50 µg	€301
Anti-Luciferase Mouse mAb (Luci)	IB		OB09	100 µg	€338
Anti-Lung Resistance Protein (LRP-56) (Ab-1) Mouse mAb (LRP-56)	FC, FS, IB, PS	human	438710	100 µl	€346
Anti-LXRα (318-329) Rabbit pAb	IB	human	440175	200 µg	€346
Anti-MacMARCKS Rabbit pAb	IC, IH	human, mouse, rabbit, rat	442708	100 µl	€338
Anti-MADD, C-Terminal (1570-1588) Rabbit pAb	ELISA, IB	human, mouse	441850	100 µg	€319
Anti-c-Maf Rabbit pAb	FS, GS, IB, IP, PS	mouse; not human	DR1021	100 µg	€384
Anti-MafA Rabbit pAb	FS, GS, IB, IP, PS	mouse; not human	DR1019	100 µg	€384
Anti-MafB Rabbit pAb	FS, GS, IB, IP, PS	mouse; not human	DR1020	100 µg	€384
Anti-Malondialdehyde Rabbit pAb	ELISA	a broad range of species	442730	50 µl	€294
Anti-Mannose 6-Phosphate Receptor Mouse mAb (2G11)	IB, IF, IH, IP	bovine, human, monkey; not hamster	444105	100 µg	€314
Anti-MAP Kinase ERK1 (339-353) (Ab-2) Rabbit pAb	IB, IP	human, mouse, rat	PC54	100 µg	€378
Anti-MAP Kinase ERK1/ERK2 (333-367) Rabbit pAb	IB, IH, IP	chicken, clam, <i>Drosophila</i> , frog, human, mouse, rat, sea star, sheep	442675	100 µg	€384
PhosphoDetect™ Anti-MAP Kinase ERK1/ERK2 (pThr ²⁰² /Tyr ²⁰⁴) Mouse mAb (12D4)	ELISA, IB, IP	human, mouse, rat	442706	1 Set	€389

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW


Anti N-Cadherin Mouse Antibody (13A9), Cat No. CA1029

Lane 1: Detection of human N-cadherin by immunoblotting. Sample: Whole cell lysate (25ug) from HeLa cells. Primary antibody: Anti N-Cadherin Mouse mAb (13A9) (Cat. No. CA1029) (1:1000) Detection: Chemiluminescence

Lane 2: Detection of human N-cadherin by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from HeLa cells. Antibody for immunoprecipitation: Anti N-Cadherin Mouse mAb (13A9) (Cat No. CA1029) (1:100/500 µg total protein). Immunoblotting conditions: same as lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-MAP Kinase ERK1/ERK2 (pThr ²⁰² /Tyr ²⁰⁴) Rabbit pAb	DB, IB	human, rat	442685	100 µl	€380
Anti-MAP Kinase ERK1/ERK2 Rabbit pAb	FC, IB, PS	hamster, human, mouse, rat, zebrafish	442704	200 µl	€390
Anti-MAP Kinase ERK2, N-Terminal (6-17) Rabbit pAb	IB	human, mouse, rat, <i>Xenopus</i> , yeast	442680	100 µg	€384
PhosphoDetect™ Anti-MARCKS (pSer ^{152/156}) Rabbit pAb	DB, IB	rat	442710	100 µl	€412
PhosphoDetect™ Anti-MARCKS (pSer ^{159,163,179}) Rabbit pAb	IB, IF, IH	human, mouse, rabbit, rat	442709	100 µl	€350
Anti-MARCKS Rabbit pAb	IB, IC	human, mouse, rabbit, rat	442707	100 µl	€338
Anti-MART-1 (Ab-1) Mouse mAb (A103)	IB, IP, PS	canine, human	OP176	100 µg	€338
Anti-Mast Cell Chymase Mouse mAb (CC1)	ELISA, IB, PS	human	444904	100 µg	€319
Anti-Mast Cell Tryptase Mouse mAb (AA1)	FS, IB, PS	canine, feline, human, monkey	444905	100 µg	€319
Anti-Matriptase/MT-SP1 Rabbit pAb	IB, IP	human	IM1014	50 µg	€189
NEW Anti-MAVS Rabbit pAb	IB, IP	human	ST1116	50 µg	€135
Anti-Maxi K ⁺ α (945-961) Rabbit pAb	IB, IC	human, mouse	444910	100 µg	€362
Anti-Mcl-1 (121-139) Rabbit pAb	ELISA, IB	canine, human, monkey, porcine, rabbit	444206	100 µg	€384
Anti-Mcl-1 (Ab-1) Mouse mAb (RC13)	IB, IP	human	AM50	100 µg	€338
Anti-MCM2 Goat pAb	IB	human	DR1007	100 µl	€397
Anti-MCM3 Goat pAb	IB	human	DR1008	100 µl	€397
Anti-MCM4 Goat pAb	IB	human	DR1009	100 µl	€397
Anti-MCM6 Goat pAb	IB	human	DR1010	100 µl	€397
Anti-MCM10 Rabbit pAb	IB, IP	human	DR1029	50 µg	€189
Anti-MDM2 (Ab-1) Mouse mAb (IF2)	FS, IB, IF, IP, PS	human; not mouse	OP46	20 µg 100 µg	€113 €338
Anti-MDM2 (Ab-1) Mouse mAb (IF2)	FS, IB, IF, IP, PS	human; not mouse, not rat	OP46T	10 µg	€70
Anti-MDM2 (Ab-2) Mouse mAb (2A10)	IB, IF, IP, PS	human	OP115	100 µg	€338
Anti-MDM2 (Ab-2) Mouse mAb (2A10)	IB, IF, IP, PS	human	OP115T	10 µg	€77
Anti-MDM2 (Ab-3) Mouse mAb (4B11)	IB, IF, IP, PS	human, mouse	OP143	10 µg 100 µg	€86 €338
Anti-MDM2 (Ab-4) Mouse mAb (2A9C1.18)	IB, IF, IP, PS	human, mouse	OP144	20 µg 100 µg	€113 €338
Anti-MDM2 (Ab-5) Mouse mAb (4B2C1.11)	IB, IF, IP, PS	human	OP145	100 µg	€338
Anti-MDM2 (Ab-6) Mouse mAb (5B10C)	IB, IF, IP, PS	human	OP146	100 µg	€338
PhosphoDetect™ Anti-MDM2 (pSer ¹⁶⁶) Rabbit pAb	IB	human, mouse, rat	DR1027	50 µl	€199

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW


Anti-Nestin Rabbit pAb, Cat No. ST1117

Lane 1: Detection of human nestin by immunoblotting. Sample: Whole cell lysate (20 µg) from SHSY-5Y cells. Primary antibody: Anti-Nestin Rabbit pAb (Cat. No. ST1117) (0.5 µg/ml). Detection: Chemiluminescence

Lane 2 and 3: Detection of nestin by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from (500 µg) from SHSY 5Y cells (lane 2) and negative control sample (lane 3). Antibody for immunoprecipitation: Anti-Nestin pAb (Cat. No. ST1117) (5 µg). Immunoblotting conditions: same as lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
MDM2 Antibody Sampler Kit-Human		human	ASK26	1 each	€202
PhosphoDetect™ Anti-MEK1 (pSer ²⁹⁸) Rabbit pAb	IB	human, mouse, rat	ST1015	10 T	€373
Anti-MEK1, N-Terminal (1-12) Rabbit pAb	IB, IH, IP	human, mouse, rat, <i>Xenopus</i>	444942	100 µg	€384
PhosphoDetect™ Anti-MEK1/2 (pThr ²⁹²) Rabbit pAb	DB, ELISA, IB, IH	human	ST1001	100 µl	€412
PhosphoDetect™ Anti-MEK1/2 (pThr ³⁸⁶) Rabbit pAb	DB, IB	human	ST1002	100 µl	€412
PhosphoDetect™ Anti-MEK1/MEK2 (pSer ^{218/222}) Rabbit pAb	DB, ELISA, IB	human	444955	100 µl	€363
Anti-MEK2 (1-13) Rabbit pAb	IB	human, mouse, rat, <i>Xenopus</i> , yeast	444949	100 µg	€319
PhosphoDetect™ Anti-MEK2 (pThr ³⁹⁴) Rabbit pAb	IB	human	ST1016	10 T	€373
Anti-MEK5 (59-74) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	444954	100 µg	€384
Anti-MEKK1, N-Terminal Rabbit pAb	ELISA, IB, IP	rat	453701	100 µg	€384
Anti-MEL-5 Mouse mAb (Ta99)	FS, IC	human	CA1005	200 µl	€338
Anti-Mer Rabbit pAb	IB	human	AP1025	50 µg	€189
PhosphoDetect™ Anti-c-Met (pTyr ^{1230/1234/1235}) Rabbit pAb	IB	human, mouse	PS1007	100 µl	€350
Anti-Metabotropic Glutamate Receptor 1α Rabbit pAb	IB, IF, IH, IP	feline, mouse, rat	445870	100 µl	€409
Anti-Metabotropic Glutamate Receptor 2/3 Rabbit pAb	IB, IH, IP	feline, mouse, rat	445871	100 µl	€409
Anti-Microphthalmia (Ab-1) Mouse mAb (C5)	GS, IB, IF, IP, PS	human, mouse	OP126L	100 µg	€338
Anti-Microtubule Associated Protein 2 Mouse mAb (M12)	ELISA, FS, IB, PS	chicken, human, mouse, rat	442695	100 µg	€339
Anti-Microtubule Associated Protein 2 Mouse mAb (SMI-31)	ELISA, FS, IB, IC, PS	mammalian, <i>Xenopus</i>	ST1110	100 µl	€283
PhosphoDetect™ Anti-MKK7 (pSer ²⁷¹ , Thr ²⁷⁵) Rabbit pAb	IB	human	ST1074	50 µl	€199
Anti-MKK7 Mouse mAb (3F5)	ELISA, IB	human	ST1031	100 µg	€397
Anti-MKP-7 (2-14) Goat pAb	IB	human, mouse	PC726	100 µg	€397
Anti-MLH1 (Ab-1) Mouse mAb (14)	IB, PS	human; not mouse	NA28	100 µg	€338
Anti-MLH1 (Ab-2) Rabbit pAb	IB	human, mouse	PC56	100 µg	€315
MMP/TIMP Family Antibody Sampler Kit I			ASK10	1 each	€461
Anti-MMP-1 (Ab-1) Mouse mAb (41-1E5)	FS, IB, PS	human; not bovine	IM35	20 µg 100 µg	€113 €338
Anti-MMP-1 (Ab-1) Mouse mAb (41-1E5)	FS, IB, PS	human; not bovine	IM35T	10 µg	€70
Anti-MMP-1 (Ab-3) Sheep pAb	IB	human	PC311	200 µg	€319
Anti-MMP-1 (Ab-6) Mouse mAb (III7)	IB, IP; NOT NS	human	IM67	100 µg	€338
Anti-MMP-1 (Ab-6) Mouse mAb (III7)	IB, IP; NOT NS	human	IM67T	10 µg	€77
Anti-MMP-1 Rabbit pAb	ELISA, IB, IF, IH, IP	human, marmoset, rabbit	444209	500 µl	€319
Anti-MMP-2 (468-483) (Ab-3) Mouse mAb (42-5D11)	FS, IB, PS	bovine, human, mouse, rat	IM33T	10 µg	€70
Anti-MMP-2 (Ab-11) Mouse mAb (6F7)	PS; NOT IB	rat	IM85	100 µg	€319
Anti-MMP-2 (Ab-3) Mouse mAb (42-5D11)	FS, IB, PS	bovine, human, mouse, rat	IM33	20 µg 100 µg	€113 €338
Anti-MMP-2 (Ab-4) Mouse mAb (75-7F7)	FS, IB, PS	bovine, human, mouse; not rat	IM51	20 µg 100 µg	€86 €338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-Myelin Basic Protein Mouse mAb (SMI-99) Cat. No. NE1019

Detection of rat myelin basic protein by staining frozen sections. Sample: Rat brain. Primary antibody: Anti-Myelin Basic Protein Mouse mAb (SMI-99) (Cat. No. NE1019) (1:1000). Detection: Fluorescence (green) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-MMP-2 (Ab-5) Sheep pAb	IB	human	PC312	200 µg	€319
Anti-MMP-2 (Ab-6) Mouse mAb (CA-4001)	FS, IB, IP, PS	human	IM53	100 µg	€338
Anti-MMP-2 (Ab-7) Rabbit pAb	IB, IP; NOT FS	human, mouse	PC342	100 µg	€319
Anti-MMP-2 (Ab-8) Mouse mAb (VB3)	IB, IF	human	IM68	100 µg	€338
Anti-MMP-2 (Ab-9) Mouse mAb (8B4)	ELISA, FS, IB, PS	chicken, human, mouse, rat	IM83	100 µg	€319
Anti-MMP-3 (Ab-1) Mouse mAb (55-2A4)	FS, IB, PS	human, porcine, rat, sheep; not bovine, not mouse, not rabbit	IM36L	100 µg	€338
Anti-MMP-3 (Ab-5) Mouse mAb (SL-1 IIIC4)	AFF PUR, FS, IB, IF, IP, PS	human	IM70	100 µg	€338
Anti-MMP-3 (Ab-5) Mouse mAb (SL-1 IIIC4)	AFF PUR, FS, IB, IF, IP, PS	human	IM70T	10 µg	€77
Anti-MMP-7 (Ab-1) Mouse mAb (141-7B2)	IB, PS	human	IM40	100 µg	€338
Anti-MMP-7 (Ab-2) Mouse mAb (176-5F12)	IB	human	IM47	100 µg	€338
Anti-MMP-7 (Ab-3) Mouse mAb (ID2)	AFF PUR, IB, IF, PS	human	IM71	100 µg	€338
Anti-MMP-7 (Ab-3) Mouse mAb (ID2)	AFF PUR, FS, IB, IF, PS	human	IM71T	10 µg	€77
Anti-MMP-7 (Ab-4) Rabbit pAb	IB	human, mouse, rat	PC492	50 µg	€253
Anti-MMP-8 (Ab-1) Mouse mAb (115-13D2)	IB, PS	human	IM38L	20 µg 100 µg	€89 €338
Anti-MMP-8 (Ab-2) Rabbit pAb	FS, IB, IP, PS	human	PC493	50 µg	€253
Anti-MMP-9 (626-644) (Ab-3) Mouse mAb (56-2A4)	FS, IB, PS	guinea pig, human, rabbit, rat; not bovine, not hamster, not mouse	IM37	100 µg	€328
Anti-MMP-9 (Ab-1) Mouse mAb (6-6B)	IB, IP, NT; NOT PS	human	IM09L	20 µg 100 µg	€109 €338
Anti-MMP-9 (Ab-1) Mouse mAb (6-6B)	IB, IP, NT; NOT PS	human	IM09T	10 µg	€70
Anti-MMP-9 (Ab-3) Mouse mAb (56-2A4)	FS, IB, PS	guinea pig, human, rabbit, rat; not bovine, not hamster, not mouse	IM37T	20 µg	€113
Anti-MMP-9 (Ab-5) Sheep pAb	IB	human	PC309	200 µg	€319
Anti-MMP-9 (Ab-7) Mouse mAb (GE-213)	FS, IP	human	IM61	100 µg	€338
Anti-MMP-9 (Ab-7) Mouse mAb (GE-213)	FS, IP	human	IM61T	10 µg	€77
Anti-MMP-9 (Ab-8) Mouse mAb (IA5)	IB, IF; NOT FS, NOT PS	human	IM72	100 µg	€338
Anti-MMP-9 (Ab-8) Mouse mAb (IA5)	FS, IB, IF, PS	human	IM72T	10 µg	€77
Anti-MMP-9 (Ab-9) Mouse mAb (ID4)	IB; NOT IF	rat	IM73	100 µg	€338
Anti-MMP-9 (Ab-10) Mouse mAb (VIIC2)	IF, NT; NOT IB	rat	IM74	100 µg	€338
Anti-MMP-9 Rabbit pAb	ELISA, IB, IF, IH, IP	human, marmoset	444236	500 µl	€319
Anti-MMP-10 (Ab-2) Mouse mAb (IB1)	IB; NOT FS, NOT PS	human	IM75	100 µg	€338
Anti-MMP-11 (Ab-2) Mouse mAb (SL3-05)	FS, IB, PS	human	IM86	100 µg	€338
Anti-MMP-12 (Ab-1) Rabbit pAb	IB; NOT FS, NOT IP, NOT PS	human	PC494	50 µg	€253
Anti-MMP-12 (Ab-2) Rabbit pAb	IB	human	IM1013	50 µl	€189
Anti-MMP-13 (Ab-1) Mouse mAb (181-15A12)	IB, PS	human	IM44L	20 µg 100 µg	€86 €338
Anti-MMP-13 (Ab-3) Mouse mAb (181-14G11)	IB, PS	human	IM64L	100 µg	€338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-Nestin Mouse mAb (2C13B9), Cat.No. ST1111

Detection of human nestin by immunocytochemistry. Sample: U251 cells fixed in methanol. Primary antibody: Anti-Nestin Mouse mAb (2C13B9) (Cat. No. ST1111) (1:200). Detection: Fluorescence (green).

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-MMP-13 (Ab-4) Mouse mAb (VIII A2)	AFF PUR, FS, IB, PS	human, mouse, rabbit	IM78T	10 µg	€77
Anti-MMP-13 (Ab-4) Mouse mAb (VIII A2)	AFF PUR, FS, IB, PS	human, mouse, rabbit	IM78	100 µg	€338
Anti-MMP-13 (Ab-6) Mouse mAb (VIIC4)	IB, IF; NOT NS	rat; not mouse	IM80	100 µg	€338
Anti-MMP-13 (Ab-7) Mouse mAb (LIPCO IID1)	IB, IF	mouse, rat	IM81	100 µg	€338
Anti-MMP-13 (Ab-8) Rabbit pAb	IB, PS	human, not mouse	PC542	100 µg	€294
Anti-MMP-20 (Ab-2) Rabbit pAb	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	PC497	50 µg	€253
Anti-c-Mos (149-177) Rabbit pAb	ELISA, IB	mouse, rat, <i>Xenopus</i>	475931	100 µg	€390
Anti-hMre-11 Rabbit pAb	IB, IF, IP	human, mouse, rat	PC388	100 µg	€319
Anti-hMre-11 Rabbit pAb	IB, IF, IP	human, mouse, rat	PC388T	10 µg	€77
Anti-MRP (QCRL-1) Mouse mAb (QCRL-1)	FC, IB, IC, IH	human	475720	500 µl	€390
Anti-MRP Mouse mAb (QCRL-3)	FC, IC, IH	human; not mouse	475723	500 µl	€390
Anti-MSH2 (Ab-1) Mouse mAb (GB12)	FS, IB, IP, PS	human; not mouse	NA26	100 µg	€338
Anti-MSH2 (Ab-1) Mouse mAb (GB12)	FS, IB, IP, PS	human; not mouse	NA26T	10 µg	€77
Anti-MSH2 (Ab-2) Mouse mAb (FE11)	FS, IB, IP, PS	human, mouse	NA27	100 µg	€338
Anti-MSH2 (Ab-2) Mouse mAb (FE11)	FS, IB, IP, PS	human, mouse	NA27T	10 µg	€77
Anti-MSH2 (Ab-3) Rabbit pAb	FS, IB, IP, PS	human, mouse	PC57	100 µg	€315
Anti-MT1-MMP (160-173) (Ab-2) Rabbit pAb	IB	human; not mouse, not rat	PC94L	100 µg	€319
Anti-MT1-MMP (319-333) (Ab-4) Mouse mAb (113-5B7)	IB, IF, IP, PS	human, mouse, rabbit, rat	IM57T	20 µg	€75
Anti-MT1-MMP (Ab-1) Mouse mAb (114-IF2)	IB, IP	human	IM39	100 µg	€324
Anti-MT1-MMP (Ab-3) Mouse mAb (114-6G6)	IB, IC, PS; NOT FS, NOT IP	human	IM42	100 µg	€338
Anti-MT1-MMP (Ab-4) Mouse mAb (113-5B7)	IB, IF, IP, PS	human, mouse, rabbit, rat	IM57	100 µg	€338
Anti-MT2-MMP (281-294) (Ab-1) Mouse mAb (162-22G5)	IB, IC, PS	human, mouse	IM48L	100 µg	€338
Anti-MT3-MMP (Ab-1) Mouse mAb (117-10C6)	IB	human	IM50L	100 µg	€338
Anti-MTA3 (400-450) Rabbit pAb	IB, IC, IP	human	IM1012	50 µg	€189
Anti-Munc-18 (58-70) Rabbit pAb	IB, IP	human, rat	475948	100 µg	€346
Anti-Muscarinic M1 Acetylcholine Receptor (227-353) Rabbit pAb	FS, IB, IP	human, rat	475923	1 Set	€206
Anti-Muscarinic M2 Acetylcholine Receptor (225-356) Rabbit pAb	FFS, FS, IB, IP; NOT PS	human, mouse, rat	475924	1 Set	€206
Anti-c-Myc (Ab-1) Mouse mAb (9E10)	FC, FS, IB, IF, IP; NOT PS	human, mouse (weakly), rat (weakly)	OP10	200 µg	€338
Anti-c-Myc (Ab-1) Mouse mAb (9E10)	FC, FS, IB, IF, IP; NOT PS	human, mouse (weakly), rat (weakly)	OP10L	200 µg	€338
Anti-c-Myc (Ab-1) Mouse mAb (9E10)	FC, FS, IB, IF, IP; NOT PS	human, mouse (weakly), rat (weakly)	OP10T	10 µg	€77
Anti-c-Myc (Ab-1) Mouse mAb (9E10) Fluorescein Conjugate	FC, IF	human, mouse (weakly), rat (weakly)	OP10F	100 µg	€346
Anti-c-Myc (Ab-2) Mouse mAb (8)	IB, IP; NOT PS	avian, human, mouse	OP30	100 µg	€338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-Pan-Axonal Neurofilament Marker Mouse mAb (SMI-312) Cat. No. NE1016

Detection of rat pan-axonal neurofilament marker by staining frozen sections. Sample: Rat brain. Primary antibody: Anti-Pan-Axonal Neurofilament Marker Mouse mAb (SMI-312) (Cat. No. NE1016) (1:1000). Detection: Fluorescence (green) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-c-Myc (Ab-2) Mouse mAb (8)	IB, IP; NOT PS	avian, human, mouse	OP30L	100 µg	€338
Anti-c-Myc (Ab-3) Mouse mAb (33)	FS, IB, IF, IP; NOT PS	avian, human, mouse	OP31	100 µg	€338
Anti-c-Myc (Ab-3) Mouse mAb (33)	FS, IB, IF, IP; NOT PS	avian, human, mouse	OP31L	100 µg	€338
PhosphoDetect™ Anti-c-Myc (pThr ⁶⁸ , Ser ⁶²) Rabbit pAb	FS, IB, IP	human, mouse	PC515	50 µl	€253
c-Myc Peptide	CT, DB		PP06	100 µg	€84
Anti-N-Myc Mouse mAb (NCM II 100)	FS, IB, IF, IP	human, mouse	OP13	100 µg	€338
Anti-N-Myc Mouse mAb (NCM II 100)	FS, IB, IF, IP	human, mouse	OP13L	100 µg	€338
Anti-Myelin Basic Protein (Ab-2) Mouse mAb (1)	ELISA, FS	bovine, human, rabbit (weakly), rat; not guinea pig	CP32	2 ml	€319
NEW Anti-Myelin Basic Protein Mouse mAb (SMI-94)	ELISA, FS, IB, IC, PS	canine, human, monkey, mouse, rabbit, rat	NE1018	100 µl	€390
NEW Anti-Myelin Basic Protein Mouse mAb (SMI-99)	ELISA, FS, IB, IC, PS	guinea pig (weakly), mammalian; not chicken, not porcine	NE1019	100 µl	€283
NEW Anti-Myelin CNPase Mouse mAb (SMI-91)	ELISA, FS, IB, IC, PS	bovine, canine, human, mouse, porcine, rat, sheep; not avian, not guinea pig	NE1020	100 µl	€283
Anti-Myelin Proteolipid Protein (Ab-2) Mouse mAb (plpc 1)	FS, IB	mammalian	NB38	100 µg	€319
Anti-Myeloperoxidase Rabbit pAb	ELISA, IEP	human	475915	1 ml	€154
Anti-MYH (33-51) Rabbit pAb	IB, IF, PS	human, rat	PC709	100 µl	€351
Anti-Myosin Rabbit pAb	FS, IB, IF, IP	guinea pig, human, mouse, rat	476126	500 µl	€319
Anti-Myosin, Smooth Muscle Heavy Chain Mouse mAb (SMMS-1)	PS	human	475985	100 µg	€319
Anti-NAK/TBK1 (712-727) Rabbit pAb	IB	human, mouse, rat	477740	100 µg	€243
Anti-NBS1 Rabbit pAb	IB, IP, NT	human	DR1033	50 µg	€189
NEW Anti-N-Cadherin Mouse mAb (13A9)	IB, IP	human	CA1029	50 µl	€135
Anti-NEIL1 Rabbit pAb	IB, IF, IP		DR1038	100 µl	€367
Anti-NEIL2 Rabbit pAb	IB	human	PC708	100 µl	€351
Anti-β-Nerve Growth Factor (Ab-2) Mouse mAb (25623.1)	ELISA, IB, NT	human, rat; not mouse	GF34L	100 µg	€319
Anti-Nerve Growth Factor Receptor Mouse mAb (192)	IF, IP; NOT IB, NOT PS	rat; not human	GR10	100 µg	€338
Anti-Nestin Mouse mAb (2C13B9)	FS, IB, IC	human	ST1111	100 µl	€135
NEW Anti-Nestin Rabbit pAb	IB, IP	human	ST1117	50 µg	€135
Anti-N-Ethylmaleimide-Sensitive Factor Mouse mAb (9G7-3)	IB, IP	bovine, canine, chicken, hamster, human, monkey, mouse, porcine, rabbit, rat, sheep	341152	100 µg	€384
Anti-N-Ethylmaleimide-Sensitive Factor Mouse mAb (NSF-1)	ELISA, IB, IP	human, rat	NB21	100 µg	€338
Anti-Netrin-1 (Ab-1) (557-568) Rabbit pAb	IB	mouse	PC344	100 µg	€338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade; **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherlon; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections


**NEW****Anti-NOD2/CARD15 Rabbit pAb, Cat No. ST1119**

Lane 1: Detection of human NOD2 by immunoblotting. Sample: Whole cell lysate (20 µg) from Jurkat cells. Primary antibody: Anti-NOD2 Rabbit pAb (Cat. No. ST1119) (2 µg/ml). Detection: Chemiluminescence

Lane 2: Detection of human NOD2 by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from Jurkat cells (lane 2) and negative control sample (lane 3). Antibody for immunoprecipitation: Anti-NOD2 Rabbit pAb (Cat. No. ST1119) (5 µg). Immunoblotting conditions: same as lane 1

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Netrin-1 (Ab-2) (497-508) Rabbit pAb	FS; NOT IB	mouse	PC364	100 µg	€338
Anti-Neural Cell Adhesion Molecule Mouse mAb (ERIC-1)	IB, IF, IP; NOT FS, NOT PS	human	CP27L	100 µg	€338
PhosphoDetect™ Anti-Neurofilament H Mouse mAb (SMI-31)	ELISA, FS, IB, IC, PS	chicken, mammalian, <i>Xenopus</i>	NE1022	100 µl	€325
Anti-Pan-Axonal Neurofilament Marker Mouse mAb	ELISA, FS, IB, IC, PS	mammalian	NE1016	100 µl	€325
Anti-Pan-Neuronal Neurofilament Marker Mouse mAb (SMI-311)	ELISA, FS, IB, IC, PS	mammalian	NE1017	100 µl	€325
Anti-Neurokinin-1 Receptor (393-407) Rabbit pAb	FS, IC	guinea pig, mouse, rat	480738	100 µl	€390
Anti-Neurokinin-3 Receptor (438-452) Rabbit pAb	FFS, FS	guinea pig, mouse, rat	480739	100 µl	€390
Anti-Neuronal Specific Enolase Rabbit pAb	FS, PS	bovine, canine, chicken, feline, guinea pig, horse, human, mouse, porcine, rabbit, sheep	480734	1 ml	€378
Anti-Neuronal Specific Enolase Rabbit pAb	FS, IEP, PS	human	PC237	6 ml	€220
Anti-Neuropeptide Y Rabbit pAb	FS	amphibian, avian, feline, fish, human, newt, rat	PC223L	100 µl	€319
Anti-Neuropilin-1 (813-827) Rabbit pAb	FS, IB	mouse, rat	PC343	100 µg	€338
Anti-Neuropilin-1 (813-827) Rabbit pAb	FS, IB	mouse, rat	PC343T	10 µg	€77
Anti-Neurotensin Rabbit pAb	FS, IF	rat	PC225L	100 µl	€259
Anti-Neutrophil Elastase Mouse mAb (121-3E10)	IB, PS; NOT FS	human	IM46	100 µg	€338
Anti-Neutrophil Elastase Rabbit pAb	ELISA, IB, IEP, IP	human	481001	1 ml	€114
Anti-NF-κB (p50) (Ab-1) Rabbit pAb	ELISA, GS, IB, IP	human; not mouse	PC136	20 µl 100 µl	€76 €330
Anti-NF-κB (p65) Mouse mAb (2A12A7)	ELISA, GS, IB	human	ST1047	100 µl	€319
Anti-NF-κB (p65, RelA) (Ab-1) Rabbit pAb	GS, IB, IP	human, mouse	PC137	20 µl 100 µl	€76 €378
Anti-NF-κB (p65, RelA) (Ab-2) Rabbit pAb	ELISA, GS, IB, IP	human	PC138	100 µl	€390
Anti-NFAT3 (887-902) Rabbit pAb	IB, IC	human, mouse	PC625	100 µg	€378
Anti-NFAT5 (1439-1455) Rabbit pAb	IB, IC	human, mouse	PC626	50 µg	€378
NF-κB Antibody Sampler Kit			ASK20	1 each	€231
Anti-NFBD1/MDC1 Rabbit pAb	IB, IP	human	DR1018	100 µg	€360
Anti-Nibrin Rabbit pAb	ELISA, IB, IP	human	PC634	50 µl	€253
Anti-Nicastrin, C-Terminal (689-709) Rabbit pAb	IB, IP, RIA	human, mouse	481905	100 µl	€397
Anti-Nicastrin, N-Terminal (62-93) Rabbit pAb	IB, IF, IP	human, mouse	481906	100 µl	€397
Anti-Nidogen Rabbit pAb	FS, IB, IF, IP	human	481978	500 µl	€319
Anti-Nitric Oxide Synthase (1131-1144), Inducible Rabbit pAb	FS, IB, IC, PS	human, mouse, rat	482728	100 µl	€413
Anti-Nitric Oxide Synthase (1411-1433), Neuronal C-Terminal Rabbit pAb	IB, IC, IF	human	482750	100 µl	€366
Anti-Nitric Oxide Synthase (2-12), Endothelial N-Terminal Rabbit pAb	IB, IC	human	482736	100 µl	€366
Nitric Oxide Synthase Antibody Set			482730	1 Set	€550
Anti-Nitric Oxide Synthase Mouse mAb (H32)	ELISA, FS, IB, IP	baboon, bovine, human, mouse, porcine, rat	482731	100 µg	€413

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ousterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-Pan Neuronal Marker Mouse mAb (SMI-311), Cat. No. NE1017

Detection of rat pan-neuronal neurofilament marker by staining frozen sections. Sample: Rat brain. Primary antibody: Anti-Pan-Neuronal Neurofilament Marker Mouse mAb (SMI-311) (Cat. No. NE1017) (1:1000). Detection: Fluorescence (green) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-Nitric Oxide Synthase, Endothelial (pSer ¹¹⁷⁷) Mouse mAb (15E2)	ELISA, IB	human, mouse	482737	100 µg	€350
Anti-Nitric Oxide Synthase, Inducible, C-Terminal (SLEMSAL) Rabbit pAb	IB, IC, PS	human; not mouse, not rat	482755	200 µl	€390
Anti-Nitrotyrosine Mouse mAb (CC22.8C7.3)	IB	a broad range of species	487923	50 µg	€191
Anti-Nitrotyrosine Rabbit pAb	IB	a broad range of species	487924	150 µg	€409
Anti-nm23 (Ab-1) Mouse mAb (NM301)	FS, IB, IP, PS	human; not rat	OP48	100 µg	€346
PhosphoDetect™ Anti-NMDA NR1 (pSer ⁸⁹⁶) Rabbit pAb	IB, IC	human, mouse, rat	454571	50 µl	€264
Anti-NMDA NR1, Splice Variant N1 Rabbit pAb	ELISA, IB	rat	454576	25 µg	€409
PhosphoDetect™ Anti-NMDA NR2B (pTyr ¹⁴⁷²) Rabbit pAb	DB, IB	rat	454583	100 µl	€412
Anti-Pan-NMDA NR1 Mouse mAb (R1JHL)	IB, IP	rat	454578	15 µg	€409
Anti-NMDAR 2A (1252-1391) Rabbit pAb	FS, IB, IP	mouse, rat	454580	10 µg	€420
Anti-NMDAR 2B (984-1104) Rabbit pAb	FS, IB, IP	fish, human, mouse, rat	454582	10 µg	€420
Anti-Pan-NMDAR 2C (25-130) Rabbit pAb	IB, IP	rat	454584	10 µg	€420
NEW Anti-NOD2/CARD15 Rabbit pAb	IB, IP	human	ST1119	50 µg	€135
Anti-Notch 1 Hamster mAb (8G10)	FC, IB, IC	mouse, rat	491010	50 µl	€231
Anti-Cleaved Notch 1 (Val ¹⁷⁴⁴) Rabbit pAb	ELISA, IB, IC, IP	human, monkey, mouse, rat	ST1028	50 µl	€278
Anti-Notch 1, Cytoplasmic Domain (2272-2488) Rabbit pAb	IB	mouse, rat	491011	50 µg	€231
Anti-Noxa Mouse mAb (114C307)	IB	human	OP180	100 µg	€351
Anti-ΔNp73 Mouse mAb (38C674)	IB	human	OP181	100 µg	€351
Anti-NRAGE Rabbit pAb	IB, IC, IP	human, mouse, rat	AP1006	100 µl	€334
Anti-NuMA (Ab-1) Mouse mAb (204-41)	IB, IF	human; not mouse, not rat	NA08	100 µg	€378
Anti-NuMA (Ab-2) Mouse mAb (107-7)	FS, IB, IF	human; not mouse, not rat	NA09L	100 µg	€378
Anti-Oct-1 (Ab-3) Mouse mAb (YL8)	GS, IB	human, mouse	NA43L	100 µg	€338
Anti-Oct-2 Mouse mAb (PT1)	GS, IB, IP	human; not mouse	NA44	100 µg	€338
Anti-Odin Rabbit pAb	IB, IP	human	ST1039	50 µg	€194
Anti-Omi/HtrA2 Mouse mAb (196C429)	IB	human	AM82	100 µg	€432
Anti-Opioid δ Receptor (3-17) Rabbit pAb	FS, IB, IF, IP	rat	PC163L	100 µl	€432
Anti-Opioid µ Receptor (384-398) Rabbit pAb	FS, IB, IF, IP	rat	PC165L	100 µl	€432
NEW Anti-ORC1 Rabbit pAb	IB, IP	human	DR1039	50 µg	€145
Anti-ORC2 (Ab-2) Mouse mAb (920 4-41)	IB, IP, PS	human	NA73	100 µg	€338
Anti-Orexin A (14-33) (Ab-1) Rabbit pAb	ELISA, FS, IF	bovine, human, mouse, rat	PC345	100 µg	€338


10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

**NEW****Anti-ORC1 Rabbit pAb Cat. No. DR1039**

Lane 1: Detection of human ORC1 by immunoblotting. Sample: Whole cell lysate (50 µg) from HeLa cells. Primary antibody: Anti-ORC1 Rabbit pAb (Cat. No. DR1039) (2 µg/ml). Detection: Chemiluminescence. Lane 2: Detection of ORC1 by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from HeLa cells. Antibody for IP: Anti-ORC1 Rabbit pAb (Cat. No. DR1039) (10 µg/mg protein). Immunoblotting conditions: same as above.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Orexin A (14-33) (Ab-2) Rabbit pAb	ELISA, FFS, FS, IF	bovine, human, mouse, rat	PC362	100 µl	€338
Anti-Orexin B (6-24) (Ab-1) Rabbit pAb	FS, IF	human, mouse, rat	PC346	25 µg	€338
Anti-Osteocalcin (1-22) Rabbit pAb	ELISA, RIA	human; not bovine, not sheep	499055	100 µl	€319
Anti-Osteonectin Rabbit pAb	ELISA, IB, PS	human	499255	100 µl	€319
Anti-Osteopontin (75-90) Rabbit pAb	ELISA, IB, PS	human	499265	100 µl	€319
Anti-Oxytocin Rabbit pAb	FS, IF	rat	PC226L	100 µl	€319
Anti-p14 ^{ARF} (Ab-1) (118-132) Rabbit pAb	IB	human; not mouse, not rat	PC409	100 µl	€294
Anti-p14 ^{ARF} (Ab-2) Mouse mAb (14P02)	IB, IF, IP	human; not mouse, not rat	NA70	100 µg	€338
Anti-p16 (Ab-1) Mouse mAb (DCS-50.1/H4)	IB	human	NA29	100 µg	€338
Anti-p19 ^{ARF} (54-75) (Ab-1) Rabbit pAb	IB	mouse	PC435	100 µl	€319
Anti-p19 ^{ARF} (54-75) Rat mAb (5-C3-1)	IB, IF, IP, PS	mouse	CB1012	50 µg	€145
Anti-p19 ^{ARF} (54-75) Rat mAb (12-A1-1)	IB, IF, IP	mouse	CB1013	50 µg	€145
Anti-p21 ^{WAF1} (Ab-3) Mouse mAb (DF10)	IB, IP	human; not mouse, not rat	OP68	100 µg	€338
Anti-p21 ^{WAF1} (58-77) (Ab-1) Mouse mAb (EA10) Fluorescein Conjugate	FC, IF	human; not mouse, not rat	OP64F	100 µg	€457
Anti-p21 ^{WAF1} (Ab-1) Mouse mAb (EA10)	FC, FS, IB, IF, IP, PS	human; not mouse, not rat	OP64	100 µg	€338
Anti-p21 ^{WAF1} (Ab-4) Mouse mAb (65)	IB, IP	mouse, rat	OP76	100 µg	€338
Anti-p21 ^{WAF1} (Ab-6) Mouse mAb (22)	IB, IP; NOT FS, NOT PS	human, mouse, rat	OP79	100 µg	€338
Anti-p27 ^{Kip-1} (182-198) Rabbit pAb	IB, IP, PS	human, rabbit, rat	506127	50 µg	€319
Anti-p27 ^{Kip-1} (Ab-2) Mouse mAb (DCS72)	IB, PS	human, mouse	NA35	100 µg	€338
PhosphoDetect™ Anti-p27 ^{Kip-1} (pThr ¹⁸⁷) (179-190) Rabbit pAb	IB	feline, human, mouse, rat	506128	50 µl	€350
Anti-p38 MAP Kinase δ Rabbit pAb	IB	human, monkey, mouse, rat	ST1077	50 µl	€176
Anti-p38 MAP Kinase (341-360) Rabbit pAb	IB, PS	human, mouse, rat	506123	200 µl	€390
PhosphoDetect™ Anti-p38 MAP Kinase (pThr ¹⁸⁰ , pTyr ¹⁸²) Rabbit pAb	IB, IC	canine, human, mouse, rat	506119	10 T	€373
Anti-p38/HOG1, C-Terminal (341-360) Rabbit pAb	IB, IH, IP	human, mouse, rat	506115	100 µg	€384
Anti-p40 (5-17) Rabbit pAb	IB, PS; NOT IF	human; not mouse	PC373	100 µl	€294
Anti-p51A Rabbit pAb	IB, IF; NOT FS, NOT IP	human, mouse	PC472	100 µl	€457
Anti-p53 (Ab-1) (Pantropic) Mouse mAb (PAb421)	FC, FS, GS, IB, IF, IP	human, monkey, mouse, rabbit, rat	OP03	20 µg 100 µg	€89 €338
Anti-p53 (Ab-1) (Pantropic) Mouse mAb (PAb421)	FC, FS, GS, IB, IF, IP	human, monkey, mouse, rabbit, rat	OP03L	100 µg	€338
Anti-p53 (Ab-1) (Pantropic) Mouse mAb (PAb421)	FC, FS, GS, IB, IF, IP	human, monkey, mouse, rabbit, rat	OP03T	10 µg	€70
Anti-p53 (Ab-1) (Pantropic) Mouse mAb (PAb421) Fluorescein Conjugate	FC, FS, IF	human, monkey, mouse, rabbit, rat	OP03F	100 µg	€346
Anti-p53 (Ab-2) (Pantropic) Mouse mAb (PAb1801)	FS, GS, IB, IP, PS	human; not mouse, not rat	OP09	100 µg 200 µg	€270 €378
Anti-p53 (Ab-2) (Pantropic) Mouse mAb (PAb1801)	FS, GS, IB, IP, PS	human; not mouse, not rat	OP09L	100 µg	€291
Anti-p53 (Ab-2) (Pantropic) Mouse mAb (PAb1801)	FS, GS, IB, IP, PS	human; not mouse, not rat	OP09T	10 µg	€77

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-p19^{ARF} (54-75) Rat mAb (5-C3-1) Cat. No. CB1012

Detection of mouse p19^{ARF} by immunoblotting. Samples: Whole cell lysate (30 mg) from MEF (lanes 1,3) and NIH3T3 (lanes 2,4) cells. Primary antibody: Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) (Cat. No. CB1013) (3 mg/ml) (lanes 1,2) or Anti-p19^{ARF} (54-75) Rat mAb (5-C3-1) (Cat. No. CB1012) (3 mg/ml) (lanes 3,4). Detection: Chemiluminescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-p53 (Ab-3) (Mutant) Mouse mAb (PAb240)	FC, FS, GS, IB, IF, IP, PS	bovine, chicken, hamster, human, mouse, rat; not <i>Xenopus</i>	OP29	20 µg 100 µg 200 µg	€86 €301 €378
Anti-p53 (Ab-3) (Mutant) Mouse mAb (PAb240)	FC, FS, GS, IB, IF, IP, PS	bovine, chicken, hamster, human, mouse, rat; not <i>Xenopus</i>	OP29L	100 µg	€291
Anti-p53 (Ab-4) (Wild type) Mouse mAb (PAb246)	IC, IP, PS; NOT IB	mouse, rat; not human	OP32	100 µg	€338
Anti-p53 (Ab-4) (Wild type) Mouse mAb (PAb246)	IC, IP, PS; NOT IB	mouse, rat; not human	OP32L	100 µg	€324
Anti-p53 (Ab-5) (Wild type) Mouse mAb (PAb1620)	FS, IF, IP, PS; NOT IB	bovine, human, mouse, primate, rat	OP33	20 µg 100 µg	€86 €338
Anti-p53 (Ab-6) (Pantropic) Mouse mAb (DO-1)	FS, IB, IC, IP, PS	feline, human; not mouse, not rat	OP43	20 µg 100 µg	€89 €344
Anti-p53 (Ab-6) (Pantropic) Mouse mAb (DO-1)	FS, GS, IB, IC, IP, PS	feline, human; not mouse, not rat	OP43L	100 µg	€344
Anti-p53 (Ab-6) (Pantropic) Mouse mAb (DO-1)	FS, IB, IC, IP, PS	feline, human; not mouse, not rat	OP43T	10 µg	€70
Anti-p53 (Ab-6) (Pantropic) Mouse mAb (DO-1) Agarose Conjugate	AC, IP	feline, human; not mouse, not rat	OP43A	0.5 ml	€432
Anti-p53 (Ab-6) (Pantropic) Mouse mAb (DO-1) Fluorescein Conjugate	FC, FS, IC, IF, PS	feline, human; not mouse, not rat	OP43F	100 µg	€346
Anti-p53 (Ab-7) (Pantropic) Sheep pAb	FS, IB, IF, IP, PS	human, mouse, rat	PC35	1 each	€397
Anti-p53 (Ab-8) (Pantropic) Mouse mAb (BP53.12)	FS, IB, IP, PS	human, primate; not mouse, not rat	OP73	100 µg	€346
Anti-p53 (Ab-11) (Pantropic) Mouse mAb (PAb1802)	IB, IP; NOT FS, NOT IF, NOT PS	human, mouse	OP104L	100 µg	€338
Anti-p53 (Ab-12) (Pantropic) Mouse mAb (DO-7)	FS, IB, IF, IP, PS	bovine, human, monkey	OP140	100 µg	€338
PhosphoDetect™ Anti-p53 (pSer ¹⁵) (Ab-3) Rabbit pAb	IB, IH	human, mouse, rat	PC386T	5 µl	€77
PhosphoDetect™ Anti-p53 (pSer ¹⁵) (Ab-6) Rabbit pAb	IB, PS	human	PC461	25 µl	€253
PhosphoDetect™ Anti-p53 (pSer ²⁰) Rabbit pAb	IB, IC, IF, PS	human, mouse	DR1023	50 µl	€199
PhosphoDetect™ Anti-p53 (pSer ³⁹²) (Ab-4) Rabbit pAb	IB	human, mouse	PC387	50 µl	€253
PhosphoDetect™ Anti-p53 (pSer ³⁹²) Mouse mAb (9F4)	ELISA, IB	human, mouse	506133	100 µg	€350
PhosphoDetect™ Anti-p53 (pSer ⁴⁶) Rabbit pAb	IB, IC, IP	human	DR1024	50 µl	€199
Anti-p53 Binding Protein 1 (Ab-1) Rabbit pAb	IB, IF, IP	human, mouse	PC712	100 µl	€378
Anti-p53 Binding Protein 1 Mouse mAb (BP13)	IB, IC, IP	human; not mouse	DR1003	50 µl	€194
Anti-p57 ^{kip-2} (Ab-2) Mouse mAb (TFE42)	IB; NOT FS, NOT IF, NOT IP, NOT PS	human	NA36	100 µg	€338

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) Cat. No. CB1013

A. Detection of mouse p19^{ARF} by immunoblotting. Samples: Whole cell lysate (25 µg) from NIH3T3 cells that have deleted the Arf gene (lane 1) and wild type mouse embryo fibroblasts (MEFs) at passage 6, which express p19^{ARF} (lane 2). Primary antibody: Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) (Cat. No. CB1013) (1 µg/ml). Detection: Chemiluminescence.

B. Detection of mouse p19^{ARF} by immunoprecipitation. Samples: Whole cell lysates (250 µg) from NIH3T3 cells that have deleted the Arf gene (lane 5) and NIH3T3 cells expressing p19^{ARF} with a HA-tag (lanes 3,4). Antibodies for IP: Normal rat IgG (negative control; lane 3) and Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) (Cat. No. CB1013). Conditions for immunoblotting: Samples: Whole cell lysate (25 µg, equivalent to 10% of above immunoprecipitations) from NIH3T3 cells that have deleted the Arf gene (lane 1), NIH3T3 cells expressing p19^{ARF} with a HA-tag (lane 2), and immunoprecipitated proteins as indicated above. Primary antibody: Anti-HA. Detection: Chemiluminescence.

C. Detection of mouse p19^{ARF} by immunocytochemistry. Samples: Wild type mouse embryo fibroblasts (MEF) at passage 5 (panels 1,2) and NIH3T3 cells, which have deleted the Arf gene (panels 3,4), fixed in methanol/acetone. Primary antibody: Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) (Cat. No. CB1013). Detection: Anti-Rat IgG secondary antibody, fluorescent conjugate (panels 1,3) and counterstained with Hoechst dye (panels 2,4).

D. Detection of mouse p19^{ARF} by immunoblotting. Samples: Whole cell lysate (30 µg) from MEF (lanes 1,3) and NIH3T3 (lanes 2,4) cells. Primary antibody: Anti-p19^{ARF} (54-75) Rat mAb (12-A1-1) (Cat. No. CB1013) (3 µg/ml) (lanes 1,2) or Anti-p19^{ARF} (54-75) Rat mAb (5-C3-1) (Cat. No. CB1012) (3 µg/ml) (lanes 3,4). Detection: Chemiluminescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-p63 (Ab-1) Mouse mAb (4A4)	IB, IF, PS	human, mouse, rat	OP132	10 µg 100 µg	€86 €338
Anti-p70S6 Kinase Rabbit pAb	IB	human, mouse, rat	ST1046	100 µg	€384
Anti-p73 (Ab-1) Mouse mAb (ER-13)	GS, IB, IP	human, monkey	OP108	100 µg	€338
Anti-p73 (Ab-2) Mouse mAb (ER-15)	GS, IB, IP	human, monkey	OP109	100 µg	€338
Anti-p73 (Ab-2) Mouse mAb (ER-15)	GS, IB, IP	human, monkey	OP109L	100 µg	€338
Anti-p73 (Ab-2) Mouse mAb (ER-15)	GS, IB, IP	human, monkey	OP109T	10 µg	€77
Anti-p73 (Ab-3) Mouse mAb (GC-15)	IB, IP	human, monkey	OP110	100 µg	€338
Anti-p73 (Ab-4) (450-636) Rabbit pAb	IB, IP, PS	human	PC385	50 µg	€294
Anti-p75 Neurotrophin Receptor Rabbit pAb	IB, IC	human	NE1024	100 µl	€368
Anti-p110 Mitochondrial Protein Mouse mAb (2G2)	IB, IF	human, monkey; not chicken, not mouse, not rat	AM67	100 µg	€319
Anti-P2X3 Receptor (383-397) Rabbit pAb	IC	human, mouse, rat	PC411	25 µg	€338
Anti-P2X4 Receptor (373-388) Rabbit pAb	FS, IB	human, mouse, rat	PC376	100 µg	€319
Anti-P2X7 Receptor (576-595) Rabbit pAb	IB	rat	506107	1 Set	€231
Anti-P2Y1 Receptor (242-258) Rabbit pAb	IB, IH	human, rat	512710	1 Set	€278
Anti-p300 Mouse mAb (NM-11)	IB, IF, IP, PS	human; not mouse, not rat	NA46	100 µg	€338
p53 Antibody Sampler Kit		a broad range of species	ASK07	1 each	€406
Anti-PAI-1 (Ab-1) Mouse mAb (2)	IB, IP	human	IM29L	100 µg	€338
Anti-PAI-1 (Ab-2) Mouse mAb (6)	FS, IB	human	IM30L	100 µg	€338
Anti-PAI-1 Mouse mAb (7)	IB, IP	human	IM31L	100 µg	€338
Anti-PAI-1 Mouse mAb (MA-33B8)	IB, NT	human, rabbit	528215	100 µg	€288
Anti-PAI-1 Mouse mAb (MA-33H1F7)	NT	human, mouse	528217	100 µg	€288
Anti-PAI-1 Rabbit pAb	ELISA, IB, PS	human	528216	100 µg	€319
PhosphoDetect™ Anti-PAK1/2/3 (pThr ⁴²³) Rabbit pAb	IB	human	ST1018	10 T	€373
Anti-Pan β-Amyloid (15-30) Rabbit pAb	DB, ELISA, IH, RIA	human	PC152	100 µg	€384
Anti-PAR (Ab-1) Mouse mAb (10H)	IB, IC	<i>Drosophila</i> , human, mouse, rat	AM80	100 µg	€360
Anti-PAR Rabbit pAb	ELISA, IB, IF, IP	a broad range of species	528815	100 µl	€384
Anti-Parathyroid Hormone Receptor Type 1 Mouse mAb (3D1.1)	FS, IB, PS, RIA	human	GF59	100 µg	€338
NEW Anti-PARC (1-763) Rabbit pAb	IB, IP	human	DR1040	50 µl	€135
Anti-PARC/H7-AP1 (900-950) Rabbit pAb	IB, IP	human; not mouse	DR1028	100 µg	€384
Anti-PARG (779-800) Rabbit pAb	ELISA, IB; NOT IF	bovine, human	AP1018	50 µl	€445
Anti-Parkin (81-98) Rabbit pAb	FS, PS; NOT IB	human, mouse, rat	PC372	25 µg	€338
Anti-PARP (215-228) Rabbit pAb	ELISA, IB, IF	human	512739	100 µl	€381
Anti-PARP-1 (197-214) Rabbit pAb	ELISA, IB, IC, IF, IP	human	512738	100 µl	€350
Anti-PARP-1 (Ab-1) Rabbit pAb	IB; NOT IP, NOT PS	human	PC100	100 µg	€319
Anti-PARP-1 (Ab-2) Mouse mAb (C-2-10)	IB, IF; NOT FS, NOT PS	hamster, human, monkey, mouse, rat; not chicken	AM30	100 µg	€338
Anti-PARP-1 (Ab-3) Mouse mAb (F1-23)	ELISA, IB, IF, IP	bovine, human; not mouse, not rat	AM68	25 µl	€248
Anti-PARP-1 Automodification Domain (509-524) Rabbit pAb	ELISA, IB, IP	bovine, hamster, human, monkey, mouse, rat	512737	100 µl	€351
Anti-PARP-1 Rabbit pAb	ELISA, IB, IC, IP, NT	bovine, gerbil, human, monkey, mouse, rat	512734	100 µl	€351
NEW Anti-PARP10 Rabbit pAb	IB	human	AP1048	50 µg	€135
Anti-Parvalbumin Rabbit pAb	FS, IF, VS	human, mouse, rat	PC255L	100 µl	€409
Anti-Pax3 Rabbit pAb	IB, IC	human, mouse	CA1010	50 µl	€361
Anti-Pax-5 (1-15) Rabbit pAb	ELISA, IB	human	PC700	100 µl	€409

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-Paxillin (pSer ¹⁷⁸) Rabbit pAb	IB	human, mouse	ST1069	50 µg	€273
PhosphoDetect™ Anti-Paxillin (pTyr ¹¹⁸) Rabbit pAb	IB, IC	chicken, human	512880	10 T	€389
Anti-Paxillin Mouse mAb (5H11)	IB, IF, IP, PS	human	CP58	100 µg	€338
Anti-P-Cadherin Mouse mAb (NCC-CAD-299)	CIA, ELISA, FS, IB; NOT PS	human	205607	100 µg	€390
Anti-PCNA (112-121) (Ab-1) Mouse mAb (PC10)	IB, IF, PS; NOT IP	human, mouse, rat, yeast	NA03T	10 µg	€77
Anti-PCNA (247-261) (Ab-5) Rabbit pAb	IP	human	PC474	5 µl 50 µl	€86 €294
Anti-PCNA (Ab-1) Mouse mAb (PC10)	IB, IF, PS; NOT IP	human, mouse, rat, yeast	NA03	200 µg	€338
Anti-PCNA (Ab-1) Mouse mAb (PC10) Fluorescein Conjugate	FC, IF	human, mouse, rat, yeast	NA03F	100 µg	€346
Anti-PDGF (Ab-1) (101-116) Rabbit pAb	ELISA, IB, IF, PS	human, rat	PC21	100 µg	€319
Anti-PDGF Receptor β (Ab-1) (425-446) Rabbit pAb	IP, PS; NOT IB	mouse, rat; not human	PC17	100 µg	€319
Anti-PDGF Receptor β (Ab-4) Mouse mAb (PR7212)	FS, IB, IF, IP, PS	primate; not mouse, not rat	GR23L	100 µg	€338
PhosphoDetect™ Anti-PDGF Receptor-β (pTyr ⁷¹⁶) Rabbit pAb	IB	human, mouse, rat	521252	50 µg	€301
NEW Anti-PDK1 (1-556) Rabbit pAb	IB	human	ST1115	50 µl	€135
Anti-PDK1 (285-559) Rabbit pAb	IB, IP	hamster, human, mouse	ST1036	50 µl	€208
PhosphoDetect™ Anti-PDK1 (pSer ²⁴¹) Rabbit pAb	IB, IC, IP	human, mouse, rat	ST1073	50 µl	€199
Anti-PEN-2, N-Terminal (1-26) Rabbit pAb	IB, IP	human, mouse	NE1008	50 µl	€278
Anti-PER1 (39-51) Rabbit pAb	FS, IB	mouse	516045	200 µl	€346
Anti-Perforin Mouse mAb (δG9)	FC, FS, IP, PS	human	516515	100 µg	€347
Anti-Perforin Mouse mAb (δG9) Fluorescein Conjugate	FC	human	516516	120 T	€413
Anti-PGC-1, C-Terminal (777-797) Rabbit pAb	ELISA, IB, IH	human, mouse	516557	100 µl	€385
Anti-P-Glycoprotein Mouse mAb (4E3.16)	FC, FS, IC, IP; NOT IB, NOT PS	human	517308	500 µl	€384
Anti-P-Glycoprotein Mouse mAb (C219)	FC, FS, IB, IC, IP, PS	hamster, human, mouse, rat	517310	500 µl	€384
Anti-P-Glycoprotein Mouse mAb (C494)	FC, FS, IB, IC, IP, PS	hamster, human	517312	500 µl	€409
Anti-PGP9.5 (187-202) Rabbit pAb	FS, IB	human, mouse, porcine, rat	NE1013	50 µl	€212
Anti-Phosphodiesterase V Rabbit pAb	IB, IP	bovine	524583	50 µl	€409
PhosphoDetect™ Anti-Phospholamban (pSer ¹⁶) Rabbit pAb	IB	bovine, human	525155	50 µg	€307
Anti-Phospholamban Mouse mAb (A1)	FS, IB, IP	bovine, canine, mouse, porcine, rat	525154	50 µg	€470
Anti-Phospholipase A ₂ , Type VI (557-576) Rabbit pAb	IB	hamster, human, mouse	525162	200 µl	€384
PhosphoDetect™ Anti-Phospholipase Cγ1 (pTyr ⁷⁸³) Rabbit pAb	IB, IC	bovine, human, mouse, rat, <i>Xenopus</i>	ST1020	10 T	€350
Anti-Phosphoserine Mouse mAb (16B4)	ELISA, IB, IP	a broad range of species	525280	1 Set	€339
Anti-Phosphoserine Mouse mAb (1C8)	ELISA, IB, IP	a broad range of species	525281	1 Set	€436
Anti-Phosphoserine Mouse mAb (4A3)	ELISA, IB, IP	a broad range of species	525283	1 Set	€436
Anti-Phosphoserine Mouse mAb (4A9)	ELISA, IB, IP	a broad range of species	525284	1 Set	€436
Anti-Phosphothreonine Mouse mAb (14B3)	ELISA, IB, IP	a broad range of species	525286	1 Set	€363
Anti-Phosphothreonine Mouse mAb (1E11)	ELISA, IB, IP	a broad range of species	525287	1 Set	€436
Anti-Phosphotyrosine (Ab-1) Mouse mAb (1G2)	FC, FS, PS; NOT IB	a broad range of species	PT01L	200 µg	€107
Anti-Phosphotyrosine (Ab-1) Mouse mAb (1G2) Agarose Conjugate	AC, IP	a broad range of species	PT02	0.25 ml	€415
Anti-Phosphotyrosine Mouse mAb (2C8)	ELISA, IB, IP	a broad range of species	525293	1 Set	€323
Anti-Phosphotyrosine Mouse mAb (PY20)	ELISA, FS, IB, IC, IP	a broad range of species	525295	100 µg	€332
Anti-Phosphotyrosine Mouse mAb (PY20) Agarose Conjugate	AC, IP	a broad range of species	525300	500 µl	€431
Anti-Phosphotyrosine Mouse mAb (PY20) Alkaline Phosphatase Conjugate	ELISA, IB, IH	a broad range of species	525305	100 µg	€339
Anti-Phosphotyrosine Mouse mAb (PY20) Biotin Conjugate	ELISA, IB, IC, IH	a broad range of species	525310	100 µg	€339
Anti-Phosphotyrosine Mouse mAb (PY20) Peroxidase Conjugate	ELISA, IB, IH	a broad range of species	525320	100 µg	€339

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ocherterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Phosphotyrosine Rabbit pAb	ELISA, IB, IP	a broad range of species	525290	100 µg	€321
Anti-PI-3-Kinase Mouse mAb (AB6)	IB, IC, IP	human, mouse, rat	528107	100 µg	€385
Anti-PI-3-Kinase p110δ, C-Terminal (1026-1044) Rabbit pAb	IB, IC	human	526553	100 µl	€343
Anti-PLD/MTA2 (652-668) Rabbit pAb	IB	human, mouse, rat	PC656	100 µg	€338
Anti-Pim1 Rabbit pAb	IB, IP	human	ST1091	50 µg	€148
Anti-Pin1 Rabbit pAb	IB, IF, IP; NOT FS, NOT PS	human	PC270	100 µg	€319
Anti-Pirh2 Rabbit pAb	IB	human	CA1012	50 µg 100 µg	€176 €319
Anti-PKC (α, β, γ) Mouse mAb (M110)	ELISA, IB, IC, IP	human, mouse, rat	PK1014	50 µl	€208
Anti-PKD2 Rabbit pAb	IB, IP	human	ST1042	50 µg	€220
Anti-PKD3 Rabbit pAb	IB	human	ST1090	50 µg	€148
Anti-PKG, C-Ter (657-671) Rabbit pAb	IB, IP	chicken, <i>E. coli</i> , guinea pig, human, mouse, porcine, rabbit, rat, <i>Xenopus</i>	PK10	100 µl	€397
PhosphoDetect™ Anti-PKR (pThr ⁴⁵¹) Rabbit pAb	IB, IH	human, mouse	527460	10 T	€373
Anti-Plasminogen Mouse mAb (10-V-1)	IB	human	528176	500 µg	€352
Anti-Pleiotrophin Goat pAb	ELISA, IB	human	PC187L	100 µg	€259
Anti-Plk1 (8-21) Rabbit pAb	IB, PS; NOT IF	human	PC382	100 µg	€319
Anti-PLZF Mouse mAb (2A9)	GS, IB, IF, IP	human, mouse, rat	OP128	10 µg 100 µg	€86 €338
Anti-PLZF Mouse mAb (2A9)	GS, IB, IF, IP	human, mouse, rat	OP128L	100 µg	€338
Anti-PMS2 Mouse mAb (9)	IB, IP	human, mouse	NA30	100 µg	€338
Anti-Polyoma Virus T Antigen (Ab-4) Mouse mAb (F4)	IF, IP; NOT IB	polyoma infected cells	DP10L	100 µg	€338
Anti-Polypyrimidine Tract Binding Protein (Ab-1) Mouse mAb (SH54)	IB, IF, PS	human; not mouse, not rat	NA63	100 µg	€338
Anti-Porin (Ab-5) (185-197) Rabbit pAb	IB	human, rat	PC548	50 µg	€270
Anti-Porin (Ab-5) Rabbit pAb	IB	human, rat	PC548T	5 µg	€77
Anti-Porin 31HL (Ab-1) Mouse mAb (89-173/016)	ELISA, IB, IC	human	529532	100 µg	€245
Anti-Porin 31HL (Ab-2) Mouse mAb (89-173/025)	ELISA, IB, IC	human	529534	100 µg	€244
Anti-Porin 31HL (Ab-3) Mouse mAb (89-173/033)	EIA, IB, IC	human	529536	100 µg	€245
PhosphoDetect™ Anti-PRAS40 (pThr ²⁴⁶) Rabbit pAb	IB	human, mouse	PS1011	10 T	€350
Anti-Presenilin 1 (14-33) (Ab-3) Goat pAb	ELISA, FS, PS	human	PC361	100 µl	€338
Anti-Presenilin 1 (303-316) (Ab-2) Rabbit pAb	IB, IF	mouse, rat	PC267	50 µg	€270
Anti-Presenilin 1 (31-46) (Ab-1) Rabbit pAb	IB, IF	mouse, rat	PC244	50 µg	€270
Anti-Presenilin 1, Loop Domain (263-407) Rabbit pAb	IB, IF, IP	human, monkey, mouse, rat	529592	100 µl	€397
Anti-Presenilin 1, N-Terminal (1-65) Rabbit pAb	IB, IF, IP	human, monkey, mouse, rat	529591	100 µl	€397
Anti-Presenilin 2 (324-335) (Ab-2) Rabbit pAb	FFS, FS, IB	human, mouse, rat	PC235	25 µg	€338
Anti-Presenilin 2 (7-24) (Ab-1) Rabbit pAb	FS, IB	human, mouse, rat	PC305	25 µg	€338
Anti-Presenilin 2, Loop Domain (269-394) Rabbit pAb	IB, IF, IP	human, mouse	529594	100 µl	€397
Anti-Prion Protein (PrP) (Ab-3) Mouse mAb (F89/160.1.5)	ELISA, FS, IB, PS	bovine, deer, elk, sheep	NB36	200 µg	€338
Anti-Prion PrP27-30, N-Terminal (79-97) Goat pAb	ELISA, PS	bovine, <i>Drosophila</i> , hamster, human, mouse, sheep	530000	100 µl	€294
Anti-Proepithelin Rabbit pAb	IB, IP	human, mouse	IM1000	100 µg	€384
Anti-Progesterone Receptor (Ab-1) Mouse mAb (PR4-12)	IB, PS	chicken, feline, human, rat	GR18	100 µg	€301
PhosphoDetect™ Anti-Progesterone Receptor (pSer ¹⁹⁰) (Ab-1) Mouse mAb (1154)	FS, IB, IP, PS	human	GR47	100 µl	€420
PhosphoDetect™ Anti-Progesterone Receptor (pSer ²⁹⁴) (Ab-2) Mouse mAb (608)	IB	human	GR48	100 µl	€426
Anti-Prohibitin Mouse mAb (II-14-10)	FS, IB, IF, PS	chicken, human, mouse, porcine, rat	CP34	100 µg	€338


10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: outherlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Prostate Apoptosis Response-4 (Ab-1) (317-332) Rabbit pAb	IB, PS	human, mouse, rat	PC426	25 µg	€338
Anti-Prostate Specific Antigen (Ab-1) Mouse mAb (EP-PR8)	FS, PS	human	OP158	100 µl	€253
Anti-20S Proteasome α1-Subunit (249-263) Rabbit pAb	IB	canine, hamster, human, mouse, rat	539145	100 µg	€290
Anti-19S Proteasome S2-Subunit Rabbit pAb	IB, IP	human	539166	100 µl	€290
Anti-19S Proteasome, S6-Subunit Rabbit pAb	IB, IP	human	539169	100 µl	€290
Anti-20S Proteasome α1, 2, 3, 5, 6, & 7-Subunits Mouse mAb (MCP231)	ELISA, IB; NOT IP	human, mouse, rabbit, rat, yeast	ST1049	100 µl	€307
Anti-20S Proteasome α3-Subunit Mouse mAb (MCP257)	IB, PS; NOT IP	human, mouse, rabbit, rat	ST1050	100 µl	€307
Anti-20S Proteasome α5-Subunit Mouse mAb (MCP196)	IB, PS; NOT IP	human, mouse, rabbit, rat	ST1051	100 µl	€307
Anti-20S Proteasome α7-Subunit Mouse mAb (MCP72)	IB, PS; NOT IP	human, rabbit, rat, yeast	ST1052	100 µl	€307
Anti-20S Proteasome β1-Subunit Mouse mAb (MCP421)	IB; NOT IP	human, rabbit	ST1054	100 µl	€307
Anti-20S Proteasome β3-Subunit Mouse mAb (MCP102)	IB; NOT IP	human, mouse, rabbit, rat	ST1055	100 µl	€307
Anti-20S Proteasome β4-Subunit (72-85) Rabbit pAb	IB	human, mouse, rat	ST1056	100 µl	€307
Anti-20S Proteasome β5i-Subunit (23-223) Rabbit pAb	IB, PS	human, mouse, rabbit	ST1057	100 µl	€307
Anti-20S Proteasome Core Subunits Rabbit pAb	IB, IP, PS	human, mouse, yeast	ST1053	100 µl	€307
Anti-26S Proteasome S4-Subunit Rabbit pAb	IB, IP	human, yeast	539167	100 µl	€290
Anti-26S Proteasome S5A-Subunit Rabbit pAb	IB, IP	human	539168	200 µl	€290
Anti-26S Proteasome, S7-Subunit Rabbit pAb	IB, IP	human	539171	100 µl	€290
Anti-Proteasome Activator 11S (PA28γ) Subunit (1-14) Rabbit pAb	IB	human	539148	100 µg	€290
Anti-Proteasome Activator PA28α Subunit (5-22) Rabbit pAb	IB	canine, hamster, human, mouse, rat	539146	100 µg	€290
Anti-Proteasome Activator PA700 Subunit 10B (278-294) Rabbit pAb	IB	canine, hamster, human, mouse, rat	539147	100 µg	€290
Anti-Protein Disulfide Isomerase Rabbit pAb	IB, IP	bovine, canine, guinea pig, hamster, human, monkey, mouse, porcine, rat, sheep, yeast	539229	50 µl	€307
Anti-Protein Kinase A Regulatory Subunit Iα Rabbit pAb	FS, IB	human, monkey, mouse, rat	539232	100 µg	€362
Anti-Protein Kinase A Regulatory Subunit Iβ Rabbit pAb	FS, IB	human, monkey, mouse, rat	539233	100 µg	€362
Anti-Protein Kinase A Regulatory Subunit IIα Rabbit pAb	IB, IH	human, monkey, mouse, rat	539234	100 µg	€362
Anti-Protein Kinase A Regulatory Subunit IIβ Rabbit pAb	FS, IB	human, monkey, mouse, rat	539235	100 µg	€362
PhosphoDetect™ Anti-Protein Kinase A Substrate Mouse mAb (15B12)	ELISA	a variety of species	539236	100 µg	€350
Anti-Protein Kinase C (Ab-1) Rabbit pAb	IB, IF	chicken, human, mouse	PC20	100 µg	€315
Anti-Protein Kinase C (Ab-2) Mouse mAb (MC5)	IB, IP, NT, PS	bovine, human, mouse, rat	OP74	100 µg	€294
Anti-Protein Kinase C _{αβ} (371-388) Rabbit pAb	ELISA, IB, IP	bovine, canine, chicken, <i>Drosophila</i> , human, monkey, mouse, ovine, porcine, rabbit, rat, <i>Xenopus</i>	539507	100 µg	€384
Anti-Protein Kinase C _α (662-672) Rabbit pAb	IB	mammalian	539601	500 µl	€347
Anti-Protein Kinase C _{βi} (661-671) Rabbit pAb	ELISA, IB, IF	rat	539523	100 µl	€384
Anti-Protein Kinase C _{βii} (661-671) Rabbit pAb	IB	mammalian	539603	100 µl	€320
Anti-Protein Kinase C _{βiii} (660-673) Rabbit pAb	ELISA, IB, IC, IF	rat	539526	100 µl	€384
Anti-Protein Kinase C _{βiii} (661-673) Rabbit pAb	IB	mammalian	539605	500 µl	€345
Anti-Protein Kinase C _γ (658-676) Rabbit pAb	IB	mammalian	539607	500 µl	€345
Anti-Protein Kinase C _δ (662-673) Rabbit pAb	ELISA, IB	rat	539532	100 µl	€384
PhosphoDetect™ Anti-Protein Kinase C _δ (pTyr ³¹¹) Rabbit pAb	IB	mouse, rat	ST1019	10 T	€373
Anti-Protein Kinase C _ε (721-737) Rabbit pAb	IB	mammalian	539609	50 µl	€324
Anti-Protein Kinase C _ε (728-737) Rabbit pAb	ELISA, IB, IF	rat	539535	100 µl	€384
Anti-Protein Kinase C _γ (306 - 318) Rabbit pAb	ELISA, IB, IF	bovine, human, rabbit, rat	539529	100 µl	€384
Anti-Protein Kinase C _γ (681-697) Rabbit pAb	IB	mammalian	539611	500 µl	€345

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchtterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Protein Kinase C _ε (580-592) Rabbit pAb	ELISA, IB, IF	rat	539538	100 µl	€384
Anti-Protein Kinase C _{Consensus} (Ac543-549) Rabbit pAb	ELISA, IB, IC, IF, IP	bovine, human, mouse, rabbit, rat	539550	100 µl	€388
Anti-Protein Kinase G Rabbit pAb	IB	bovine	370661	50 µl	€409
PhosphoDetect™ Anti-Protein Kinase G Substrate Mouse mAb (7H1)	ELISA	a variety of species	539725	100 µg	€350
Anti-Protein Phosphatase 1α, C-Terminal Rabbit pAb	IB, IP	human, mouse, rabbit, rat	539517	100 µl	€409
Anti-Protein Phosphatase 1β, C-Terminal (318-327) Rabbit pAb	IB, IP	human, rabbit, rat	539537	100 µl	€384
Anti-Protein Phosphatase 1γ1, C-Terminal (311-323) Rabbit pAb	IB, IP	bovine, human, mouse, rat, <i>Xenopus</i>	539543	100 µl	€384
Anti-Protein Phosphatase 2A/A (7-19) (PR65) Rabbit pAb	IB, IC	mammalian	539509	100 µl	€409
Anti-Protein Phosphatase 2A/Bα (14-27) (PR55) Rabbit pAb	IB, IC	mammalian	539521	100 µl	€384
Anti-Protein Phosphatase 2A/Bγ, N-Terminal (53-66) Rabbit pAb	IB, IP	rat	539546	100 µl	€384
Anti-Protein Phosphatase 2A/C (298-309) Rabbit pAb	IB, IC	mammalian	539525	100 µl	€409
Anti-Protein Phosphatase 2Bα (482-494) Rabbit pAb	IB, IC	mammalian	539530	100 µl	€409
Anti-Protein Phosphatase 2Cγ Mouse mAb (7-53)	IB, IF, IP; NOT FS, NOT PS	human; not mouse, not rat	PH04	100 µg	€338
Anti-Protein Phosphatase Inhibitor-2 Rabbit pAb	IB	human	539642	100 µg	€266
Anti-Protein Tyrosine Phosphatase μ (Ab-2) Mouse mAb (SK15)	IB, IP; NOT FS, NOT PS	human; not mouse, not rat	PH07	100 µg	€338
Anti-Protein Tyrosine Phosphatase-Pest Mouse mAb (Ag 25 2-82)	IB	human; not mouse, not rat	PH06	100 µl	€338
Anti-Protein Tyrosine Phosphatase, T Cell (Ab-1) Mouse mAb (CF4-1D)	IB, IC, IP	human; not mouse	PH03L	100 µg	€338
Anti-Prothrombin (Ab-1) Mouse mAb (GMA-095)	ELISA, IB	human	GF52	100 µg	€248
Anti-PSD95 (Ab-2) Mouse mAb (6G6-1C9)	IB, IF	bovine, mouse, rat	529541	100 µg	€397
Anti-PSD-95 (Ab-2) Mouse mAb (7E3-1B8)	IB, IF, IP	rat	CP35	100 µl	€338
Anti-Pso4 Rabbit pAb	IB, IP	human	DR1022	50 µg	€189
Anti-PSTAIR (42-57) Rabbit pAb	ELISA, IB, IP	a broad range of species	541330	100 µg	€384
PhosphoDetect™ Anti-PTEN (pSer ¹⁸⁰) Rabbit pAb	IB, IC, IP, PS	human, mouse, rat	ST1072	50 µl	€199
Anti-PTHLP (Ab-1) Mouse mAb (212-10.7)	FC, FS, IB, IF, IP, PS	human	GF08	100 µg	€338
Anti-PTHLP (Ab-2) (34-53) Rabbit pAb	FC, IB, IF, IP, PS	human	PC09	100 µg	€319
Anti-PTPase 1B (Ab-1) Mouse mAb (FG6-1G)	IB, IF, IP	avian, hamster, human, mouse (weakly), rat	PH01	100 µg	€338
Anti-PTPase 1B (Ab-2) Mouse mAb (AE4-2J)	IB, IF, IP	human; not mouse	PH02	100 µg	€338
Anti-PUMA (Ab-1) (2-16) Rabbit pAb	IB	human	PC686	100 µg	€350
PhosphoDetect™ Anti-Pyk2 (pTyr ⁴⁰²) Mouse mAb (14F6)	ELISA, IB	human, mouse, rat	549902	100 µg	€350
PhosphoDetect™ Anti-Pyk2 (pTyr ⁴⁰²) Rabbit pAb	IB, IC	chicken, human, rat	549900	10 T	€373

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOP: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections




NEW

Anti-Rat Blood-Brain Barrier Mouse mAb (SMI-71)

Cat. No. NE1026 Detection of rat blood-brain barrier protein by staining frozen sections. Sample: Rat brain. Primary antibody: Anti-Rat Blood-Brain Barrier Mouse mAb (SMI-71) (Cat. No. NE1026) (1:1000). Detection: Fluorescence (red) with Hoechst 33342 counterstain.

Product	Applications	Species Reactivity	Cat No.	Size	Price
PhosphoDetect™ Anti-Pyk2 (pTyr ^{579/580}) Rabbit pAb	IB	chicken, human, rat	ST1023	10 T	€373
PhosphoDetect™ Anti-Pyk2 (pTyr ⁵⁸⁰) Rabbit pAb	IB, IC	human, rat	ST1021	10 T	€373
PhosphoDetect™ Anti-Pyk2 (pTyr ⁸⁸¹) Rabbit pAb	IB	human, rat	ST1022	10 T	€373
Anti-Rab3A (1-18) Rabbit pAb	IB, IP	canine, human, mouse, rat	552097	100 µg	€362
Anti-Rab4 Rabbit pAb	ELISA, IB, IF, IP	canine, Chinese hamster, human, mouse, rat	552104	100 µl	€384
Anti-Rab5 Rabbit pAb	ELISA, IB, IH	human	552108	50 µl	€390
Anti-Rab6 Rabbit pAb	ELISA, IB	human	552109	50 µl	€390
Anti-Rab9 Mouse mAb (mab9)	IB, IF, IH, IP	bovine, canine, hamster, human, monkey, rat	552101	100 µl	€338
Anti-Rab18 Rabbit pAb	ELISA, IB	rat	552126	50 µl	€390
Anti-Rad9 Mouse mAb (93A535)	FS, IB, IP; NOT IF, NOT PS	human	NA72	100 µg	€346
Anti-RAD17 (Ab-1) (2-15) Rabbit pAb	IB	human	PC614	100 µg	€378
PhosphoDetect™ Anti-Rad17 (pSer ⁶⁴⁵) Rabbit pAb	IB	human	DR1013	100 µl	€373
Anti-Rad50 Mouse mAb (2C6)	IB, IP	human	552140	50 µg	€301
Anti-Rad51 (Ab-1) Rabbit pAb	FS, IB, IP, PS	human	PC130	100 µl	€319
Anti-Rad51 (Ab-1) Rabbit pAb	FS, IB, IP, PS	human	PC130T	10 µl	€77
Anti-Rad51 (Ab-2) Mouse mAb (51RAD01)	IB, IP, PS	human	NA71	100 µg	€338
PhosphoDetect™ Anti-Raf (pSer ²⁵⁹) Rabbit pAb	IB, IP, PS	human, mouse, rat	ST1076	50 µl	€199
Anti-A-Raf Protein Kinase, C-Terminal (584-597) Rabbit pAb	ELISA, IB, IP	human, mouse, rat	553001	100 µg	€384
PhosphoDetect™ Anti-c-Raf (pSer ⁶²¹) Mouse mAb (6B4)	ELISA, IB	human	553010	100 µg	€350
PhosphoDetect™ Anti-c-Raf (pTyr ³⁴⁰ , Tyr ³⁴¹) Rabbit pAb	IB, IC	human, mouse, rat	553009	10 T	€373
Anti-Raf1 (Ab-2) (380-403) Rabbit pAb	IB	human, mouse	PC91	100 µg	€319
Anti-RAGE (42-59) Goat pAb	PS	human	ST1025	100 µl	€319
Anti-RANK Mouse mAb (9A725)	IB	human	AM56	100 µg	€338
Anti-RANKL Mouse mAb (12A668)	IB	human, mouse	AM57	100 µg	€338
Anti-RAP Mouse mAb (7F1)	FC, IB, IH	human; not mouse, not rabbit	553507	100 µg	€288
Anti-Raptor Rabbit pAb	IB	human, mouse, rat	ST1048	100 µg	€384
Anti-Ras (31-43) Rabbit pAb	ELISA, IB, IH, IP	human	553571	100 µg	€384
Anti-c-H-Ras (Ab-1) Mouse mAb (F235-1.7.1)	FS, IB, IF, IP, PS	human, mouse, rat	OP23	100 µg 200 µg	€291 €378
Anti-c-H-Ras (Ab-1) Mouse mAb (F235-1.7.1)	FS, IB, IF, IP, PS	human, mouse, rat	OP23T	10 µg	€77
Anti-c-K-Ras (Ab-1) Mouse mAb (234-4.2)	FS, IB, IP, PS	human, mouse, rat	OP24	100 µg	€291

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections




NEW

Anti-TNPO1 Rabbit pAb, Cat No. AP1051 Detection of human TNPO1 by immunoblotting: Sample: Whole cell lysate (30 µg) from HeLa cells. Primary antibody: Anti-TNPO1 Rabbit pAb (Cat. No. AP1051) (1:500). Detection: Chemiluminescence

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-c-K-Ras (Ab-1) Mouse mAb (234-4.2)	FS, IB, IP, PS	human, mouse, rat	OP24T	10 µg	€77
Anti-c-N-ras (Ab-1) Mouse mAb (F155-227)	FS, IB, IP, PS	human, mouse, rat	OP25	100 µg	€324
Anti-v-H-Ras (Ab-1) Rat mAb (Y13-259)	FS, IB, IF, IP, NT, PS	human, mouse, rat	OP01	100 µg	€338
Anti-v-H-Ras (Ab-1) Rat mAb (Y13-259)	FS, IB, IF, IP, NT, PS	human, mouse, rat	OP01L	100 µg	€338
Anti-v-H-Ras (Ab-1) Rat mAb (Y13-259) Agarose Conjugate	IP	human, mouse, rat	OP01A	.5 ml	€457
Anti-v-H-Ras (Ab-2) Rat mAb (Y13-238)	IF, IP; NOT IB, NOT NS	human, mouse, rat	OP04	100 µg	€338
Anti-Pan-Ras (Ab-1) Mouse mAb (F132-62)	FS, IF, IP, PS; NOT IB	human, mouse, rat	OP21	100 µg	€338
Anti-Pan-Ras (Ab-2) Mouse mAb (F111-85)	IB, IP	human, mouse, rat	OP22	100 µg	€338
Anti-Pan-Ras (Ab-3) Mouse mAb (RAS 10)	FC, FS, IB, IF, IP, PS	human, mouse, rat	OP40	100 µg	€338
Anti-Pan-Ras (Ab-4) Mouse mAb (RAS 11)	IB, IP	human, mouse, rat	OP41	100 µg	€338
Anti-Pan-Ras Asp ¹² (Ab-1) Rabbit pAb	IB	human, mouse, rat	PC10L	100 µg	€319
Anti-Pan-Ras Val ¹² (Ab-1) Mouse mAb (DWP)	ELISA, IB, IP	human, mouse, rat	OP38	100 µg	€338
Anti-Rat Blood-Brain Barrier Mouse mAb (SMI-71)	ELISA, FS, PS; NOT IB	rat	NE1026	100 µl	€283
Anti-Rb (Ab-1) Mouse mAb (C36)	IP	human, mouse	OP28	100 µg	€248
Anti-Rb (Ab-1) Mouse mAb (C36) Agarose Conjugate	IP	human, mouse	OP28A	0.5 ml	€457
Anti-Rb (Ab-2) Rabbit pAb	IB	human, mouse	PC13	100 µg	€248
Anti-Rb (Ab-5) Mouse mAb (LM95.1)	IB, PS	human	OP66	100 µg	€248
Anti-Rb (Ab-6) Mouse mAb (AF11)	IB; NOT FS, NOT PS	human, mouse, rat, <i>Xenopus</i>	OP77	100 µg	€248
Anti-Rb (Ab-8) Mouse mAb (XZ55)	IB, IP; NOT FS, NOT IF, NOT PS	chicken, human, mouse, <i>Xenopus</i>	OP99L	100 µg	€248
Anti-Rb (Ab-11) Mouse mAb (1F8)	IB, IF, IP, PS	human; not mouse, not rat	OP136	100 µg	€338
PhosphoDetect™ Anti-Rb (pSer ²⁴⁹ , Thr ²⁵²) Rabbit pAb	IB	human, mouse	PC640	50 µg	€384
PhosphoDetect™ Anti-Rb (pSer ⁶¹²) Rabbit pAb	IB	human	CA1006	10 T	€373
PhosphoDetect™ Anti-Rb (pSer ⁷⁸⁰) Rabbit pAb	IB	human	CA1007	10 T	€373
PhosphoDetect™ Anti-Rb (pThr ³⁷³) Rabbit pAb	IB	human	PC504	25 µg	€282
PhosphoDetect™ Anti-Rb (pThr ⁸²⁶) Rabbit pAb	IB	human	PC644	50 µg	€384
Anti-RCAS1 (Ab-1) Mouse mAb (5E4)	ELISA, FC, IB	human	AM75	100 µg	€360
Anti-C-Reactive Protein Rabbit pAb	ELISA, IP, RIA	human	235752	1 ml	€89
Anti-Reelin (164-496) Mouse mAb (G10)	ELISA, FS, IB, IF, IP, PS	rodent	553731	50 µl	€352
Anti-Reelin (40-189) Mouse mAb (142)	ELISA, IB, IF, IP, PS	a broad range of species	553730	100 µl	€338
Anti-19S Regulator ATPase Subunit Rpt1 Mouse mAb (MSS1-104)	IB; NOT IP	human, mouse	ST1058	100 µl	€307
Anti-19S Regulator non-ATPase Subunit Rpn10 Mouse mAb (S5a-18)	IB, IP	human	ST1060	100 µl	€307

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti Tip60 Rabbit pAb, Cat. No. DR1041

Lane 1: Detection of human TIP60 by immunoblotting. Sample: Whole cell lysate (25 µg) from Saos-2 cells. Primary antibody: Anti-TIP60 Rabbit pAb (Cat. No. DR1041) (1:3000). Detection: Chemiluminescence.

Lanes 2 and 3: Detection of TIP60 by immunoprecipitation. Samples: Whole cell lysate (500 µg) from Saos-2 cells (lane2). Antibody for immunoprecipitation: anti-TIP60 Rabbit pAb (Cat. No. DR1041) (1:100). Conditions for immunoblotting: same as lane 1.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-cRel Rabbit pAb	ELISA, GS, IB, IP	human; not mouse	PC139	20 µl 100 µl	€76 €330
Anti-Relaxin 2 Rabbit pAb	ELISA, IH, RIA	human	553850	100 µg	€319
Anti-RENT1 Goat pAb	IB	human	DR1015	100 µg	€384
Anti-Replication Protein A (Ab-1) Mouse mAb (RPA70-9)	FS, IB, IF, IP, PS	human, yeast	NA13	100 µg	€338
Anti-Replication Protein A (Ab-2) Mouse mAb (RPA34-19)	IB, IF, IP, PS	hamster, human	NA18	100 µg	€338
Anti-Replication Protein A (Ab-3) Mouse mAb (RPA34-20)	IB, IF, IP, PS	human, mouse, yeast	NA19L	100 µg	€338
Anti-Resistin (83-91) Rabbit pAb	ELISA, IB	mouse	554101	100 µl	€409
Anti-Retinoic Acid Receptor α (Ab-1) Rabbit pAb	IB	human; not mouse, not rat	PC92L	100 µg	€319
Anti-Retinoic Acid Receptor β Mouse mAb (336)	IB, IH	human	554722	100 µg	€319
PhosphoDetect™ Anti-RhoA (pSer ¹⁸⁰) Rabbit pAb	IB, IP	bovine, human	ST1035	50 µg	€181
Anti-Ribonucleotide Reductase R2 (165-174) Rabbit pAb	IB, PS	human	PC540	100 µg	€294
Anti-RIP (Ab-1) (52-63) Rabbit pAb	PS; NOT IB, NOT IF	human; not mouse	PC533	100 µg	€294
Anti-RIP2 (275-286) Rabbit pAb	IB, PS; NOT IF	human; not mouse	PC534	100 µg	€294
Anti-RKIP Rabbit pAb	IB, IC	human, monkey, mouse, rat	ST1041	50 µl	€189
Anti-Ryanodine Receptor Mouse mAb (C3-33)	FS, IB, IP, PS	amphibian, avian, canine, fish, rabbit (weakly), rat	NR07	100 µg	€319
Anti-Ryanodine Receptor Mouse mAb (RYR.1)	ELISA, FC, IB, IH, IP, RIA	bovine, canine, human, mouse, porcine, rat	559279	50 µg	€390
Anti-S100 Protein (Ab-1) Mouse mAb (4C4.9)	FS, PS	bovine, human	OP159	100 µl	€253
PhosphoDetect™ Anti-SAPK/JNK (pThr ¹⁸³ , Tyr ¹⁸⁵) Rabbit pAb	IB, IH	<i>C. elegans</i> , chicken, <i>Drosophila</i> , human, mouse, rat	559309	10 T	€373
Anti-SAPK/JNK Rabbit pAb	IB, IC	human, mouse, rat	559304	200 µl	€390
Anti-p38/SAPK 2a Mouse mAb (20B11)	ELISA, IB	human, mouse	559322	100 µg	€409
Anti-PL-Scramblase (309-317) Rabbit pAb	IB, PS	human, mouse	PC367	100 µg	€338
Anti-rSec6 Mouse mAb (9H5)	IB, IF	bovine, canine, chicken, hamster, human, monkey, mouse, porcine, rabbit, rat	559240	100 µg	€413
Anti-Secretory Leukocyte Protease Inhibitor Rabbit pAb	IB, IP	human, mouse	IM1001	100 µg	€384
PhosphoDetect™ Anti-SEK1/MKK4 (pThr ²⁶¹) Rabbit pAb	IB	human, mouse, rat	564970	100 µl	€350
Anti-SERCA2 ATPase Mouse mAb (1ID8)	IB, IF, IH	canine, guinea pig, human, porcine, rabbit, rat	564702	100 µl	€384
Anti-Serotonin (5-HT) Rabbit pAb	FS, IF	rat	PC228L	100 µl	€319
Anti-Serotonin (5-HT) Transporter (602-622) Rabbit pAb	FS, IB, IF	rat	PC177L	100 µl	€445
Anti-Serotonin 2A (5-HT _{2A}) Receptor (22-41) Rabbit pAb			PC176	100 µl	€381
Anti-Serotonin 3 (5-HT ₃) Receptor (444-457) Rabbit pAb	FS; NOT IB	mouse, rat	PC347	25 µg	€338
Anti-Serotonin 3 (5-HT ₃) Receptor (444-457) Rabbit pAb	FS; NOT IB	mouse, rat	PC347T	5 µg	€77
Anti-Serotonin 4A (5-HT _{4A}) Receptor (364-380) Rabbit pAb	ELISA, IB, IC, IP, PS	mouse, rat	NE1010	50 µl	€194
Anti-Serotonin 7 (5-HT ₇) Receptor (8-23) Rabbit pAb	FS	rat	PC249L	100 µl	€445
Anti-Serum Amyloid P Component Rabbit pAb	ELISA, IEP, IF, IH, OT, RID	human; not mouse	565191	1 ml	€192
PhosphoDetect™ Anti-Shc (pTyr ^{239/240}) Mouse mAb (1E3)	ELISA, IB	human, mouse	566808	1 Set	€350
Anti-Shc Rabbit pAb	ELISA, IC, IH, IP	human, mouse, rat	ST1033	50 µg	€278
PhosphoDetect™ Anti-Shc/p66 (pSer ³⁶) Mouse mAb (6E10)	ELISA, IB	human, mouse	566807	1 Set	€350
PhosphoDetect™ Anti-SHIP1 (pTyr ¹⁰²⁰) Rabbit pAb	IB	human, mouse	ST1081	50 µl	€199
Anti-SHP-1 Mouse mAb (1SH01)	IB		566806	100 µg	€319
PhosphoDetect™ Anti-SHP-2 (pTyr ⁵⁴²) Rabbit pAb	IB, IP	human, mouse, rat	ST1082	50 µl	€199
Anti-SHP-2 Rabbit pAb	IB, IP, PS	human, mouse, rat	ST1083	50 µl	€189
Anti-Sialyl Lewis a Mouse mAb (KM231)	AB, ELISA, FC, IB, PS, RIA	human	565942	100 µg	€432


10T: 10 tests by immunoblotting; **AB:** adhesion blockade; **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Sialyl Lewis x Mouse mAb (KM93)	AB, FC, PS	human	565953	100 µg	€432
Anti-α2,3-(N)-Sialyltransferase Rabbit pAb	IB	rat	566219	100 µl	€189
Anti-SIRPα1, C-Terminal (487-503) Rabbit pAb	IB	human	566310	100 µg	€319
Anti-SLAP-2 Rabbit pAb	IB, IP	human, mouse	ST1037	50 µg	€194
Anti-SLP-76 Mouse mAb (AS55)	IB, IP	human	567350	50 µg	€338
Anti-Smac/DIABLO (222-237) Rabbit pAb	IB	human, mouse, rat	567365	100 µg	€301
Anti-Smac/DIABLO (Ab-1) (199-212; 226-239) Rabbit pAb	IB	human	PC574	50 µg	€248
Anti-Smac/DIABLO (Ab-3) (188-198) Rabbit pAb	PS; NOT IB	human	PC605	100 µl	€338
Anti-Smac/DIABLO Rat mAb (10G7)	IB, IC, IP	human, mouse	567360	100 µg	€482
Anti-SMAD1 Rabbit pAb	IB	human	566416	100 µg	€264
PhosphoDetect™ Anti-SMAD1 (pSer ^{463/465}) Rabbit pAb	IB	human, mouse	566411	50 µg	€307
PhosphoDetect™ Anti-SMAD2 (pSer ^{465/467}) Rabbit pAb	IB	human, mink, mouse	566415	50 µl	€262
Anti-SMAD3 (206-220) Rabbit pAb	IB	human, mouse, rat	566414	50 µg	€231
Anti-SMC1 Rabbit pAb	IB, IF, IP	human	PC737	100 µl	€384
Anti-SMC3 Rabbit pAb	IB, IF, IP	human	PC740	100 µl	€384
Anti-SMC5 Rabbit pAb	IB, IP	human	DR1030	50 µg	€189
Anti-SMC6 Rabbit pAb	IB, IP	human	DR1031	50 µg	€189
Anti-SNAP-23A (2-19) Rabbit pAb	IB	bovine, canine, hamster, human	487913	100 µg	€319
Anti-SNAP-25 Rabbit pAb	IB, IH, IP	crayfish, mouse, rat, <i>Torpedo</i> sp. (Electric Ray)	487912	1 Set	€409
Anti-Sp1 (520-534) Rabbit pAb	ELISA, IB	human	PC701	100 µl	€445
Anti-β1-Spectrin Mouse mAb (4C3)	IB, IF	human, mouse, rat	CP33	100 µl	€315
Anti-Sphingosine Kinase 1 (371-384) Goat pAb	IB	human, mouse	PC727	100 µg	€397
Anti-Spi2A (278-294; 406-423) Rabbit pAb	IB	mouse	AP1022	50 µg	€220
Anti-Sprouty (1-207) Rabbit pAb	IB, IP	human	567890	50 µg	€231
PhosphoDetect™ Anti-Src (pTyr ⁴¹⁸) Mouse mAb (9A6)	ELISA, IB, IC	human, mouse, rat	569372	100 µg	€350
PhosphoDetect™ Anti-Src (pTyr ⁴¹⁸) Rabbit pAb	FS, IB, PS	chicken, human, mouse	569373	10 T	€408
PhosphoDetect™ Anti-Src (pTyr ⁶²⁹) Rabbit pAb	IB, IH	chicken, human, mouse, rat	569374	10 T	€408
Anti-v-Src (Ab-1) Mouse mAb (327)	IB, IP; NOT FS, NOT IF, NOT PS	human, mouse, rat	OP07	100 µg	€338
Anti-v-Src (Ab-1) Mouse mAb (327)	IB, IP; NOT FS, NOT IF, NOT PS	human, mouse, rat	OP07L	100 µg	€338
Anti-v-Src (Ab-1) Mouse mAb (327)	IB, IP; NOT FS, NOT IF, NOT PS	human, mouse, rat	OP07T	10 µg	€77
Anti-v-Src (Ab-1) Mouse mAb (327) Agarose Conjugate	IB, IP	human, mouse, rat	OP07A	0.5 ml	€413
Anti-SRC-3 (590-851) Rabbit pAb	IB, IF, IP, PS	human, mouse	567822	100 µl	€409
Anti-STAM1 Rabbit pAb	IB, IP	human, mouse	ST1040	50 µg	€164
Anti-STAM2 Rabbit pAb	IB, IP	human, mouse	ST1038	50 µg	€164
Anti-STAR Rabbit pAb	IB	rat	569420	100 µg	€384
PhosphoDetect™ Anti-STAT1 (pSer ⁷²⁷) Rabbit pAb	IB	human, mouse, rat	569383	10 T	€373
PhosphoDetect™ Anti-STAT1 (pTyr ⁷⁰¹) Rabbit pAb	IB	human	569382	10 T	€373
Anti-STAT1 Rabbit pAb	FC, IB, IP	human, mouse, rat	569386	10 T	€390
PhosphoDetect™ Anti-STAT3 (pSer ⁷²⁷) Rabbit pAb	IB	human, mouse, rat	569392	10 T	€373
PhosphoDetect™ Anti-STAT3 (pTyr ⁷⁰⁵) Mouse mAb (9E12)	ELISA, IB, IP	human, mouse	569384	1 Set	€350
Anti-STAT3 Rabbit pAb	IB, IC, IP	human, mouse, rat	569388	100 µl	€340
PhosphoDetect™ Anti-STAT5 (pTyr ⁶⁹⁴) Rabbit pAb	IB, IC	human, mouse	575142	100 µl	€350
PhosphoDetect™ Anti-STAT6 (pTyr ⁶⁴¹) Rabbit pAb	FC, IB, IF	human, mouse	575143	100 µl	€350
Anti-Stathmin Rabbit pAb	IB, IH, IP	mammalian	569391	100 µl	€445
Anti-Substance P Rabbit pAb	FS, IF	rat	PC232L	100 µl	€259

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-SUMO-1/Sentrin (78-89) Rabbit pAb	PS; NOT IB, NOT IF	human	PC603	100 µl	€294
Anti-Superoxide Dismutase (Cu/Zn Enzyme) Sheep pAb	FFS, FS, IB, IEP, PS	human, monkey, mouse, rat	574597	1 ml	€192
Anti-Superoxide Dismutase (Mn Enzyme) Sheep pAb	FFS, FS, IB, PS	hamster (weakly), human, mouse, rabbit, rat	574596	1 ml	€192
Anti-Survivin (Ab-1) (1-12) Rabbit pAb	IB	human	PC527	100 µg	€338
Anti-Survivin Rabbit pAb	ELISA, IB, IF, IH	human, mouse	574660	50 µg	€319
Anti-SV40 T Antigen (Ab-1) Mouse mAb (PAb419)	IB, IF, IP	SV40 infected cells	DP01	100 µg	€338
Anti-SV40 T Antigen (Ab-1) Mouse mAb (PAb419)	IB, IF, IP	SV40 infected cells	DP01L	100 µg	€338
Anti-SV40 T Antigen (Ab-2) Mouse mAb (PAb416)	IB, IF, IP	SV40 infected cells	DP02	200 µg	€338
Anti-SV40 T Antigen (Ab-2) Mouse mAb (PAb416) Agarose Conjugate	AFF PUR, IP	SV40 infected cells	DP02A	0.5 ml	€445
Anti-SV40 T Antigen (Ab-3) Mouse mAb (PAb280)	FS, IB, IP, PS	SV40 infected cells	DP14L	100 µg	€338
Anti-Synapse-Associated Protein 102 Mouse mAb (7D3)	FS, IB, IF	mouse, rat	573902	100 µg	€384
PhosphoDetect™ Anti-Synapsin I (pSer ⁶⁰³) Rabbit pAb	DB, IB	rat	574774	100 µl	€373
PhosphoDetect™ Anti-Synapsin I (pSer ⁹) Rabbit pAb	DB, IB	rat	574773	100 µl	€444
Anti-Synapsin I Rabbit pAb	ELISA, FS, IB, IC, IP, PS	bovine, human, mouse, rat	574777	10 µg	€278
Anti-Synapsin I Rabbit pAb	ELISA, FS, IB, IC, IP, PS	bovine, human, mouse, rat	574778	50 µl	€141
Anti-Synapsin II (52-65) Rabbit pAb	IB	mouse, rat	574779	100 µl	€384
Anti-Synaptobrevin (Ab-2) Mouse mAb (SP11)	ELISA, IB, PS	human, porcine; not rat	NB07	200 µg	€319
Anti-Synaptophysin (Ab-1) Mouse mAb (SP15)	ELISA, IB, PS	hamster, human, porcine, rat	NB06	200 µg	€319
Anti-Synaptophysin Mouse mAb (SP15)	FS, IB, IC, IP	bovine, feline, hamster, human, mouse, porcine, rat	573822	100 µl	€371
Anti-Synaptotagmin (Ab-1) Mouse mAb (ASV30)	IB, IH, IP	bovine, chicken, frog, mouse, rabbit, rat, shark	573824	200 µg	€384
PhosphoDetect™ Anti-Synaptotagmin (pSer ³⁰⁹) Rabbit pAb	DB, IB	human, mouse, rat	PS1002	100 µl	€385
PhosphoDetect™ Anti-Synaptotagmin (pThr ²⁰²) Rabbit pAb	DB, IB	human, mouse, rat	PS1006	100 µl	€385
PhosphoDetect™ Anti-Synaptotagmin (pThr ²⁰²) Rabbit pAb	DB, IB	human, mouse, rat	PS1006	100 µl	€385
Anti-Synaptotagmin I (1-23) Rabbit pAb	IB, IC, IP	rat, <i>Xenopus</i>	573828	100 µl	€371
Anti-Syntaxin 1 Rabbit pAb	IB, IC, IH	crayfish, mouse, rat, <i>Torpedo</i> sp. (Electric Ray)	574784	200 µg	€413
Anti-Syntaxin 2 (1-19) Rabbit pAb	ELISA, IB, IP	bovine, canine, chicken, hamster, human, monkey, mouse, porcine, rabbit, rat, sheep	574786	100 µg	€384

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-TSC1 (Tuberous Sclerosis 1) Rabbit pAb Cat. No. AP1032


A: Detection of human TSC1 by immunoblotting. Samples: Whole cell lysate (50 µg) from HeLa cells. Primary antibody: Anti-TSC1 (Tuberous Sclerosis 1) Rabbit pAb (Cat. No. AP1032) (0.2 µg/ml). Detection: Chemiluminescence B: Detection of human TSC1 by immunoprecipitation followed by immunoblotting. Sample: Whole cell lysate from HeLa cells. Antibody for immunoprecipitation: Anti-TSC1 (Tuberous Sclerosis 1) Rabbit pAb (Cat. No. AP1032) (1 µg/mg lysate). Immunoblotting conditions same conditions as lane A.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Syntaxin 3 (2-16) Rabbit pAb	IB	human, mouse, rat	574787	1 Set	€231
Anti-Syntaxin 4 (2-23) Rabbit pAb	IB	mouse, rat	574788	1 Set	€231
Anti-Syntaxin Mouse mAb (SP6)	IB	bovine, feline, hamster, human, mouse, porcine, rat	573831	100 µg	€384
Anti-α-Synuclein (Ab-1) (123-140) Guinea pig pAb	FS, IB, IF	bovine, human, mouse, rat	PC325	50 µl	€338
Anti-α-Synuclein, C-Terminal (116-131) Goat pAb	ELISA, PS	human, mouse, rat	575000	100 µl	€313
Anti-TACE (807 - 823) Rabbit pAb	IB	human, mouse, rat	PC491	100 µg	€338
Anti-Tal-1 (7-21) Rabbit pAb	ELISA, IB	human, mouse	PC702	100 µl	€409
Anti-Tankyrase Mouse mAb (19A520)	IB	human	575400	100 µg	€338
PhosphoDetect™ Anti-TAO2 (pSer ¹⁸¹) Rabbit pAb	DB, IB	human, mouse, rat, <i>Xenopus</i>	PS1004	100 µl	€385
Anti-TAP1 (735-748) Rabbit pAb	IB, IC, IP	human	579110	100 µl	€319
Anti-Tau (Ab-1) Mouse mAb (TAU-2)	FS, IB, IF; NOT IP, NOT PS	bovine, chicken, feline, human, monkey	NB16	100 µg	€338
PhosphoDetect™ Anti-Tau (pSer ¹⁹⁹ , Ser ²⁰²) Rabbit pAb	IB	human	577802	10 T	€373
PhosphoDetect™ Anti-Tau (pSer ²¹⁴) Rabbit pAb	IB	human	577811	10 T	€373
PhosphoDetect™ Anti-Tau (pSer ²⁶²) Rabbit pAb	IB	human, mouse, rat	577814	10 T	€373
PhosphoDetect™ Anti-Tau (pSer ³⁵⁶) Rabbit pAb	IB	human	NE1005	10 T	€373
PhosphoDetect™ Anti-Tau (pSer ³⁹⁶) Rabbit pAb	IB	human	577815	10 T	€373
PhosphoDetect™ Anti-Tau (pSer ⁴²²) Rabbit pAb	IB	human	577817	10 T	€373
PhosphoDetect™ Anti-Tau (pThr ¹⁸¹) Rabbit pAb	IB	human	577804	10 T	€373
PhosphoDetect™ Anti-Tau (pThr ²¹²) Rabbit pAb	ELISA, IB, IH	human	577810	10 T	€373
PhosphoDetect™ Anti-Tau (pThr ²³¹) Rabbit pAb	ELISA, IB, IH	human, mouse, rat	577813	10 T	€373
Anti-Tau Mouse mAb (TAU-5)	FS, IB, IP, PS	human, mouse, rat, sheep	577801	100 µg	€378
Anti-TBK1 Mouse mAb (108A429)	IB	human	AM76	100 µg	€360
Anti-TCP-1α Rat mAb (91a)	IB, IF, IP	<i>C. elegans</i> , canine, <i>Drosophila</i> , guinea pig, hamster, human, monkey, mouse, porcine, rabbit, rat	580558	50 µg	€290
Anti-Telomerase (Ab-1) Rabbit pAb	FS, IB, IP	human, mouse	582000	100 µl	€417
Anti-Telomerase (Ab-2) Rabbit pAb	FS, IB, IP, PS	human, mouse	582005	100 µl	€436
Anti-Telomeric Repeat Binding Factor-1 Rabbit pAb	ELISA, IB	human	581420	50 µg	€384
Anti-hTERT (348-358) Rabbit pAb	PS; NOT IB, NOT IF	human; not mouse	PC563	100 µg	€294
Anti-TFF3 Mouse mAb (15C6)	IB, IP	human	585350	100 µg	€409
Anti-Thrombospondin (Ab-1) Mouse mAb (46.4)	FC, FS, IF, IP, PS; NOT IB	human	BA18	100 µg	€338
Anti-Thrombospondin (Ab-3) Mouse mAb (A6.1)	FC, FS, IB, IF, NT; NOT PS	bovine, canine, horse, human, porcine, rat, sheep	BA24	100 µg	€338
Anti-Thrombospondin Rabbit pAb	ELISA, IEP	human	605230	1 ml	€201
Anti-Thymidine Phosphorylase Mouse mAb (P-GF.44C)	FS, IB, NT, PS	human	GF40	100 µg	€338
Anti-Thymosin β4 Rabbit pAb	ELISA, FS	human	609270	100 µl	€319
Anti-Tiam1 Rabbit pAb	IB	human, mouse	ST1070	50 µg	€189
Anti-Tie1, C-Terminal Rabbit pAb	ELISA, IB, IF, IH	bovine, human, mouse, rat	610200	50 µg	€360
PhosphoDetect™ Anti-Tie2 (Ab-1) (pTyr ^{1102/1108}) Rabbit pAb	ELISA, PS	human	PC449	25 µl	€282
PhosphoDetect™ Anti-Tie2 (pTyr ^{1108/1117}) (Ab-2) Rabbit pAb	ELISA, PS	human	PC450	25 µl	€282
Anti-Tie2 Mouse mAb (Ab33)	FS, IB, IP; NOT PS	human, mouse, porcine, rat	610205	50 µg	€231
Anti-Tie2, N-Terminal Rabbit pAb	ELISA, IB	bovine, human, mouse, rat	610204	50 µg	€360
Anti-TIMP-1 (Ab-1) Mouse mAb (7-6C1)	IB, IC, PS	bovine, human, rat; not mouse	IM32	20 µg 100 µg	€113 €338
Anti-TIMP-1 (Ab-2) Mouse mAb (147-6D11)	IB, PS	human	IM41	100 µg	€338
Anti-TIMP-1 (Ab-4) Mouse mAb (102D1)	FS, IB, PS	human	IM63	100 µg	€338

10T: 10 tests by immunoblotting; **AB**: adhesion blockade **AC**: affinity chromatography; **AFF PUR**: affinity purification; **APOP**: apoptosis studies; **DB**: dot blot; **CIA**: cell inhibition assays; **CT**: competition studies; **ELISA**: enzyme-linked immunosorbent assay; **FC**: flow cytometry; **FS**: frozen sections; **GS**: gel shift; **IB**: immunoblotting; **IC**: immunocytochemistry; **IEP**: immunoelectrophoresis; **IF**: immunofluorescence; **IH**: immunohistochemistry; **IP**: immunoprecipitation; **NT**: neutralization; **OT**: ousterlony; **PS**: paraffin sections; **RIA**: radioimmunoassay; **RID**: radial immunodiffusion; **TCELL**: t cell co-stimulation; **VS**: vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-TIMP-1/TIMP-2 (Ab-3) Mouse mAb (102B1)	IB, IH; NOT PS	human, rat	IM62	100 µg	€338
Anti-TIMP-2 (Ab-1) Mouse mAb (T2-101)	FS, IB, IP; NOT PS	human; not mouse, not rat	IM11L	100 µg	€338
Anti-TIMP-2 (Ab-2) Mouse mAb (67-4H11)	FS, IB, PS	bovine, guinea pig, human, mouse, rabbit, rat	IM56	100 µg	€328
Anti-TIMP-2 (Ab-2) Mouse mAb (67-4H11)	FS, IB, PS	bovine, guinea pig, human, mouse, rabbit, rat	IM56T	20 µg	€113
Anti-TIMP-2 (Ab-4) Mouse mAb (T2-N IC3)	AFF PUR, FS, IB, PS	human	IM82	100 µg	€338
Anti-TIMP-3 Mouse mAb (136-13H4)	FS, IB, PS; NOT IF, NOT IP	baboon, human, rabbit; not bovine, not chicken, not porcine, not rat	IM43L	100 µg	€338
Anti-TIMP-3 Mouse mAb (136-13H4)	FS, IB, PS; NOT IF, NOT IP	baboon, human, rabbit; not bovine, not chicken, not porcine, not rat	IM43T	10 µg 20 µg	€77 €113
Anti-TIP60 Rabbit pAb	IB, IP	human	DR1041	50 µl	€135
Anti-Tissue Factor Mouse mAb (TF9-10H10)	DB, ELISA, FC, FS, IB, IF, PS; NOT NS	human, primate; not rabbit	612161	500 µg	€307
Anti-Tissue Plasminogen Activator (Ab-1) Mouse mAb (GMA-043)	ELISA, IB	human	GF53	100 µg	€248
Anti-Tissue Transglutaminase (676-687) Rabbit pAb	IB	human; not mouse, not rat	PC446	100 µg	€294
Anti-Tn Antigen, Sialylated Mouse mAb (B35.1)	ELISA, IB, IH	human	613802	500 µl	€314
Anti-Toll-Like Receptor 2 (179-204) Goat pAb	ELISA, IC	human	PC674	100 µg	€409
Anti-Toll-Like Receptor 3 Mouse mAb (40C1285)	IB	human	OP182	100 µg	€351
Anti-Toll-Like Receptor 5 (151-181) Goat pAb	ELISA, IC	human	PC675	100 µg	€409
Anti-Toll-Like Receptor 6 Mouse mAb (86B1153)	IB	human	OP183	100 µg	€351
Anti-Toll-Like Receptor 8 Mouse mAb (44C143)	FC, IB	human	OP184	100 µg	€351
Anti-Toll-Like Receptor 9 (Ab-1) Mouse mAb (26C593)	FC, IB	human	OP185	100 µg	€351
Anti-TopBP1 Rabbit pAb	IB, IP	human	PC743	100 µl	€384
Anti-Topoisomerase II (Ab-1) Mouse mAb (SWT3D1)	IB, PS; NOT IP	human; not mouse, not rat	NA14	100 µg	€338
PhosphoDetect™ Anti-mTOR (pSer ²⁴⁴⁹) Rabbit pAb	IB	human	PS1020	50 µg	€183
Anti-Toxin A Goat pAb	IB, RIA		616360	500 µl	€220
Anti-TRAF1 Mouse mAb (3D4)	IB	mouse	AP1002	100 µg	€397
Anti-TRAF2 Rabbit (1-501) pAb	IB	human, mouse	AP1040	50 µg	€176
Anti-TRAF6 (1-305) Rabbit pAb	IB, IP	mouse	AP1004	100 µg	€397
Anti-TRAIL Receptor-1 (DR4) (Ab-2) Rabbit pAb	IB	human	PC398	100 µg	€319
Anti-TRAIL Receptor-2 (DR5) (Ab-1) Rabbit pAb	IB, PS	human	PC392	200 µl	€319
Anti-TRAIL Receptor-2, N-Terminal (68-94) Goat pAb	ELISA, IB, IC, PS	human, mouse	616380	100 µl	€409
Anti-TRAIL Receptor-3, N-Terminal (73-103) Goat pAb	ELISA, IB, IC	human, mouse	616383	100 µl	€409
Anti-TRAIL Receptor-4 (249-263) Rabbit pAb	IB	human	616396	100 µg	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade; **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections



NEW

Anti-WASP Rabbit pAb Cat. No. ST1113

Detection of human WASP by immunoblotting. Sample: Molecular weight markers (lane 1) and whole cell lysate (50 µg) from U937 cells (lane 2). Primary antibody: Anti-WASP Rabbit pAb (Cat. No. ST1113) (0.3 µg/ml). Detection: Chemiluminescence.

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-TRAIL, C-Terminal (261-277) Rabbit pAb	ELISA, IB	human	616378	100 µg	€319
Anti-Transferrin Goat pAb	IEP, IP, RID	human	616423	1 ml	€75
Anti-Transferrin Receptor (Ab-1) Mouse mAb (T56/14)	FC, IF, IP; NOT IB	human	GR08L	100 µg	€338
Anti-Transferrin Receptor (Ab-2) Mouse mAb (42/6)	IF, IP, NT	human, rabbit, rat	GR09L	100 µg	€338
Anti-Transforming Growth Factor-α (Ab-1) Mouse mAb (134A-2B3)	IB	human, rat	GF06	100 µg	€248
Anti-Transforming Growth Factor-α (Ab-2) Mouse mAb (213-4.4)	FS, IC, PS; NOT IB	canine, hamster, human, mouse, rat	GF10	100 µg	€248
Anti-Transforming Growth Factor-α (Ab-3) Mouse mAb (189-2130.1)	NT	human, rat	GF15	100 µg	€248
Anti-Transforming Growth Factor-α (Ab-3) Mouse mAb (189-2130.1)	NT	human, rat	GF15L	100 µg	€248
Anti-Transforming Growth Factor-β1 (Ab-1) Mouse mAb (9016.2)	ELISA, IB, NT	human	GF33L	100 µg	€248
Anti-Transforming Growth Factor-β3 (Ab-1) Mouse mAb (236-5.2)	FS, IB, IP, PS	human, mouse, rat	GF16	200 µg	€248
Anti-Transglutaminase II Rabbit pAb	FS, IB, IP, PS	human, mouse, rat	616465	500 µl	€350
Anti-TRB3 (1-145) Rabbit pAb	IB, IC, IP	human, mouse, rat	ST1032	100 µl	€361
Anti-TRF2 Mouse mAb (4A794)	IB, IC, IP	human	OP129	100 µg	€338
Anti-c-Trk (777-790) Rabbit pAb	IB, IF, IP, PS	human, primate, rat	PC31	100 µg	€283
PhosphoDetect™ Anti-TrkA (pTyr ⁴⁹⁰) Rabbit pAb	IB, IP	human, rat	625725	100 µl	€374
Anti-TrkB (807-822) Rabbit pAb	IB, IF	human, rat	PC86	100 µg	€319
Anti-Trypsin Rabbit pAb	ELISA, IB, IEP	human	650297	1 ml	€384
Anti-Tryptase Rabbit pAb	ELISA, IB	human	650367	500 µg	€407
Anti-Tryptophan Hydroxylase (Ab-1) Mouse mAb (PH8)	FS, IB, PS	human, rat	OP71L	100 µg	€324
Anti-TSC1 (Tuberous Sclerosis 1) Rabbit pAb	IB, IP	human, mouse	AP1032	50 µg	€183
PhosphoDetect™ Anti-Tuberin/TSC2 (pThr ¹⁴⁶²) Rabbit pAb	IB	human, mouse	ST1084	50 µl	€231
Anti-α-Tubulin Mouse mAb (DM1A)	IB, IF	chicken, gerbil, human, mouse, rat	CP06	100 µg	€324
Anti-β-Tubulin Mouse mAb (DM1B)	IB, IF	chicken, gerbil, human, mouse, rat	CP07	200 µg	€338
Anti-Tumor Necrosis Factor Receptor (Ab-1) Mouse mAb (16803.1)	ELISA, IB, NT	human	GR28L	100 µg	€248
Anti-Tumor Necrosis Factor Receptor-1 (29-44) Rabbit pAb	ELISA, IB	canine, human, monkey, rabbit, rat	654216	100 µg	€413
Anti-Tumor Necrosis Factor-α (Ab-1) Mouse mAb (1825.121)	ELISA, IB, NT	human; not mouse	GF31	100 µg	€248
Anti-Tumor Necrosis Factor-α Rabbit pAb	ELISA, FS, IB, NT, PS, RIA	human, mouse (weakly)	654250	1 mg	€413
Anti-Tumor Necrosis Factor-α Rabbit pAb	ELISA, IB, NT	mouse, rat	654300	1 ml	€413
PhosphoDetect™ Anti-Tyk2 (pTyr ^{1054/1055}) Rabbit pAb	IB	human	656000	100 µl	€350
Anti-Tyro 3 Rabbit pAb	IB, IP	human	ST1086	50 µg	€189
Anti-Tyrosinase Mouse mAb (T311)	FS, IB, PS	human	OP160	100 µl	€253
PhosphoDetect™ Anti-Tyrosine Hydroxylase (pSer ¹⁹) Rabbit pAb	IB	rat	NE1000	100 µl	€412
PhosphoDetect™ Anti-Tyrosine Hydroxylase (pSer ³¹) Rabbit pAb	DB, IB, IF, IH	rat	NE1001	100 µl	€412
PhosphoDetect™ Anti-Tyrosine Hydroxylase (pSer ⁴⁰) Rabbit pAb	IB	rat	657016	100 µl	€363
Anti-Tyrosine Hydroxylase Mouse mAb (6D7)	IB, IH	a broad range of species	657010	100 µg	€343
Anti-Tyrosine Hydroxylase Rabbit pAb	DB, IB, IF, IH	a broad range of species	657012	100 µl	€343
Anti-Ubc9 Rabbit pAb	IB; NOT FS, NOT PS	human, mouse; not yeast	PC131	100 µg	€255
Anti-Ubiquilin (2-18) Rabbit pAb	IB	mouse	662052	100 µg	€384
Anti-Ubiquitin (Ab-1) Mouse mAb (6C1)	IB	bovine, human	CC37	100 µl	€396
Anti-Ubiquitin Conjugating Enzyme E2 N (2-19; 131-148) Rabbit pAb	IB	human	662118	100 µg	€384

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-Ubiquitin Mouse mAb (2C5)	IB	bovine, human	662097	100 µg	€457
Anti-Ubiquitin Rabbit pAb	IB, IH, IP	bovine, canine, chicken, <i>Drosophila</i> , <i>E. coli</i> , frog, guinea pig, hamster, human, monkey, mouse, rabbit, rat, sheep, yeast	662099	50 µl	€271
Anti-Ubiquitin-Activating Enzyme E1A Rabbit pAb	IB, IH	human	662102	200 µl	€307
Anti-Ubiquitin-Activating Enzyme E1A/E1B Rabbit pAb	IB	human	662106	200 µl	€307
Anti-Uncoupling Protein-1 (UCP-1) (145-159) Rabbit pAb	IB	human, mouse, rat	662045	100 µl	€377
Anti-Uncoupling Protein-2 (UCP-2) (144-157) Rabbit pAb	ELISA, IB, IC	human, mouse, rat	662047	100 µl	€409
Anti-Uncoupling Protein-3 (UCP-3) (147-159) Rabbit pAb	ELISA, IB, IC	mouse, rat; not human	662048	100 µl	€409
Anti-Universal Nitric Oxide Synthase (Ab-2) Mouse mAb (NOS-3F7-B11-B5)	FS, IB	bovine, mouse, rat	NB01	200 µl	€319
Anti-uPA (Ab-1) Mouse mAb (5)	FS, IB, IP, PS	human	IM13L	100 µg	€338
Anti-uPA (Ab-3) Mouse mAb (16)	IB, IP	human	IM15L	100 µg	€338
Anti-USF-1 (1-14) Rabbit pAb	ELISA, IB	human	PC703	100 µl	€409
Anti-VAMP2 (36-56) Rabbit pAb	IB, IC, IH, IP	bovine, mouse, rat, <i>Xenopus</i>	627724	100 µl	€413
Anti-Vanilloid Receptor Like Protein-1 (744-761) Rabbit pAb	FS, IB, IC	mouse, rat	PC421	10 µg 100 µg	€86 €319
Anti-Vasoactive Intestinal Peptide Receptor 1 Mouse mAb (AS58)	ELISA, IB, IC, IP	human, rat	676386	100 µg	€384
Anti-Vasopressin Rabbit pAb	FS, IF	rat	PC234L	100 µl	€319
PhosphoDetect™ Anti-VASP (pSer ¹⁵⁷) Mouse mAb (5C6)	ELISA, IB, IC, IP	human, mouse	676604	100 µg	€389
PhosphoDetect™ Anti-VASP (pSer ²³⁹) Mouse mAb (16C2)	ELISA, IB, IC, IP	human, mouse	676602	100 µg	€363
Anti-VASP Rabbit pAb	ELISA, IB, IC, IP	canine, chicken, human, marsupial, mouse, porcine, rat	676600	100 µl	€409
Anti-Vav-2, SH2 Domain (578-878) Sheep pAb	IB, IP	human	676652	250 µg	€319
Anti-VEGF (Ab-2) Rabbit pAb	IB, PS	human, mouse	PC37	100 µg	€319
Anti-VEGF (Ab-3) Mouse mAb (14-124)	IB, PS	human	GF25	100 µg	€338
Anti-VEGF (Ab-4) Rabbit pAb	IB, IP, PS	human, mouse, rat	PC315	20 µg 500 µg	€86 €338
Anti-VEGF (Ab-5) Mouse mAb (JH121)	FS, IP, NT, PS; NOT IB	human, rabbit	GF41	100 µg	€338
Anti-VEGF (Ab-6) Mouse mAb (VG76e)	IB	bovine, human, porcine	GF46	50 µg	€338
PhosphoDetect™ Anti-VEGF Receptor 1 (Ab-2) (pTyr ¹²¹³) Rabbit pAb	ELISA, PS	human	PC459	25 µl	€282
Anti-VEGF Receptor 1 Goat pAb	ELISA, IB, PS	human	PC322L	100 µg	€527
Anti-VEGF Receptor 2 (931-997) Rabbit pAb	IB, IH, IP	human, mouse	676488	50 µl	€231
Anti-VEGF Receptor 2 Mouse mAb (CH-11)	IB	human, mouse	676486	50 µg	€301
PhosphoDetect™ Anti-VEGF Receptor 2/3 (pAb-1) Rabbit pAb	IB, PS	human	PC460	25 µl	€282
PhosphoDetect™ Anti-VEGFR2 (pTyr ^{1054/1059}) Rabbit pAb	IB	mouse	PS1013	10 T	€350
PhosphoDetect™ Anti-VEGFR2 (pTyr ¹²¹⁴) Rabbit pAb	IB	human, mouse	PS1012	10 T	€350
Anti-VEGF-Receptor 3 Rabbit pAb	ELISA, IB	mouse, rat	343009	50 µg	€360
Anti-Vesicular Monoamine Transporter (VMAT-2) (496-515) Guinea pig pAb	FS, IF	rat	PC301	50 µl	€324
Anti-VGAT (511-525) Guinea Pig pAb	IH	mouse, rat	676780	100 µl	€362
Anti-Vimentin/LN6 (Ab-1) Mouse mAb (V-9)	FS, IB, IF, PS	chicken, human, porcine, rat; not mouse	IF01	100 µg	€338
Anti-Vinculin Mouse mAb (V284 or VLN01)	IB, IP, PS	human	CP74	100 µg	€338
Anti-Vitamin D Receptor (Ab-1) Rat mAb (9A7γE10.4)	FS, GS, IB, IP, PS	amphibian, chicken, hamster, human, mouse, porcine, rat	GR37	100 µg	€338
Anti-Vitronectin Rabbit pAb	ELISA, IB, OT	human	681125	1 ml	€319
Anti-von Hippel-Lindau (VHL) Mouse mAb (Ig33)	FS, IB, IP; NOT PS	human	OP102	100 µg	€338

10T: 10 tests by immunoblotting; AB: adhesion blockade AC: affinity chromatography; AFF PUR: affinity purification; APOF: apoptosis studies; DB: dot blot; CIA: cell inhibition assays; CT: competition studies; ELISA: enzyme-linked immunosorbent assay; FC: flow cytometry; FS: frozen sections; GS: gel shift; IB: immunoblotting; IC: immunocytochemistry; IEP: immunoelectrophoresis; IF: immunofluorescence; IH: immunohistochemistry; IP: immunoprecipitation; NT: neutralization; OT: ouchterlony; PS: paraffin sections; RIA: radioimmunoassay; RID: radial immunodiffusion; TCELL: t cell co-stimulation; VS: vibratome sections

NEW

Product	Applications	Species Reactivity	Cat No.	Size	Price
Anti-WASL/N-WASP Rabbit pAb	IB	human	AP1049	50 µg	€135
Anti-WASP Rabbit pAb	IB	human	ST1113	50 µg	€199
PhosphoDetect™ Anti-WASP (pSer ^{483/484}) Rabbit pAb	IB	human	ST1087	50 µg	€199
Anti-Wilm's Tumor Protein Mouse mAb (WLM04)	IF, PS	human	OP139	100 µg	€338
Anti-Wilm's Tumor Protein Rabbit pAb	IB, IP, PS	human, mouse, rat	CA1026	50 µl	€150
Anti-WRN (400-450) Rabbit pAb	IB, IP	human	DR1032	50 µg	€189
Anti-WWOX (28-42) Rabbit pAb	IB, IC, IP	human, mouse, rat	AP1008	100 µl	€338
Anti-WWOX (89-107) Rabbit pAb	IB, IC	human, mouse, rat	AP1007	100 µl	€338
PhosphoDetect™ Anti-WWOX (pTyr ³³) Rabbit pAb	IB, PS	human, mouse, rat	AP1009	100 µl	€284
Anti-XIAP (Ab-2) Mouse mAb (117318)	NT	human	AM81	100 µg	€378
Anti-XPA (Ab-1) Mouse mAb (5A2)	IB	human	NA66	50 µg	€338
Anti-XPA (Ab-2) Mouse mAb (2A4)	IP; NOT IB	human	NA69	100 µg	€338
Anti-XPG Mouse mAb (8H7)	IB, IP	human	NA67	50 µg	€338
Anti-XRCC3 Rabbit pAb	IB	human	PC691	100 µl	€409
Anti-ZIP Kinase (279-298) Rabbit pAb	IB	human	691700	100 µg	€319

10T: 10 tests by immunoblotting; **AB:** adhesion blockade **AC:** affinity chromatography; **AFF PUR:** affinity purification; **APOP:** apoptosis studies; **DB:** dot blot; **CIA:** cell inhibition assays; **CT:** competition studies; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **GS:** gel shift; **IB:** immunoblotting; **IC:** immunocytochemistry; **IEP:** immunoelectrophoresis; **IF:** immunofluorescence; **IH:** immunohistochemistry; **IP:** immunoprecipitation; **NT:** neutralization; **OT:** ouchterlony; **PS:** paraffin sections; **RIA:** radioimmunoassay; **RID:** radial immunodiffusion; **TCELL:** t cell co-stimulation; **VS:** vibratome sections

Secondary Antibodies & Conjugates

Anti-Human Antibodies & Conjugates

Product	Cat No.	Size	Price
Goat Anti-Human			
Goat Anti-Human IgA, α -Chain Specific	401125	5 ml	€147
Goat Anti-Human IgA, α -Chain Specific Alkaline Phosphatase Conjugate	401132	1 ml	€206
Goat Anti-Human IgA, α -Chain Specific Fluorescein Conjugate	401134	2 ml	€120
Goat Anti-Human IgA, α -Chain Specific Peroxidase Conjugate	401135	2 ml	€153
Goat Anti-Human IgE, ϵ -Chain Specific	411520	1 ml	€89
Goat Anti-Human IgG, γ -Chain Specific	401441	1 ml	€107
Goat Anti-Human IgG, γ -Chain Specific Alkaline Phosphatase Conjugate	401442	1 ml	€206
Goat Anti-Human IgG, γ -Chain Specific Peroxidase Conjugate	401445	2 ml	€153
Goat Anti-Human IgG, Fc Fragment Specific	401439	5 ml	€89
Goat Anti-Human IgG, H & L Chain F(ab') ₂ Fragment Peroxidase Conjugate	401432	500 μ l	€162
Goat Anti-Human IgG, H & L Chain Specific	401954	1 ml	€89
Goat Anti-Human IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401453	1 ml	€206
Goat Anti-Human IgG, H & L Chain Specific Fluorescein Conjugate	401434	2 ml	€113
Goat Anti-Human IgG, H & L Chain Specific Peroxidase Conjugate	401455	2 ml	€147
Goat Anti-Human IgM, μ -Chain Specific Alkaline Phosphatase Conjugate	401902	1 ml	€206
Goat Anti-Human IgM, μ -Chain Specific Peroxidase Conjugate	401905	2 ml	€153
Mouse Anti-Human			
Mouse Anti-Human IgG, Fc Fragment Specific (HP6043) Peroxidase Conjugate	411550	500 μ g	€252
Mouse Anti-Human IgA, Fab Fragment Specific (HP6123)	411420	400 μ g	€192
Mouse Anti-Human IgA, Secretory (HP6141)	411423	400 μ g	€192
Mouse Anti-Human IgE, Fc Fragment Specific (HP6029)	411507	200 μ l	€120
Mouse Anti-Human IgG ₁ , Fc Fragment Specific (HP6070)	411455	400 μ g	€185
Mouse Anti-Human IgG ₁ , Fc Fragment Specific (HP6069A)	411452	200 μ l	€117
Mouse Anti-Human IgG ₁ , Fc Fragment Specific (HP6069B) Biotin Conjugate	411543	200 μ g	€153
Mouse Anti-Human IgG ₁ , Fc Fragment Specific (HP6069P)	411451	400 μ g	€192
Mouse Anti-Human IgG ₁ , Fc Fragment Specific (HP6070)	411454	200 μ l	€137
Mouse Anti-Human IgG ₁ , Fc Fragment Specific Alkaline Phosphatase Conjugate	401459	250 μ g	€334
Mouse Anti-Human IgG ₂ , Fc Fragment Specific Alkaline Phosphatase Conjugate	401462	250 μ g	€334
Mouse Anti-Human IgG ₂ , Fc Fragment Specific (HP6002) Biotin Conjugate	411462	200 μ g	€153
Mouse Anti-Human IgG ₂ , Fd Fragment Specific (HP6014)	411466	400 μ g	€192
Mouse Anti-Human IgG ₂ , IgG ₃ , IgG ₄ (HP6030)	411427	400 μ g	€192
Mouse Anti-Human IgG ₃ , Fc Fragment Specific Alkaline Phosphatase Conjugate	401465	250 μ g	€334
Mouse Anti-Human IgG ₃ , Hinge Specific (HP6047) Biotin Conjugate	411483	200 μ g	€153

Product	Cat No.	Size	Price
Mouse Anti-Human IgG ₁ Fc Fragment Specific (HP6023B) Biotin Conjugate	411493	200 µg	€153
Mouse Anti-Human IgG ₁ Fc Fragment Specific (HP6025)	411492	400 µg	€100
Mouse Anti-Human IgG ₁ Fc Fragment Specific (HP6025) Alkaline Phosphatase Conjugate	401468	250 µg	€334

Rabbit Anti-Human

Rabbit Anti-Human IgM, µ-Chain Specific	401852	1 ml	€113
---	--------	------	------

Anti-Mouse Antibodies & Conjugates

Immobilized Anti-Mouse IgG	121937	2 ml	€206
Hybridoma Subtyping Kit, Mouse	386445	1 kit	€338

Goat Anti-Mouse

Goat Anti-Mouse IgG	0502	200 µg	€66
Goat Anti-Mouse IgG Alkaline Phosphatase Conjugate	DC05L	200 µg	€84
Goat Anti-Mouse IgG Fluorescein Conjugate	DC13L	100 µg	€66
Goat Anti-Mouse IgG Peroxidase Conjugate	DC02L	200 µg	€84
Goat Anti-Mouse IgG TEXAS RED® Conjugate	DC18L	100 µg	€84
Goat Anti-Mouse IgG, γ-Chain Specific	401223	1 ml	€107
Goat Anti-Mouse IgG, H & L Chain β-Galactosidase Conjugate	401607	1 ml	€206
Goat Anti-Mouse IgG, H & L Chain Specific	401211	1 ml	€107
Goat Anti-Mouse IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401212	1 ml	€206
Goat Anti-Mouse IgG, H & L Chain Specific Alkaline Phosphatase Conjugate, Adsorbed	401252	1 ml	€260
Goat Anti-Mouse IgG, H & L Chain Specific Biotin Conjugate	401213	2 ml	€114
Goat Anti-Mouse IgG, H & L Chain Specific Biotin Conjugate, Adsorbed	401216	2 ml	€176
Goat Anti-Mouse IgG, H & L Chain Specific Fluorescein Conjugate	401214	2 ml	€114
Goat Anti-Mouse IgG, H & L Chain Specific Fluorescein Conjugate, Adsorbed	401234	1 ml	€206
Goat Anti-Mouse IgG, H & L Chain Specific Fluorescein Conjugate	401244	1 ml	€206
Goat Anti-Mouse IgG, H & L Chain Specific Phycoerythrin Conjugate	401239	1 ml	€270
Goat Anti-Mouse IgG, H & L Chain Specific Rhodamine Conjugate	401217	2 ml	€120
Goat Anti-Mouse IgG, H & L Chain Specific Rhodamine Conjugate, Adsorbed	401245	2 ml	€176
Goat Anti-Mouse IgG, H & L Chain Specific TEXAS RED® Conjugate	401230	2 ml	€161
Goat Anti-Mouse IgG, H & L Chain Specific TEXAS RED ^{trp} Conjugate, Adsorbed	401246	2 ml	€192
Goat Anti-Mouse IgG, H & L Chain Specific, F(ab') ₂ Fragment Peroxidase Conjugate	401233	1 ml	€233
Goat Anti-Mouse IgG, H & L Chain Specific Peroxidase Conjugate	401215	2 ml	€159
Goat Anti-Mouse IgM, µ-Chain Specific TEXAS RED® Conjugate	401296	2 ml	€130
Goat Anti-Mouse IgM, H & L Chain Specific Alkaline Phosphatase Conjugate	401295	1 ml	€206
Goat Anti-Mouse IgM, H & L Chain Specific Fluorescein Conjugate	401224	2 ml	€114
Goat Anti-Mouse IgM, H & L Chain Specific Peroxidase Conjugate	401225	2 ml	€159
Goat Anti-Mouse Total Ig Peroxidase Conjugate	DC08L	200 µg	€84
Goat Anti-Mouse IgG, γ-Chain Specific Peroxidase Conjugate	401207	2 ml	€201
Goat Anti-Mouse IgG, H&L Chain Specific Peroxidase Conjugate	401253	2 ml	€196

Rabbit Anti-Mouse

Rabbit Anti-Mouse IgG, H & L Chain Specific	402334	1 ml	€154
Rabbit Anti-Mouse IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401262	1 ml	€206
Rabbit Anti-Mouse IgG, H & L Chain Specific Fluorescein Conjugate	401219	2 ml	€120
Rabbit Anti-Mouse IgG, H & L Chain Specific Peroxidase Conjugate	402335	2 ml	€161

Anti-Rabbit Antibodies & Conjugates

Goat Anti-Rabbit IgG Alkaline Phosphatase Conjugate	DC06L	200 µg	€84
---	-------	--------	-----

Product	Cat No.	Size	Price
Goat Anti-Rabbit IgG Biotin Conjugate	OS03B	200 µg	€74
Goat Anti-Rabbit IgG Fluorescein Conjugate	DC15L	100 µg	€66
Goat Anti-Rabbit IgG Peroxidase Conjugate	DC03L	200 µg	€84
Goat Anti-Rabbit IgG TEXAS RED® Conjugate	DC20L	100 µg	€84
Goat Anti-Rabbit IgG, Fc Fragment Specific Alkaline Phosphatase Conjugate	401371	1 ml	€206
Goat Anti-Rabbit IgG, H & L Chain Specific	401311	1 ml	€100
Goat Anti-Rabbit IgG, H & L Chain Specific, Adsorbed	401962	1 ml	€147
Goat Anti-Rabbit IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401312	1 ml	€206
Goat Anti-Rabbit IgG, H & L Chain Specific Alkaline Phosphatase Conjugate, Adsorbed	401352	1 ml	€229
Goat Anti-Rabbit IgG, H & L Chain Specific Biotin Conjugate	401313	2 ml	€107
Goat Anti-Rabbit IgG, H & L Chain Specific Fluorescein Conjugate	401314	2 ml	€113
Goat Anti-Rabbit IgG, H & L Chain Specific Fluorescein Conjugate, Adsorbed	401319	2 ml	€206
Goat Anti-Rabbit IgG, H & L Chain Specific Peroxidase Conjugate	401315	2 ml	€148
Goat Anti-Rabbit IgG, H & L Chain Specific Peroxidase Conjugate	401353	2 ml	€206
Goat Anti-Rabbit IgG, H & L Chain Specific Peroxidase Conjugate	401393	2 ml	€254
Goat Anti-Rabbit IgG, H & L Chain Specific Phycoerythrin Conjugate	401299	1 ml	€261
Goat Anti-Rabbit IgG, H & L Chain Specific Rhodamine Conjugate	401321	2 ml	€113
Goat Anti-Rabbit IgG, H & L Chain Specific TEXAS RED® Conjugate	401355	2 ml	€120

Anti-Rat Antibodies & Conjugates

Goat Anti-Rat IgG Peroxidase Conjugate	DC01L	200 µg	€84
Goat Anti-Rat IgG, H & L Chain Specific	401411	1 ml	€100
Goat Anti-Rat IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401409	1 ml	€206
Goat Anti-Rat IgG, H & L Chain Specific Fluorescein Conjugate	401414	2 ml	€113
Goat Anti-Rat IgG, H & L Chain Specific Peroxidase Conjugate	401416	2 ml	€159

Other Secondary Antibodies & Conjugates

Anti-Goat Antibodies & Conjugates

Rabbit Anti-Goat IgG, Fc Fragment Specific Alkaline Phosphatase Conjugate	401502	1 ml	€206
Rabbit Anti-Goat IgG, Fc Fragment Specific Peroxidase Conjugate	401504	2 ml	€192
Rabbit Anti-Goat IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	401512	1 ml	€207
Rabbit Anti-Goat IgG, H & L Chain Specific Fluorescein Conjugate	401514	2 ml	€122
Rabbit Anti-Goat IgG, H & L Chain Specific Peroxidase Conjugate	401515	2 ml	€159


Anti-Sheep Antibodies & Conjugates

Rabbit Anti-Sheep IgG, H & L Chain Specific	401966	1 ml	€107
Rabbit Anti-Sheep IgG, H & L Chain Specific Alkaline Phosphatase Conjugate	402102	1 ml	€206
Rabbit Anti-Sheep IgG, H & L Chain Specific Fluorescein Conjugate	402104	2 ml	€120
Rabbit Anti-Sheep IgG, H & L Chain Specific Peroxidase Conjugate	402100	2 ml	€159

Calbiochem®

Diabetes and Obesity

Including unique research reagents and kits



Request your FREE copy today!
www.calbiochem.com/diabetes-obesity

Request
a copy!



Calbiochem®

New Research Tools for Signal Transduction and Life Science Research

Highlighting new research
reagents and kits

Request
a copy!

Call us or visit our website:
www.calbiochem.com/newtools

Antibody Support Products

Product	Application	Cat No.	Size	Price
Adjuvants & Immunization Tools				
N-Acetyl-D-glucosaminyl-(β 1,4)-N-acetylmuramyl-L-alanyl-D-isoglutamine		110150	5 mg	€216
N-Acetylmuramyl-L-alanyl-D-isoglutamine		110125	1 mg	€55
DNP-Hemocyanin Conjugate, Keyhole Limpet		324121	100 mg	€102
Freund's Complete Adjuvant, Modified, <i>Mycobacterium butyricum</i>		344289	1 Set	€126
Freund's Incomplete Adjuvant, Modified		344291	10 ml	€32
Hemocyanin, <i>Concholepas concholepas</i> , Glycerol		374803	1 g	€422
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i>		374805	50 mg 250 mg 1 g	€53 €212 €662
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i> , Ammonium Sulfate		374811	50 mg 250 mg 1 g	€46 €184 €590
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i> , BES/Mg ²⁺ Buffer		374807	250 mg	€204
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i> , Glycerol		374817	50 mg 250 mg 1 g	€53 €202 €606
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i> , High Purity, Endotoxin-Free, Glycerol		374813	5 mg 25 mg	€64 €238
Hemocyanin, Keyhole Limpet, <i>Megathura crenulata</i> , High Purity, Endotoxin-Free, Sterile-Liquid		374825	25 mg	€179
Pam ₃ Cys-Ser-(Lys) ₄ , Hydrochloride		506350	2 mg	€269

Avidin, Streptavidin, Biotin & Conjugates

Avidin, Alkaline Phosphatase Conjugate	EIA, IH	189726	2 ml	€183
Avidin, Egg White		189725	10 mg 100 mg	€64 €390
Avidin, Fluorescein Conjugate	IF, IH	189727	1 ml	€147
Avidin, Peroxidase Conjugate	ELISA, IH	189728	2 ml	€154
Biotin, Fluorescein Conjugate	EIA, IH	343222	1 mg	€114
Biotin, Peroxidase Conjugate		203194	1 ml	€102
Biotin-NHS		203112	25 mg 250 mg 1 g	€34 €231 €589
Biotin-NHS, Water-Soluble		203118	50 mg	€136
Biotin-X-NHS		203188	25 mg 50 mg	€63 €110
Biotin-X-NHS Kit		203187	1 kit	€224
Biotin-X-NHS, Water-Soluble		203189	100 mg	€164
Biotin-XX-Hydrazide		203110	50 mg	€185

ABPUR: antibody purification; **EIA:** enzyme immunoassay; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **IB:** immunoblotting; **IC:** immunocytochemistry; **IH:** Immunohistochemistry; **IP:** immunoprecipitation; **IP:** immunoprecipitation; **NT:** neutralization; **PS:** paraffin sections; **RIA:** radioimmunoassay

Product	Application	Cat No.	Size	Price
Biotin-XX-NHS		203114	25 mg 100 mg	€108 €352
D-(+)-Biotin		2031	1 g	€52
MBB		442631	25 mg	€122
Streptavidin		189730	1 mg 10 mg	€97 €605
Streptavidin, β -Galactosidase Conjugate	EIA, IH	569404	2 ml	€206
Streptavidin, Alkaline Phosphatase Conjugate	ELISA, IB, IH	189732	2 ml	€300
Streptavidin, Alkaline Phosphatase Conjugate	ELISA, IB, IH	OR04L	200 μ g	€74
Streptavidin, Fluorescein Conjugate	IF, IH	189734	1 ml	€147
Streptavidin, Horseradish Peroxidase Conjugate	ELISA, IB, IH	OR03L	200 μ g	€74
Streptavidin, Peroxidase Conjugate	EIA	189733	1 ml	€184
Streptavidin, Phycoerythrin Conjugate	EIA, FC	189737	1 ml	€223
Streptavidin, TEXAS RED® Conjugate	ELISA, FC, IH	189738	1 ml	€154

Extracts, Lysates, Cell Pellets & Blots

HL-60 Cell Extract, Induced	IB	383500	200 μ l	€91
INSTA-Blot™, Human Cell Lines	IB	407210	1 blot	€330
INSTA-Blot™, Human Tissues	IB	407211	1 blot	€330
INSTA-Blot™, Rat Tissues	IB	407213	1 blot	€330
p21 H-Ras (Gly ¹²) Western Blot Standard	IB	WA01	2 μ g	€133
K-ras Val ¹² Western Blot Standard	IB	WA04	2 μ g	€133
A431 Cell Pellet		WB36	10 lane	€133
HeLa Cell Extract, Normal		WB29	200 μ g	€105
HeLa Cell Pellet		WB37	10 lane	€133
Jurkat Cell Pellet		WB43	10 lane	€133
MCF7 Cell Pellet		WB39	10 lane	€133
SK-BR-3 Cell Pellet		WB38	10 lane	€133
HeLa Cell Nuclear Extract		WB64	200 μ g	€175

Markers

Perfect Protein™ Markers, 15–150 kDa	69149	100 lanes	€110
Perfect Protein™ Markers, 10–225 kDa	69079	100 lanes	€110
Strep•Tag® II Perfect Protein™ Markers, 16–100 kDa, Lyophilized	71614	100 lanes	€162
Perfect Protein™ Western Markers	69959	25 lanes	€59
Perfect Protein™ AP Western Blot Kit	69965	25 blots	€158
Perfect Protein™ HRP Western Blot Kit	69078	25 blots	€158
Trail Mix™ Protein Markers	70980	100 lanes	€135
Trail Mix™ Western Markers	70982	25 lanes	€69
Trail Mix™ AP Western Blot Kit	71047	25 blots	€169
Trail Mix™ HRP Western Blot Kit	71048	25 blots	€169
4X SDS Sample Buffer	70607	2 m	€29

Immunodetection Reagents

5% Alkali-soluble Casein	70955	225 ml	€66
AntigenPlus™ Buffer Set	71292	1 Set	€128
5X AntigenPlus™ Buffer (pH 6)	71289	100 ml	€53
5X AntigenPlus™ Buffer (pH 7.4)	71291	100 ml	€53

ABPUR: antibody purification; **EIA:** enzyme immunoassay; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **IB:** immunoblotting; **IC:** immunocytochemistry; **IH:** Immunohistochemistry; **IP:** immunoprecipitation; **IP:** immunoprecipitation; **NT:** neutralization; **PS:** paraffin sections; **RIA:** radioimmunoassay

Product	Application	Cat No.	Size	Price
5X AntigenPlus™ Buffer (pH 10)		71290	100 ml	€53
Biotin/Streptavidin Immunostaining Kit, Anti-Mouse IgG, Peroxidase		568108	1 kit	€308
Biotin/Streptavidin Immunostaining Kit, Anti-Rabbit IgG, Peroxidase		568112	1 kit	€308
BLOT-FastStain™ Reagent	IB	WB56	1 each	€292
BLOT-QuickBlocker™ Reagent	IB	WB57	175 g	€69
BrdU Immunohistochemistry System		HCS30	1 each	€515
ChromoSave™ Reagent		238540	30 ml	€94
Fluorescent Mounting Media, Aqueous		HC08	8 ml	€69
FluorSave™ Reagent		345789	20 ml	€77
ImmunoPen™		402176	1 each	€78
MOWIOL® 4-88 Reagent		475904	100 g	€102
Mouse UniTect™ ABC Kit	IC, IH, PS	XHC01	1 each	€170
Rabbit UniTect™ ABC Kit	IC, IH, PS	XHC02	1 each	€182
Western-Re-Probe™ Reagent	IB	WB59	100 ml	€167

Protein A & Protein G Reagents

PANSORBIN® Cells	IP	507858	1 g 5 g 10 g 25 g 50 g	€77 €306 €486 €965 €1700
PANSORBIN® Cells, Lyophilized	EIA, IP, RIA	507862	1 g 5 g	€76 €291
PANSORBIN® Cells, Standardized	IP	507861	25 ml 50 ml	€159 €279
Protein A, <i>Staphylococcus aureus</i>	AB PUR, IH	539202	5 mg	€155
Protein A, Alkaline Phosphatase Conjugate	ELISA, IB	539251	1 mg	€178
Protein A, Biotin Conjugate		203195	2 mg	€107
Protein A, Peroxidase Conjugate	ELISA, IB	539253	1 mg	€120
Protein G, <i>Streptococcus</i> sp.		539302	1 mg	€163
Protein G, Alkaline Phosphatase Conjugate		539305	500 µg	€137
Protein G, Fluorescein Conjugate	IC	539326	1 mg	€189
Protein G, Peroxidase Conjugate	ELISA, IB	539322	500 µg	€137
Protein G, Recombinant, <i>E. coli</i>	AB PUR	539303	1 mg 5 mg	€94 €320
TACHISORB Immunoabsorbent, Rabbit	RIA	575544	55 ml	€104

Substrates for AP & HRP

AEC	IB, IH	152224	10 ml	€89
4-Chloro-1-naphthol		220480	5 g	€82
ABTS™ Chromophore	ELISA	194434	100 ml	€45
ABTS™ Chromophore, Diammonium Salt	ELISA	194430	1 g	€39
ABTS™ Chromophore, Diammonium Salt		194437	1 each	€53
AEC Substrate Buffer	IB, IH	152226	50 ml	€53
AP Detection Reagent Kit		69264	1 ea 5 ea	€89 €370
5-Bromo-4-chloro-3-indolyl- α -D-galactopyranoside		203780	25 mg	€106
5-Bromo-4-chloro-3-indolyl- α -D-glucopyranoside		203781	25 mg	€102

ABPUR: antibody purification; **EIA:** enzyme immunoassay; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **IB:** immunoblotting; **IC:** immunocytochemistry; **IH:** Immunohistochemistry; **IP:** immunoprecipitation; **IP:** immunoprecipitation; **NT:** neutralization; **PS:** paraffin sections; **RIA:** radioimmunoassay

Product	Application	Cat No.	Size	Price
5-Bromo-4-chloro-3-indolylphosphate, Disodium Salt		203788	100 mg 1 g	€41 €235
BCIP/NBT		203790	100 ml	€52
BCIP/TNBT		203795	100 ml	€78
CDP-Star® AP Substrate		69086	40 ml	€164
Chlorophenolred-β-D-galactopyranoside, Sodium Salt		220588	250 mg	€85
DAB Substrate Buffer		281753	50 ml	€53
DAB, Tetrahydrochloride		281751	10 ml	€89
o-Phenylenediamine, Dihydrochloride, 10-mg Tablets		523121	1 each	€113
o-Phenylenediamine, Dihydrochloride, 2-mg Tablets		523120	1 each	€102
Phenylphosphate, Disodium Salt		516798	25 g	€39
Streptavidin AP LumiBlot™ Kit		70554	25 blots	€304
SuperSignal® HRP Substrate		69059	2 x 25 ml	€82
TMB, Insoluble		613548	100 ml	€71
TMB, Soluble		613544	100 ml	€53
TMB Solution		CL07	100 ml 1000 ml	€69 €455

ABPUR: antibody purification; **EIA:** enzyme immunoassay; **ELISA:** enzyme-linked immunosorbent assay; **FC:** flow cytometry; **FS:** frozen sections; **IB:** immunoblotting; **IC:** immunocytochemistry; **IH:** Immunohistochemistry; **IP:** immunoprecipitation; **IP:** immunoprecipitation; **NT:** neutralization; **PS:** paraffin sections; **RIA:** radioimmunoassay

Protein A/G Affinity Resins

Protein A/G Affinities for IgG from Various Species

Species	Affinity for Protein A	Affinity for Protein G
Human	++++	++++
Equine	++	++++
Bovine	++	++++
Porcine	+++	+++
Sheep	+/-	++
Goat	-	++
Rabbit	++++	+++
Chicken	-	+
Hamster	+	++
Guinea Pig	++++	++
Rat	+/-	++
Mouse	++	++

Protein A/G Affinities for Various Antibody Subclasses

Species	Affinity for Protein A	Affinity for Protein G
Human IgG ₁	++++	++++
Human IgG ₂	++++	++++
Human IgG ₃	-	++++
Human IgG ₄	++++	++++
Rat IgG ₁	-	+
Rat IgG _{2a}	-	++++
Rat IgG _{2b}	-	++
Rat IgG _{2c}	+	++
Mouse IgG ₁	+	++++
Mouse IgG _{2a}	++++	++++
Mouse IgG _{2b}	+++	+++
Mouse IgG ₃	++	+++

Prices and availability are subject to change. ©Copyright 2007 EMD Chemicals, an affiliate of Merck KGaA, Darmstadt, Germany. All rights reserved. Each product is sold with a limited warranty which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. EMD Chemicals products may not be resold, modified for resale, or used to manufacture commercial products without written approval of EMD Chemicals. CALBIOCHEM® and PANSORBIN® are registered trademarks of EMD Chemicals in the United States and in certain other jurisdictions. AntigenPlus™, ChromoSave™, FluorSave™, ImmunoPen™, Interactive Pathways™, LumiBlot™, Perfect Protein™, PhosphoDetect™, SourceBook™, Trail Mix™, and UniTect™ are trademarks of EMD Chemicals. ABTS™ is a trademark of Roche Diagnostics, BLOT-FastStain™, BLOT-QuickBlocker™, and Western-Re-Probe™ are trademarks of Geno Technology, Inc., CDP-Star® is a registered trademark of Applera Corporation INSTA-Blot™ is a trademark of Imgenex, MOWIOL® is a registered trademark of Clariant GmbH, Strep-Tag® is a registered trademark of IBA GmbH, SuperSignal® is a registered trademark of Pierce Biotechnology, Inc., TEXAS RED® is a registered trademark of Molecular Probes, Inc.

Thank you
for purchasing
Merck products!



Calbiochem®

Accelerate Your Cell Signaling
Research with our
Interactive Pathways™ Resource

More than **20** Interactive Pathways
to choose from, backed by 30 years of
experience in signal transduction.

Over 10,000 products
to choose from

Select a Pathway, choose your
target, and click to find related:

- antibodies
- assay kits
- enzymes and substrates
- inhibitors
- proteins




Visit calbiochem.com/pathways
to view pathways for:

- Akt/PKB
- Alzheimer's Disease
- Apoptosis
- GSK-3
- Insulin/IGF
- NEW** • mTOR
- MAPK Family
- Nitric Oxide
- PKC
- Tyrosine Kinases
- Ubiquitin-Proteasome and more

To order a poster please contact
your local Merck office or go to
www.merckbio.eu/literature

For more information or to place an order,
contact your local distributor (see back cover).

INTERACTIVE PATHWAYS



Austria

VWR International GmbH
Tel: +43 1 97 002 0
Fax: +43 1 97 002 600
E-mail: info@at.vwr.com
www.vwr.com

Belgium

VWR International bvba/sprl
Tel: +32 16 385 011
Fax: +32 16 385 385
E-mail: customerservice@be.vwr.com
www.vwr.com

Bulgaria

AQUACHIM
Merck KGaA Representative
Tel: +359 2 98 33 760
Fax: +359 2 97 12 092
E-mail: merck-aquachim@bgnet.bg

Croatia

Merck d.o.o.
Tel: +385 1 4864 111
Fax: +385 1 4864 191
E-mail: kemija@merck.t-com.hr

Czech Republic

Merck Spol s.r.o.
Tel: +420 323 619 323
Fax: +420 323 619 303
E-mail: infochemie@merck.cz
www.merck.cz

Denmark

VWR International ApS
Tel: +45 43 86 87 88
Fax: +45 43 86 87 90
E-mail: info@dk.vwr.com
www.vwr.com

Finland

VWR International Oy
Tel: +358 9 804 551
Fax: +358 9 804 55200
E-mail: info@fi.vwr.com
www.vwr.com

France

VWR International S.A.S
Tel: 0 825 02 30 30 (0.15€TTC/min)
Fax: 0 825 02 30 35 (0.15€TTC/min)
E-mail: info@fr.vwr.com
www.vwr.com

Greece

Merck E.P.E
Tel: +302 10 9885 300, 349, 350
Fax: +302 10 9885 400
E-mail: merckhel@otenet.gr

Hungary

Merck Kft.
Tel: +36 1 463 8100
Fax: +36 1 463 8101
E-mail: kemia@merck.hu
www.merck.hu

Italy

VWR International s.r.l.
Tel: +39 02 3320311 centralino
Tel: +39 02 332031323 hot line
Fax: 800 152999
E-mail: biotech@it.vwr.com
www.vwr.com

Macedonia

VARUS
Merck KGaA Representative
Tel: +389 2 3115 122
Fax: +389 2 3137 013
E-mail: varus@mt.net.mk

Netherlands

VWR International B.V.
Tel: +31 20 4808 400
Fax: +31 20 4808 480
E-mail: info@nl.vwr.com
www.vwr.com

Norway

VWR International AS
Tel: +47 02290
Fax: +47 22 90 00 40
E-mail: info@no.vwr.com
www.vwr.com

Poland

Merck Sp. z.o.o.
Tel: +48 22 53 59 700
Fax: +48 22 53 59 703
E-mail: karina.blachnio@merck.pl
www.merck.pl

Portugal

VWR International Lda.
Tel: +351 213600 770
Fax: +351 213600 798/9
E-mail: info@pt.vwr.com
www.vwr.com

Romania

Merck Romania S.R.L.
Tel: +40 21 319 8850
Fax: +40 21 319 8848
E-mail: info.merck@merck.ro
www.merck.ro

Russia

Merck KGaA Representative
Office Moscow
Tel: +7 (495) 935 70 46
Fax: +7 (495) 935 70 44
E-mail: yulia.vaulina@merck.de
www.merck.ru

Serbia and Montenegro

NOVOS d.o.o.
Merck KGaA Representative
Tel: +381 11 316 4339
Fax: +381 11 316 3180
E-mail: novos@beocity.net

Slovakia

Merck spol. s.r.o.
Tel: +421 2 4926 7111
Fax: +421 2 4926 7777
E-mail: merck@merck.sk
www.merck.sk

Slovenia

Merck d.o.o.
Tel: +386 1 560 3 800
Fax: +386 1 560 3 830
E-mail: info.chem@merck.si

South Africa

Merck Laboratory Supplies Division
Tel: +27 11 345 9000
National Contact No. 08600 63725
Fax: +27 11 345 9251
E-mail: labsupply@merck.co.za
www.merck.co.za

Spain

VWR International Eurolab S.L.
Tel: +34 902 19 56 44
Fax: +34 902 19 59 39
E-mail: info@es.vwr.com
www.vwr.com

Sweden

VWR International AB
Tel: +46 8 621 34 00
Fax: +46 8 621 34 66
E-mail: info@se.vwr.com
www.vwr.com

Switzerland

VWR International AG
Tel: +41 44 745 1313
Fax: +41 44 745 1310
E-mail: info@ch.vwr.com
www.vwr.com

Turkey

ALFRED PALUKA & Co.
Merck KGaA Representative
Istanbul
Tel: +90 212 245 1246
Fax: +90 212 252 4458
E-mail: info@alfredpaluka.com

