

1.14428

MQuant®

Ammonium Test

Ammonium-Test

Test Ammonium

Test Amonio

Color card

Farbkarte

Carte colorimétrique

Tarjeta colorimétrica

Supelco®


www.sigmaaldrich.com

NH₄⁺




© 2024

1.



2.



7.69144.2805-XXXXXX 07/2024

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.sigmaaldrich.com/colorimetry

EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321

MilliporeSigma Canada Ltd., 2149 Winston Park Dr,
Oakville, Ontario, L6H 6J8, Canada,
Phone: +1 800-565-1400

0 0.025 0.050 0.075 0.10 0.15 0.20 0.25 0.30 0.40 mg/l NH₄⁺



20 ml



2 ml NH₄-1



1 x NH₄-2




5 min



2 x NH₄-3



7 min !



0 0.02 0.04 0.06 0.08 0.12 0.16 0.19 0.23 0.31 mg/l NH₄-N

