

XXXXXX

7.76067.0007 10/2022

Supelco

sigmaaldrich.com

50 Tests

1.16993

Reflectoquant®

**Calcium Test**

for RQflex® 10 and RQflex® 10 plus

**Calcium-Test**

für RQflex® 10 und RQflex® 10 plus



For IfU please scan code with  
My M Safety App or go to  
sigmaaldrich.com

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

2.5 - 45.0 mg/l Ca

Reaction time: 90 sec



5 ml



2 s



1 x Ca-1



10 x Ca-2



80 s



©2022

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany  
EMD Millipore Corporation, Burlington MA 01803, USA  
Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd. 2149  
Winston Park, Dr. Oakville, Ontario, L6H 6J8