



QUALITY CONTROL RANGES

Human Insulin ELISA

Catalog #: EZHI-14K

QC Catalog #: E6000-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	EUV-129	8.0 - 17	$\mu\text{U/mL}$
QC 2	EUV-229	23 - 47	$\mu\text{U/mL}$

E6000-K: EUV-129.229 (11 February 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.