



## QUALITY CONTROL RANGES

### RAT/MOUSE FGF-21 ELISA CATALOG # EZRMFGF21-26K

#### STANDARD, LYOPHILIZED CATALOG # E8026-K

STANDARD	RAT/MOUSE FGF-21 VALUE
LOT NUMBER: 826-13K	11,411 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6026-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RMFGF-113	142 - 294	pg/mL
QC 2	RMFGF-213	1003 - 2082	pg/mL

E6026-K: RMFGF-113.213 (16 APRIL 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA