

Product Information

UBIQUITIN, N-TERMINAL FLAG-TAGGED

Lyophilized

U5382**Product Description**

The ubiquitin-proteasome pathway is a controlled ATP dependent cellular protein degradation system. This system plays a central role in degradation of short-lived and regulatory proteins important in a variety of basic cellular processes. These processes include amongst others regulation of the cell cycle, modulation of cell surface receptors and ion channels and antigen presentation. The ubiquitin proteasome pathway involves an enzymatic cascade through which multiple ubiquitin molecules are covalently attached to the protein substrate to form a poly-ubiquitin-protein conjugate, which is then degraded by the 26S proteasome complex.^{1,2} Recently it was suggested that exogenous ubiquitin can induce apoptosis in various cell lines.³ Ubiquitin, N-terminal FLAG-tagged (FLAG-ubiquitin, Molecular Weight: approximately 10 kDa) can replace the ubiquitin in formation of poly-ubiquitin-protein conjugates. The FLAG tag enables separation and enrichment of the protein conjugates on anti-FLAG affinity columns and detection of conjugates in western blots using anti-FLAG antibodies.⁴

Reagents

Purity: minimum 95% by SDS-PAGE

Preparation Instructions

Solubility: 10 mg/mL in 50 mM Tris, pH 7.4

Storage/Stability

Store at -20 °C.

References

1. Hershko, A. and Ciechanover, A., Annu. Rev. Biochem, 61, 761-807 (1992).
2. Schwartz, A. and Ciechanover, A, Annu. Rev. Med., 50, 21645-21649 (1999).
3. Matsumura, D.H., et al., Blood, 95, 2577-2585 (2000).
4. Brizzard, B.L., et al., Biotechniques, 16, 730-735 (1994).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.