



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sia.com
sigma-aldrich.com

Product Information

CHLORAMPHENICOL SELECTIVE SUPPLEMENT

Product Number **C 7470**

Product Description

An antibiotic supplement recommended for the selective isolation of yeasts and molds.

Components

(per vial, sufficient for 500 ml of medium)

Item

Chloramphenicol	50.00 mg
-----------------	----------

Preparation Instructions

Rehydrate the contents of 1 vial aseptically with 2 mls of ethanol. Mix well and aseptically add it to 500 mls of sterile, molten (R 1273) Rose Bengal Agar Base. Mix well and pour into sterile petri plates.

Storage Temperature

Store at 2-8 °C.

Precautions and Disclaimer For laboratory use only.
Not for drug, household or other uses.

C7470 6/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.