

Product No. E-9769
Epithelial Neutrophil Activating Peptide-78 (ENA-78)
Recombinant, Human
Expressed in *E. coli*

Description

Epithelial Neutrophil Activating Peptide-78 (ENA-78) was initially discovered from the conditioned medium of human pulmonary epithelial cells (A549) which had been stimulated with TNF- α or IL-1 β .¹ ENA-78 is an 8.3 kD protein with 78 amino acids containing 4 cysteines positioned identically to those of IL-8.1 ENA-78 is a member of the alpha (C-X-C) supergene family which includes platelet factor-4, macrophage inflammatory protein-2α and macrophage inflammatory protein-2β. ENA-78 shares 53% sequence homology with NAP-2 and 52% sequence homology with GROα.1 ENA-78 shares several properties of neutrophil activation with NAP-2 and IL-8. ENA-78 induces chemotactic activity in neutrophils, as well as release of elastase from cytochalasin B-pretreated neutrophils and the induction of cytosolic calcium release.² In response to ENA-78, neutrophils migrate into inflamed joints of patients with rheumatoid arthritis.3

Performance Characteristics

The biological activity of recombinant, human ENA-78 is measured by its ability to induce myeloperoxidase release from cytochalasin B-treated neutrophils. The EC_{50} is defined as the effective concentration of growth factor that elicits a 50% increase in myeloperoxidase release from neutrophils in a cell based bioassay.

Product Information

Expressed in *E. coli* Molecular Weight: 8.3 kD

Purity: ≥ 97% as determined by SDS-PAGE

EC₅₀: 1 - 10 μg/ml Package Size: 10 μg/vial Formulation: Lyophilized from a $0.2\ \mu m$ -filtered

solution of PBS, pH 7.4.

Carrier Protein: 500 µg bovine serum albumin (BSA)

Sterility: 0.2 µm-filtered, aseptic fill Endotoxin: ≤0.1 ng/µg ENA-78

Reconstitution and Use

Reconstitute the contents of the vial using 0.2 μ m-filtered PBS containing 0.1% HSA or BSA to a concentration of not less than 10 μ g/ml.

Storage

Prior to reconstitution, store at -20° C for no more than 6 months. After reconstitution, store at 2-8°C for a maximum of one months. For extended storage, freeze in working aliquots at -70° C or -20° C. Repeated freezing and thawing is not recommended.

References

- 1. Walz, A., et al., J. Exp. Med., **174**, 1355 (1991).
- Walz, A., et al., in Chemotactic Cytokines, Plenum Publishing Corp. (In press).
- 3. Koch, A., et al., J. of Clin. Invest., **94**, 1012 (1994).
- 4. Schröeder, J., et al., J. Immunol., **139**, 3474 (1987)

Sigma warrants that its products conform to the information contained in this and other Sigma publications. Purchaser must determine the suitability of the products for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale. Issued 05/95.