

02287 Park and Sanders Selective Supplement II, modified

An antibiotic supplement recommended for the selective isolation and enumeration of thermotolerant *Campylobacter* species .

Composition
(per vial sufficient for 1000 ml medium)
Cefoperazone 32mg
Amphotericin B 5mg

Directions
Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add it to 1000 ml of sterile *Campylobacter* Selective Agar, Base (Cat. No. 21378) containing 50- 70 ml sterile defibrinated blood and Park and Sanders Selective Supplement I (Cat. No. 17191). Mix well and dispense.

