



QUALITY CONTROL RANGES

High Sensitivity Human Amyloid β 40 ELISA Kit
Catalog# EZHS40

Human Amyloid β 40 Brain ELISA Kit
Catalog#EZBRAIN40

Control Catalog# 0QC

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	0QC-113	24 - 51	pg/mL
QC 2	0QC-213	94 - 196	pg/mL

0QC: 0QC-113.213 (02 October 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA