



RABBIT ANTI- CANNABINOID RECEPTOR 1 AFFINITY PURIFIED POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB9415
LOT NUMBER:	
QUANTITY:	50 µg
CONCENTRATION:	1.0 mg/mL
SPECIFICITY:	Recognizes Cannabinoid Receptor 1 [CB1].
IMMUNOGEN:	Synthetic peptide from the 3 rd cytoplasmic domain of human Cannabinoid Receptor 1.
APPLICATIONS:	Immunohistochemistry on formalin fixed paraffin embedded tissue sections. Immunocytochemistry on acetone fixed cells. Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITIES:	Human, rat and monkey. Other species have not been tested. The immunogen sequence has 100% homology with mouse.
FORMAT:	Affinity purified immunoglobulin.
PRESENTATION:	Liquid in PBS, pH 7.7 with 0.01% sodium azide.
STORAGE/HANDLING:	Maintain at -70°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com

AB9415/8-3-2005/CB