

QUALITY CONTROL RANGES



MILLIPLEX[®] MAP Mouse Kidney Injury Magnetic Bead Panel 2

Kit Catalog # MKI2MAG-94K

Control Catalog # MKI2-6094-2
Lot# MKI2-104 and MKI2-204

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Cystatin C	Control 1	0.09 – 0.19	ng/mL
	Control 2	0.89 – 1.9	ng/mL
EGF	Control 1	0.32 – 0.66	ng/mL
	Control 2	2.8 – 5.7	ng/mL
Lipocalin-2 /NGAL	Control 1	0.03 – 0.07	ng/mL
	Control 2	0.33 – 0.68	ng/mL
Clusterin	Control 1	1.8 – 3.6	ng/mL
	Control 2	16 – 33	ng/mL
OPN	Control 1	0.07 – 0.15	ng/mL
	Control 2	0.65 – 1.4	ng/mL

MKI2-6094-2: MKI2-104.204 (29 April 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.