



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sia.com
sigma-aldrich.com

Product Information

STEAROYL CHLORIDE

Product Number **S 5253**

The melting point of this product is approximately 22 °C. Prior to opening ampule, liquify contents by incubating at 37 °C.

Heating in hot water baths (100 °C) is not recommended as pressure may build and burst the ampule.

KMR 03/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.