

MOUSE ANTI INFLUENZA A H1N1 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8261	QUANTITY:	100 µg
LOT NUMBER:		CONCENTRATION:	1 mg/mL
CLONE NAME:	9B3.2	HOST/ISOTYPE:	Ms/ IgG _{2a}
SPECIFICITY:	Influenza A H1N1 antigen. No reactivity shown to Influenza B strains.		
IMMUNOGEN:	Influenza blend		
APPLICATIONS:	IF Optimal dilutions must be determined by end user		
SPECIES REACTIVITY:	Reacts strongly with the following H1N1 strain: Beijing. A/ Texas /36/91, A/Berkeley/1/98, A/HongKong/503/97, A/Nanchang /16A/98, A/PR8/34 and New Caledonia strain.		
FORMAT:	Purified immunoglobulin		
PRESENTATION:	Liquid in 0.02 M Phosphate Buffer + 0.25 M Sodium Chloride + 0.1% Sodium Azide, pH = 7.6		
STORAGE/HANDLING:	Store at 2° to 8°C for up to 12 months from date of receipt.		

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.