



MOUSE ANTI-RESPIRATORY SYNCYTIAL VIRUS (RSV) MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB858-4
LOT NUMBER:	
VOLUME:	100 µL
SPECIFICITY:	Blend of all RSV specifications of MAB858-1, 2 and 3. Recognizes the fusion, G, and nuclear protein of both A and B strains.
IMMUNOGEN:	A2
ISOTYPE:	IgG
CLONE NAME:	Blend of MAB858-1, 2 and 3 Clones 133-1H, 131-2G, and 130-12H
APPLICATIONS:	ELISA at 1:800+ Indirect Immunofluorescence at 1:100-200+ (fresh frozen tissue sections) Final working dilutions must be determined by end user.
FORMAT:	Ascites fluid.
PRESENTATION:	Liquid
STORAGE/HANDLING:	Maintain at -20°C in aliquots for up to 12 months. Short term storage, 3 weeks, at +2-8°C is acceptable, with the addition of 0.1% sodium azide suggested. Do not store in diluted format. Avoid repeated freeze/thaw cycles.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.