



## CHICKEN ANTI-LDL RECEPTOR AFFINITY PURIFIED POLYCLONAL ANTIBODY

<b>CATALOG NUMBER:</b>	AB3461
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	2 x 50 µg
<b>CONCENTRATION:</b>	1 mg/mL
<b>SPECIFICITY:</b>	LDL Receptor
<b>IMMUNOGEN:</b>	Human LDL receptor, amino acids 29-205.
<b>APPLICATIONS:</b>	Western blot: 1:500-1:2,000 (colormetric)  It is expected that the antibody will also work for immunocytochemistry and immunohistochemistry. Specific reactivity has not been confirmed. Optimal working dilutions must be determined by the end user.
<b>SPECIES REACTIVITY:</b>	Human. The antibody may cross react with mouse and rat.
<b>FORMAT:</b>	Affinity purified immunoglobulin.
<b>PRESENTATION:</b>	Liquid in PBS, containing 0.02 % sodium azide, pH 7.2
<b>STORAGE/HANDLING:</b>	Maintain at -20°C in undiluted aliquots up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
<b>BACKGROUND REFERENCE:</b>	Agnello, V., et al., <i>PNAS</i> (1999) <b>96</b> :12766-12771.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002-2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209  
Australia +61 3 9839 2000  
www.millipore.com