

## QUALITY CONTROL RANGES

### Glucagon Chemiluminescent ELISA Catalog #: EZGLU-30K

### GLUCAGON CHEMILUMINESCENT STANDARD, LYOPHILIZED CATALOG#: E8030-K

STANDARD	Glucagon Concentration
Lot Number: E8030-10K	2.0 ng/mL

The exact value of lyophilized standard is shown above.  
Reconstitute with **2.0mL** distilled or deionized water.

### Control Catalog#: E6030-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GLU-110	0.08 – 0.17	ng/mL
QC 2	GLU-210	0.67 – 1.4	ng/mL

E6030-K: GLU-110.210 (18 December 2025)

**SigmaAldrich.com/support • SigmaAldrich.com**

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.