



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Rat/Mouse Neuropeptide Magnetic Bead Panel
Kit Catalog # RMNPMAG-83K
Control Catalog # RMNP-6083
Lot# RMNPM-108 and RMNPM-208

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
α -MSH	Control I	400 – 831	pg/mL
	Control II	2084 – 4328	pg/mL
β -Endorphin	Control I	660 – 1447	pg/mL
	Control II	3991 – 8289	pg/mL
Neurotensin	Control I	221 – 458	pg/mL
	Control II	1290 – 2680	pg/mL
Oxytocin	Control I	122 – 253	pg/mL
	Control II	658 – 1366	pg/mL
Substance P	Control I	23 – 47	pg/mL
	Control II	135 – 281	pg/mL

RMNP-6083: RMNPM-108.208 (05 February 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA