

38390 Preston Selective Supplement

This supplement is recommended for selective isolation of *Campylobacter* species.

Composition

(Per vial sufficient for 500 ml medium)

Polymyxin B 2500IU

Rifampicin 5mg

Trimethoprim lactate 5mg

Amphotericin B 5mg

Directions

Rehydrate the contents of 1 vial aseptically with 2.5 ml of ethanol 95%. Mix well (few particulates may be observed) and aseptically add it to 500 ml sterile, molten, cooled (45-50°C) Preston Broth Base (74173). Mix well and dispense into sterile tubes/ petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

