

67728 *C. perfringens* Millichrome™ plus Supplement 1

Selective supplement for detection and enumeration of *C. perfringens*

Package information:

1 vial à 10g (with screw cap) is sufficient for 5L

Directions

- In a transparent vessel, add 2 g of Millichrome™ *C. perfringens* supplement 1 in 20 mL of purified water.
- Swirl well until complete dissolution.
- Filter to sterilize at 0.45 µm.
- Add 20 mL/L of the supplement 1 solution to the melted *C. perfringens* Millichrome™ plus Agar base at 45/50 °C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

