



MOUSE ANTI-INFLUENZA A MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8254
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Influenza A H3 antigen. Reacts strongly with the following H3N2 strains: A/Wuhan/396/95; A/Johannesburg/33/94; A/Shandong/93; A/Shanghai/16/89; All CA Strains 1986 to Present. No reactivity shown with the following H1N1 strains: A/Taiwan/1/86; A/Texas/36/91; A/USSR/90/77; All CA strains 1986 to present. No reactivity shown to Influenza B strains.
IMMUNOGEN:	Influenza blend.
ISOTYPE:	IgG _{2a}
FORMAT:	Purified immunoglobulin from Protein A Chromatography.
STORAGE/HANDLING:	Maintain at 2-8°C for up to 12 months.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209
www.millipore.com