

Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

Interleukin-2 from mouse, recombinant expressed in *E. coli*

Product Number **I0523** Storage Temperature –20 °C

Synonyms: IL-2; T Cell Growth Factor; Aldesleukin

Product Description

Interleukin-2 is a T cell-derived cytokine, which was first described as a T cell growth factor (TCGF). In addition to its potent T cell growth stimulatory activity, multiple biological functions of IL-2 also include B cell growth and differentiation, and generation of LAK cells and NK cells. Natural mouse IL-2 is produced from T cells. Recombinant, mouse IL-2 is produced from a DNA sequence encoding the mature mouse IL-2. Mouse IL-2 contains 149 amino acids with disulfide bonds between Cys⁷² and Cys¹²⁰. The disulfide bond is essential for its biological activity. The amino acid sequence of mouse IL-2 has 63% homology with human IL-2. Human IL-2 is active on both human and mouse lymphocytes. Mouse IL-2 is active on mouse lymphocytes.

Molecular Weight: 17.2 kDa

Each vial contains 20 μ g of Interleukin-2 with 1 mg of bovine serum albumin lyophilized from a 0.2 μ m-filtered solution of 10 mM sodium citrate, pH 4.0. Sterility: 0.2 μ m-filtered, aseptic fill

Purity: ≥97% (SDS-PAGE)

Endotoxin: \leq 0.1 ng/ μ g IL-2 (1 EU/ μ g)

ED₅₀: =0.2 ng/ml (corresponds to $\geq 5 \times 10^6$ units/mg) The biological activity of recombinant, mouse IL-2 is measured by its ability to stimulate the proliferation of the IL-2 dependent mouse cytotoxic T cell line, CTLL-2. The ED₅₀ is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell based bioassay.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Centrifuge vial prior to opening. Reconstitute the contents of the vial in water to a concentration range of 0.1–1.0 mg/ml.

Storage/Stability

Store the lyophilized powder at $-20\,^{\circ}$ C. After reconstitution, store at 2–8 $^{\circ}$ C for a maximum of one week. For extended storage, freeze in working aliquots at $-70\,^{\circ}$ C or $-20\,^{\circ}$ C for a maximum of 1 year. Repeated freezing and thawing is not recommended.

References

- 1. Hatakeyama, M., *et al.*, Peptide Growth Factors and their Receptors I, Sporn, M., *et al.*, eds., Springer-Verlag (New York, NY: 1991) p. 523.
- 2. Kashima, N., et al., Nature, 313, 402 (1985).
- 3. Yokata, T., et al., Proc. Natl. Acad. Sci. USA, **82**, 68 (1986).
- Gearing, A., et al., Lymphokines and Interferons, A Practical Approach, Clemens, M., et al., eds., IRL Press (Oxford, UK: 1987) p. 296.

LCM,MAM 12/05-1