

XXXXXX

7.76048.0005 01/2025

Supelco

sigmaaldrich.com

50 Tests

1.16982


Reflectoquant®

Iron Test for RQflex® 20

Eisen-Test für RQflex® 20

0.5 - 20.0 mg/l Fe²⁺

Reaction time: 15 sec



Calculation of the iron(III) content:

$$\text{mg/l Fe}^{3+} = \text{mg/l total iron} - \text{mg/l Fe}^{2+}$$

Berechnung des Eisen(III)-Gehalts:

$$\text{mg/l Fe}^{3+} = \text{mg/l Gesamteisen} - \text{mg/l Fe}^{2+}$$

Please consult Instruction for Use for further information and observe the RQflex® manual.



For IfU please scan code with My M Safety App or go to sigmaaldrich.com



©2025

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany
 EMD Millipore Corporation, Burlington MA 01803, USA
 MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
 Ontario, L6H 6J8, Canada