3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

Anti-Netrin-G1a

produced in goat, affinity isolated antibody

Catalog Number N4037

Product Description

Anti-Netrin-G1a (Mouse) was produced in goat using purified insect cell line *Sf* 21-derived, recombinant mouse Netrin-G1a (rmNetrin-G1a) as the immunogen. . The antibody was affinity isolated.

Anti-Netrin-G1a (Mouse) recognizes mouse Netrin-G1a in direct ELISAs and western blots. This antibody shows less than 2% cross-reactivity with rmNetrin-1, rchNetrin-2 and rhNetrin-4.

Semaphorins, neuropilins and netrins are among a number of molecules and their receptors that regulate the developing nervous system to guide the development of neural circuits.¹ Although first identified as axon guidance cues,^{2,3} it is now apparent that many of these same factors are not limited to the guidance of growing axons, but have roles in a range of processes from the guidance of cell migration to the regulation of the immune response, angiogenesis, lung branching morphogenesis, nervous system regeneration, and cancer.⁴⁻⁹

The semaphorins make up the largest family of axon guidance cues. They are characterized by the presence of an approximately 500 amino acid N-terminal semaphorin (Sema) domain. Semaphorins function mainly as chemorepellents that direct axons away from tissues.³ Semaphorin 3A (Sema3A) has been shown to be repellent to cortical axons and to inhibit axon branching.¹⁰ The transmembrane protein semaphorin 6A has been shown to repel embryonic sympathetic axons.¹¹ The actions of the various semaphorin 3A has been found to inhibit tumor development whereas semaphorin 6A may contribute to tumor progression.⁹

Neuropilins are the ligand binding moieties in the class 3 Semaphorin receptor complexes that subsequently activate signaling through associated plexins. Two types have been identified so far: Neuropilin-1 (Npn-1) and Neuropilin-2 (Npn-2) receptors. At the amino acid sequence level, Npn-1and Npn-2 share 44% identity. Npn-1 and Npn-2 show different expression patterns in developing neurons of the central and peripheral nervous systems, and show different binding specificities for different members of the semaphorin family. Both also function as receptors for some forms of vascular endothelial growth factor (VEGF).¹²

Netrins are a family of laminin-related small proteins that are involved in axon guidance and neurite outgrowth. Netrin-1 has been shown to attract cortical growth cones and promote axon branching.¹⁰ Netrin-4 (first named β -netrin) was found to promote neurite elongation from olfactory bulb explants.¹³ Netrin-G1 consists of at least six isoforms of which five are predominantly anchored to the plasma membrane via glycosyl phosphatidyl-inositol linkages, and lack appreciable affinity to receptors for classical netrins.^{14,15}

Reagent

The antibody is supplied lyophilized from a 0.2 μm filtered solution in phosphate-buffered saline (PBS) with 5% trehalose.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reconstitute with sterile PBS. If 1 mL of PBS is used, the antibody concentration will be 100 μ g/mL.

Storage/Stability

Lyophilized powder can be stored intact at room temperature for several weeks. For extended storage, it should be stored at -20 °C or below. The reconstituted solution can be stored at 2-8 °C for up to 2 weeks. For longer storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

The recommended working dilutions are 0.1-0.2 μ g/ml for immunoblotting and 0.5-1 μ g/ml for Direct ELISA.

Note: In order to obtain best results and assay sensitivities of different techniques and preparations, determination of optimal working dilutions by titration test is recommended.

References

- Tsim, T.Y., et al., Expression of axon guidance molecules and their related genes during development and sexual differentiation of the olfactory bulb in rats, *Neuroscience* **123**, 951-965 (2004).
- de Castro F., Chemotropic molecules: guides for axonal pathfinding and cell migration during CNS development, *News Physiol. Sci.* 18, 130-136 (2003).
- 3. Fiore, R. and Puschel, A.W., The function of semaphorins during nervous system development, *Front. Biosci.* **1**, s484-499 (2003).
- Tamagnone, L. and Comoglio, P.M., To move or not to move? Semaphorin signalling in cell migration, *EMBO Rep.* 5, 356-361(2004).
- Eichmann, A., et al., Guidance of vascular and neural network formation, *Curr. Opin. Neurobiol.* 15, 108-15 (2005).
- Goshima, Y., et al., Semaphorins as signals for cell repulsion and invasion, *J. Clin. Invest.* **109**, 993– 998 (2002).
- De Winter, F., et al., Neuropilin and class 3 semaphorins in nervous system regeneration, *Adv. Exp. Med. Biol.* **515**, 115-139 (2002).
- de Wit, J. and Verhaagen, J., Role of semaphorins in the adult nervous system, *Prog. Neurobiol.* **71**, 249-267 (2003).

- 9. Neufeld, G., et al., Semaphorins in cancer, *Front. Biosci.* **10**, 751-760 (2005).
- Dent, E.W., et al., Netrin-1 and Semaphorin 3A Promote or Inhibit Cortical Axon Branching, Respectively, by Reorganization of the Cytoskeleton, *J. Neurosci.* 24, 3002-3012 (2004).
- 11. Xiao-Mei, X., et al., The Transmembrane Protein Semaphorin 6A Repels Embryonic Sympathetic Axons, *J. Neurosci.* **20**, 2638-2648 (2000).
- Neufeld, G., et al., The interaction of Neuropilin-1 and Neuropilin-2 with tyrosine-kinase receptors for VEGF, *Adv. Exp. Med. Biol.* 515, 81-90 (2002).
- 13. Koch, M., et al., A novel member of the netrin family, β -netrin, shares homology with the β chain of laminin: identification, expression, and functional characterization, *J. Cell Biol.* **151**, 221-234 (2000).
- Nakashiba, T., et al., Netrin-G1: a Novel Glycosyl Phosphatidylinositol-Linked Mammalian Netrin That Is Functionally Divergent from Classical Netrins, *J. Neurosci.* 20, 6540-6550 (2000).
- 15. Nakashiba, T., et al., Complementary expression and neurite outgrowth activity of netrin-G subfamily members, *Mech. Dev.* **111**, 47-60 (2002).

AH,PHC 09/10-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.