

## Certificate of Analysis

### Anti-Caspase-8 (N-terminus), clone E6, rabbit monoclonal

(Rabbit monoclonal IgG)

Catalog # 04-574

Lot #

**Immunogen:** A synthetic peptide corresponding to N-terminal residues of human Caspase-8. The antibody should recognize all splice isoforms of Caspase-8.

**Specificity:** Recognizes human Caspase-8 (N-Term).

**Molecular Weight:** 55 kDa.

**UniProt ID:** Q14790

**Species Cross-reactivity:** Reacts with human. Not believed to cross-react with mouse or rat.

**Formulation:** 100  $\mu$ L of rabbit monoclonal IgG in 50mM Tris-Glycine (pH 7.4), 0.15 M NaCl, 0.01% sodium azide, 0.05% BSA and 40% glycerol.

**Storage and Stability:** Stable for 2 years at -20°C from date of shipment.

**Handling Recommendations:** Upon first thaw, and prior to removing the cap, centrifuge the vial and gently mix the solution. Aliquot into microcentrifuge tubes and store at -20°C. **Avoid repeated freeze/thaw cycles, which may damage IgG and affect product performance.** Note: Variability in freezer temperatures below -20°C may cause glycerol-containing solutions to become frozen during storage.

**FOR RESEARCH USE ONLY  
NOT FOR USE IN HUMANS**


### Applications:

Western blot: A 1:500-1:1000 dilution

Immunoprecipitation: 1:50 dilution

Immunohistochemistry: 1:100 dilution

Flow cytometry: 1:50 dilution



### Quality Control Testing:

**Western blot Analysis:** Lysate from Jurkat cells was resolved by electrophoresis, transferred to PVDF and probed with anti-caspase-8 (1:1,000). Proteins were visualized using goat anti-rabbit secondary antibody conjugated to HRP and chemiluminescence detection. Arrow indicates caspase-8, (55 kDa).

### General References:

1. Medema, J.P. et al. (1997). FLICE is activated by association with the CD95 death-inducing signaling complex (DISC). *EMBO J.* **16**(10): 2794–804.
2. Muzio, M. et al. (1996). FLICE, a novel FADD-homologous ICE/CED-3-like protease, is recruited to the CD95 (Fas/APO-1) death-inducing signaling complex. *Cell* **85**: 817–827.
3. Boldin, M.P. et al. (1996). Involvement of MACH, a novel MORT1/FADD-interacting protease, in Fas/APO-1 and TNF receptor-induced cell death. *Cell* **85**: 803–815.
4. Srinivasula, S.M. et al. (1996). Molecular ordering of the Fas-apoptotic pathway: the Fas/APO-1 protease Mch5 is a CrmA-inhibitable protease that activates multiple Ced-3/ICE-like cysteine proteases. *Proc. Natl. Acad. Sci. U S A* **93**(25): 14486–91.

### Immunoblot Protocol

1. Perform SDS-polyacrylamide gel electrophoresis (SDS-PAGE) on cell lysate samples (cell lysis buffer: 50 mM Tris-HCl, pH 7.4; 1% NP-40; 0.25% sodium deoxycholate; 150 mM NaCl; 1 mM EDTA; 1 mM PMSF; 1  $\mu$ g/ml each aprotinin, leupeptin, pepstatin; 1 mM  $\text{Na}_3\text{VO}_4$ , 1 mM NaF) and transfer the proteins to PVDF. Wash the blotted PVDF twice with TBST.
2. Block the blotted PVDF in freshly prepared 5% BSA in TBS with 0.05% Tween<sup>®</sup>-20 for 1 hour at room temperature with constant agitation.
3. Incubate the PVDF with 1:1,000 dilution of **anti-caspase-8** diluted in freshly prepared 5% BSA in TBST for 2 hours at room temperature or overnight with agitation at 4°C.
4. Wash the PVDF 3 times with TBST.
5. Incubate the PVDF in the secondary reagent of choice (a goat anti-rabbit HRP conjugated IgG, Catalog # 12-348, 1:5000 dilution was used) in TBST/5% BSA for 1 hour with agitation at room temperature.
6. Wash the PVDF 3-5 times with TBST.
7. Use detection method of choice (enhanced chemiluminescence was used).



**Rabbit Monoclonals Produced Using Technology from Epitomics, Inc. Under Patent No. 5,675,063**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

