

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

### **GOAT SERUM**

Product Number G 9023

### **Product Description**

Goat Serum is a pooled preparation obtained from a normal donor herd.

Immunoelectrophoresis against anti-goat serum results in multiple arcs of precipitation. The presence of endogenous antibodies has not been assayed.

Normal goat serum is often used as a blocking agent or negative control in immunoassays.

## Reagents

This product is provided as a liquid containing preservative.

### **Precautions and Disclaimer**

This product is for research use only. Consult the MSDS for information regarding hazards and safe handling practices.

### Storage/Stability

Store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended.

JWM,DAA,SM 08/20-1