

Product Information Sheet

Interleukin-2 (IL-2)

Human, Recombinant

SRP3085

Product Description

Recombinant Human Interleukin-2 is a 15.5 kDa protein containing 134 amino acid residues including one intrachain disulfide bond.

Interleukin-2 (IL-2), also known as T Cell Growth Factor (TCGF) and aldesleukin, is an immunomodulatory lymphokine produced by T cells in response to antigenic or mitogenic stimulation.¹ IL-2/IL-2 receptor signaling is required for T cell proliferation and other functions, which are essential for immune response. Interleukin-2 promotes growth and differentiation of B cells, NK cells, lymphokine-activated killer (LAK) cells, monocytes, macrophages, and oligodendrocytes. It has been shown to affect the activation and proliferation of NK cells, induce γ -interferon and B cell growth factor secretion,^{2,3} and modulate the expression of the IL-2 receptor (IL-2R).⁴ Interleukin-2 has been isolated from a number of cell types^{5,6} and has been reproduced by recombinant DNA technology.⁷

Reagent

Recombinant Human Interleukin-2 is supplied as lyophilized with no additives. It is sterile filtered through a 0.2 μ m filter.

Storage/Stability

Prior to reconstitution, store the lyophilized protein at -20 °C. It is stable for up to a few weeks at room temperature, but is best stored at -20 °C. For extended storage, after reconstitution, store in working aliquots at -20 °C. Repeated freezing and thawing is not recommended.

Preparation Instructions

Reconstitute the contents of the vial using 100 mM acetic acid to a concentration of 0.1-1.0 mg/mL. This solution can then be diluted into other aqueous buffers and stored at 2-8 °C for up to 1 week. For extended storage, freeze in working aliquots. It may be advisable to centrifuge the vial prior to reconstitution.

Product Profile

The biological activity of Recombinant Human Interleukin-2 is measured by the dose-dependent stimulation of the proliferation of the IL-2 dependent mouse cytotoxic T-cell line, CTLL-2.⁸ The EC₅₀ is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell-based bioassay.

References

1. Smith, K. A., Science, 240, 1169 (1988).
2. Morgan, D. A., et al., Science, 193, 1007 (1976).
3. Inaba, K., et al., J. Exp. Med., 158, 2040 (1983)
4. Smith, K. A., et al., Proc. Natl. Acad. Sci. USA, 82, 864 (1985).
5. Kniep, E. M., et al., Eur. J. Biochem., 143, 199 (1984).
6. Robb, R. J., et al., Proc. Natl. Acad. Sci., USA, 80, 5990 (1983)
7. Taniguchi, T., et al., Nature, 302, 305 (1983)
8. Gearing, A., et al., In: Lymphokines and Interferons, A Practical Approach, Clemens, M., et al. (eds.), 296 (IRL Press, Oxford, 1987).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.