

43963 Moxalactam Supplement

An antibiotic supplement recommended for selective isolation and cultivation of *Listeria monocytogenes*.

Composition

(per vial, sufficient for 1000 ml medium)

Moxalactam 20.0 mg

Directions

Rehydrate the contents of 1 vial with 10 ml of sterile distilled water.

Mix well and add it to 1000 ml of LPM Agar (Cat. No. 62653).

Mix well and pour
into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

