

80704 VCC Selective Supplement (Vancomycin Cefixime Cefsulodin Selective Supplement; ECO157:H7 Selective Supplement)

An antibiotic supplement recommended for selective isolation of Escherichia coli O157:H7.

Composition

(per vial, sufficient for 1000 ml medium)

Vancomycin 8.0 mg

Cefixime 0.05 mg

Cefsulodin 10.0 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix gently to dissolve the contents completely and aseptically add the rehydrated contents to 1000 ml of sterile, cooled (45-50°C) Modified tryptone soy broth (Cat. No 08069).

Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

