

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Monoclonal Anti-Vasoactive Intestinal Peptide Receptor 2 (VIPR2), Clone AS69 Purified Mouse Immunoglobulin

Product Number V 6257

Product Description

Monoclonal Anti-Vasoactive Intestinal Peptide Receptor 2 (VIPR2) (mouse IgG1 isotype) is derived from the hybridoma produced by the fusion of mouse myeloma cells and splenocytes from a BALB/c mouse immunized with a unique peptide corresponding to a portion of rat Vasoactive Intestinal Peptide Receptor 2 (VIPR2). The antibody was purified from tissue culture supernatant using immobilized protein A.

Monoclonal Anti-Vasoactive Intestinal Peptide Receptor 2 (VIPR2) recognizes VIPR2 protein from rat tissue by immunoblotting and immunohistochemistry. The antibody does not recognize VIPR1 protein.

Vasoactive Intestinal Peptide (VIP) is a 28 amino acid neuropeptide with a broad range of biological activities. VIP is widely distributed throughout the central nervous system (CNS), peripheral nervous system and nonneuronal tissue such as mast cells and leukocytes.¹⁻⁶ In the CNS, VIP acts as a neurotransmitter or neuromodulator.^{1,2} It effects behavior, secretion, metabolism, electrophysiology, and secretion. It has trophic and mitogenic activity on neural tissue during embryogenesis, and inhibits tumors.^{7,8} VIP also has important functions in the immune⁹, cardiovascular¹⁰, reproductive¹¹, pulmonary¹², and gastrointestinal systems. General effects include immunosuppression, vasodilation, hormone secretion, bronchodilation, and increased gastric motility.

VIP binding sites in tissues can be divided into two types based on their sensitivity to GTP. Two GTP sensitive VIP receptors have been cloned and characterized, VIPR1 and VIPR2. These have been found to be G protein coupled receptors linked to adenylate cyclase. VIPR1 and VIPR2 are also referred to as VPAC1 and VPAC2 respectively since they respond similarly to VIP and pituitary adenylate cyclase activating polypeptide (PACAP) in stimulating cAMP production.¹³ Although structurally related, VIPR1 and VIPR2 exhibit differences in expression and function and VIP has a 3-10 fold preference for VIPR1 over VIPR2 receptors.

VIPR1 is expressed throughout the central nervous system (predominantly in the cerebral cortex and hippocampus), in peripheral tissues including liver, lung and intestine and in T lymphocytes.¹⁴ VIPR1 mediates suppression of chemotaxis and matrix metalloprotein-ase expression elicited by some cytokines and chemo-kines, tumor cell migration induced by VIP, and vaso-dilation.

VIPR2 is expressed throughout the central nervous system, but to varying degrees. The highest expression levels are in the thalamus and suprachiasmatic nucleus, but VIPR2 is also present in the hippocampus, brainstem, spinal cord and dorsal root ganglia.¹⁵ VIPR2 is also expressed in several peripheral tissues including pancreas, skeletal muscle, heart, kidney, adipose tissue, testis and stomach.¹³

Reagents

Monoclonal Anti-Vasoactive Intestinal Peptide Receptor 2 (VIPR2) is supplied as 100 μ g of purified mouse immunoglobulin (IgG1) at 1.0 mg/ml in 0.2 μ m sterile-filtered phosphate buffered saline with 0.1% sodium azide.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling

Storage/Stability

Store at –20°C. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

The recommended working concentration is $0.5-2 \mu g/ml$ for immunoblotting.¹⁶ This antibody has been used at $0.25-2 \mu g/ml$ along with the appropriate secondary antibodies to stain rat VIPR2 expressed on rat hippocampal cells fixed with paraformaldehyde followed by 0.1% Triton X-100.¹⁶

Note: In order to obtain best results and assay sensitivities of different techniques and preparations, we recommend determining optimal working dilutions by titration test.

References

- Sundler, F., et al., In: "Ann. N.Y. Acad. Sci." Vol. 52: Vasoactive Intestinal Peptide and Related Peptides, Said, S., and Chutt, V., (Eds.), New York Academy of Science, New York, p. 143 (1988).
- Abrams, G., et al., In: "Hanbook of Chemical Neuroanatomy" Vol. 4: GABA and Neuropeptides in the CNS, Part I, Björklund, A., and Hökfelt, T., (Eds.), Elsevier Science Publishers B.V., Amsterdam, p. 335 (1985).

- 3. Hökfelt, T., et al., Neuroscience, 2, 885 (1977).
- 4. Lorén, I., et al., Neuroscience, 4, 1953 (1979).
- 5. Sims, K., et al., Brain Research, 186, 165 (1980).
- Shimatsu, A., et al., Endocrinology, **108**, 395 (1981).
- 7. Gressens, P., et al., Nature, 362, 155-158 (1993).
- 8. Alle, J., et al., J. Neurooncology., 3, 197 (1985).
- 9. De la Fuente, M., et al., Adv. Neuroimmunol, **6**, 75-91 (1996).
- Ganong, W.F., et al., J Hypertens Suppl., 2, S15-S23 (1984).
- 11. Fahrenkrug, J., et al., Ann. NY Acad. Sci., **527**, 393-404 (1988).
- 12. Lundberg, J.M., et al., Ann. NY Acad. Sci., **629**, 332-339 (1991).
- 13. Harmar, A.J., et al., Pharmacol. Rev., **50**, 265-270 (1998).
- 14. Ishihara, T. et al., Neuron, 8, 811-819 (1992).
- 15. Sheward, W.J. et al., Neuroscience, **67**, 409-418 (1995).
- Goetzl, E.J. et al., Mol. Cell. Neurosci., 5, 145-152 (1994).

mje 8/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.