

44051 Thayer Martin Supplement I

An antibiotic supplement recommended for the selective isolation of *Neisseria gonorrhoeae* and *Neisseria meningitidis*.

Composition

(per vial, sufficient for 500 ml medium):

Vancomycin	1.5mg
Colistin methane sulphonate	3.75mg
Trimethoprim	2.5mg
Nyastin	6'250 units

Directions

Rehydrate the contents of one vial aseptically with 5 ml sterile water. Mix well and add to 240 ml sterile, molten Thayer-Martin Medium Base (Prod No 88296)/Chocolate Agar Base along with the rehydrated contents of one bottle of Thayer-Martin Supplement II (Prod No 52083) and 250 ml of sterile 2% hemoglobin (Prod No 08449) solution. Mix gently and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

