



Quality Control Ranges

High Sensitivity Human Amyloid β 42 ELISA Kit Catalog# EZHS42

Human Amyloid β 42 Brain ELISA Kit Catalog#EZBRAIN42

STANDARD, LYOPHILIZED CATALOG # 2STM

Reconstitute standard with **0.25 mL** distilled or deionized water.

STANDARD	Amyloid β 1-42 Standard Value
LOT NUMBER: 2STM-11K	793 ng/mL

Standard	Concentration of Amyloid β 1-42	Volume of Standard and Sample Diluent to Add	Volume of Standard to Add
Stock (S)	25,000 pg/mL	968 μ L	32 μ L of reconstituted Standard

CONTROL CATALOG# 2QC

Control	Lot number	Expected range	Units
QC 1	2QC-111	19 - 40	pg/mL
QC 2	2QC-211	82 - 171	pg/mL

2QC: 2QC-111.211 (28 June 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) •
emdmillipore.com