



QUALITY CONTROL RANGES
Milliplex Rat Vascular Injury Panel 1
Catalog # RV1MAG-26K

Control Catalog # RV1-6026, Lot # RV1-105 and RV1-205

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CAV-1	Control 1	8.2 – 17	ng/mL
	Control 2	63 – 131	ng/mL
CTGF	Control 1	1.6 – 3.3	ng/mL
	Control 2	12 – 26	ng/mL
IL-6	Control 1	0.90 – 1.9	ng/mL
	Control 2	6.7 – 14	ng/mL
MCP-1	Control 1	0.92 – 1.9	ng/mL
	Control 2	7.0 – 15	ng/mL
TNF- α	Control 1	0.16 – 0.33	ng/mL
	Control 2	1.3 – 2.6	ng/mL
GRO/KC/CINC-1	Control 1	0.89 – 1.9	ng/mL
	Control 2	6.9 – 14	ng/mL
TIMP-1	Control 1	1.4 – 2.9	ng/mL
	Control 2	11 – 22	ng/mL
tPAI-1	Control 1	0.85 – 1.8	ng/mL
	Control 2	6.2 – 13	ng/mL
VEGF	Control 1	0.13 – 0.26	ng/mL
	Control 2	1.0 – 2.0	ng/mL

Note: The Quality Control Ranges are generated using serum matrix provided in the kit.
Quality Control Ranges in culture media are not tested.

RV1-6026: RV1-105.205 (17 April 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA