

Product Information

Tissue Section Slide, Multiple Organs from human adult normal neural tissues

Product Number **H 2661**
Store at Room Temperature

Product Description

This slide contains multiple paraffin embedded sections prepared from different human normal tissues. The tissue sections are fixed in paraformaldehyde after excision and then sections approximately 5 µm thick are embedded in paraffin. The tissue slide is suitable for most immunohistochemistry applications and as an immunohistochemistry control slide. The slide may also be used for *in situ* hybridization and *in situ* PCR studies, including tissue specific protein and mRNA expression comparing the location of novel genes among different tissues.

Component

The major organ multiple tissue slides are offered as individual slides.

Storage/Stability

It is recommended to store the slide at room temperature.

Precautions and Disclaimer

The slides contain human source material. Refer to Material Safety Data Sheet for hazard information. Lot-specific HBsAg and HIV test results are available on the product labels and/or Certificate of Analysis. All human source products should be handled using the same safety precautions employed when handling any potentially infectious material.

LCM/MAM 7/02

Table 1.

Section location for Tissue Section Slide, Multiple Organs from human adult normal neural tissues

Medulla oblongata	Corpus collosum	Spinal cord	Peripheral nerve	Slide Identification
Cerebral cortex	Mesencephalon	Cerebellum	Pone	

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.