

## Product Information

# Anti-Transportin 1 Antibody, Mouse Monoclonal

Clone D45, purified from hybridoma cell culture

**T0825**

## Product Description

Monoclonal Anti-Transportin 1 (mouse IgG1 isotype) is derived from the hybridoma D45 produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with recombinant human transportin 1.<sup>1</sup> The isotype is determined using ImmunoType™ Kit (T0825 ISO-1) and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents (T0825 ISO-2).

Monoclonal Anti-Transportin1 recognizes human,<sup>1</sup> bovine, canine, rat, and mouse transportin 1 (approx. 90 kDa). The epitope recognized by the antibody resides in the second quarter of transportin 1.<sup>1</sup> The antibody can be used in various immunochemical techniques including immunoblotting,<sup>1</sup> immunoprecipitation,<sup>1</sup> and immunocytochemistry.<sup>1</sup>

Nucleocytoplasmic transport, which consists of entry and exit of molecules to and from the nucleus, is crucial for cell function. Both the entry and exit processes play an important role in the regulation of diverse cellular processes including growth factor-mediated signaling, stress responses, cell cycle control, and gene transcription and translation.<sup>2</sup> Eukaryotic cells are equipped with machinery charged with the responsibility of transporting a vast number of molecules in and out of the nucleus in a rapid, accurate, and often regulated manner. The cargoes for this machinery are diverse, comprising proteins as well as RNA-protein complexes (RNPs, ribonucleoproteins). Proteins and RNAs are imported and exported through nuclear pore complexes (NPCs), which perforate the double bilayer of the nuclear envelope.<sup>3</sup> The nucleocytoplasmic transport is mediated by specific soluble receptors. The import receptors are called importins and the export receptors known as exportins. The family of import receptors consists of the importin  $\beta$  family that bind to the importin  $\alpha$  adaptor protein. Transportin belongs to the importin  $\beta$  family of receptors and is responsible for importing hnRNP A1 back to the nucleus after it exports mRNA. Transportin recognizes a 38 amino acid motif called M9 domain that is enriched with aromatic residues and glycine. This motif is found in hnRNP A1, which is an abundant mRNA binding protein.<sup>1, 5, 6</sup>

Transportin 1 is localized in the cytoplasm, nucleoplasm, and nuclear rim, similar to the localization of importin  $\beta$ . This suggests that transportin 1 may interact with NPC during translocation. Monoclonal antibodies specific for Transportin 1 are an important tool for studying import and export of proteins to and from the nucleus.

## Reagent

Monoclonal Anti-Transportin 1 is supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 1% bovine serum albumin and 15 mM sodium azide.

Antibody Concentration: Approx. 0.5 mg/mL.

## Precautions and Disclaimer

Due to the sodium azide content, a safety data sheet (SDS) for this product has been sent to the attention of the safety officer of your institution. Consult the SDS for information regarding hazards and safe handling practices.

## Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in frost-free freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilutions should be discarded if not used within 12 hours.

## Product Profile

By immunoblotting, a working antibody concentration of 0.1-0.2 µg/mL is recommended using HeLa nuclear extract.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining the optimal working concentration by titration.

## References

1. Siomi, M.C., et al., J. Cell Biol., **138**, 1181-1192 (1997).
2. Mattaj, I.W., and Englmeier, L., Annu. Rev. Biochem., **67**, 265-306 (1998).
3. Ryan, K.J., and Wentz, S., Curr. Opin. Cell Biol., **12**, 361-371 (2000).
4. Weis, K., et al., Cell, **112**, 441-451 (2003).
5. Siomi, M.C., et al., Mol. Cell. Biol., **18**, 4141-4148 (1998).
6. Pollard, V.W., et al., Cell, **86**, 985-994 (1996).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

T0825 dat Rev 04/22

**MILLIPORE**  
**SIGMA**