

SIGMA QUALITY CONTROL TEST PROCEDURE

Product Information**Enzymatic Assay of ASPARAGINASE
(EC 3.5.1.1)****PRINCIPLE:**

L-Asparagine + H₂O $\xrightarrow{\text{Asparaginase}}$ L-Aspartate + NH₃

CONDITIONS: T = 37°C, pH = 8.6, A_{436nm}, Light path = 1 cm

METHOD: Spectrophotometric Stop Rate Determination

REAGENTS:

- A. 50 mM Tris Buffer, pH 8.6 at 37°C
(Prepare 100 ml in deionized water using Trizma Base, Prod. No. T-1503. Adjust to pH 8.6 at 37°C with 1 M HCl.)
- B. 189 mM L-Asparagine Solution
(Prepare 10 ml in deionized water using L-Asparagine, Anhydrous, Prod. No. A-0884.)
- C. 6 mM Ammonium Sulfate Standard Solution ((NH₄)₂SO₄ Std)
(Prepare 100 ml deionized water using Ammonium Sulfate, Grade I, Prod. No. A-5132.)
- D. 1.5 M Trichloroacetic Acid (TCA)
(Prepare 10 ml in deionized water using Trichloroacetic Acid, 6.1 N Solution, Stock No. 490-10.)
- E. Ammonia Color Reagent
(Use Nessler's Reagent, Aldrich Stock No. 34,518-8.)
- F. Asparaginase Enzyme Solution
(Immediately before use, prepare a solution containing 2.0 - 4.0 units/ml of Asparaginase in cold deionized water.)

Enzymatic Assay of ASPARAGINASE (EC 3.5.1.1)

PROCEDURES:

Step 1:

Pipette (in milliliters) the following reagents into suitable tubes:

	<u>Test</u>	<u>Blank</u>	<u>Std 1</u>	<u>Std 2</u>	<u>Std 3</u>	<u>Std Blank</u>
Reagent A (Buffer)	1.00	1.00	1.00	1.00	1.00	1.00
Reagent B (L-Asparagine Soln)	0.10	0.10	----	----	----	----
Reagent C ((NH ₄) ₂ SO ₄ Std)	----	----	0.25	0.50	1.00	----
Deionized Water	0.90	0.90	0.85	0.60	0.10	1.10

Equilibrate to 37°C. Then add:

Reagent F (Enzyme Solution)	0.10	----	----	----	----	----
--------------------------------	------	------	------	------	------	------

Immediately mix by inversion and incubate at 37°C for 30 minutes. Then add:

Reagent D (TCA)	0.10	0.10	0.10	0.10	0.10	0.10
Reagent F (Enzyme Solution)	----	0.10	----	----	----	----

Mix by inversion. Centrifuge for 2 minutes to clarify.

Step 2:

Pipette (in milliliters) the following reagents into suitable containers:

Deionized Water	4.30	4.30	4.30	4.30	4.30	4.30
Supernatant (from Step 1)	0.20	0.20	0.20	0.20	0.20	0.20
Reagent E (Ammonia Color Reagent)	0.50	0.50	0.50	0.50	0.50	0.50

Immediately mix by inversion and after 1 minute record the A_{436nm} for Standards, Tests, and Blanks.

Enzymatic Assay of ASPARAGINASE (EC 3.5.1.1)

CALCULATIONS:

Standard Curve:

$$\Delta A_{436\text{nm}} \text{ Standard} = A_{436\text{nm}} \text{ Standard} - A_{436\text{nm}} \text{ Standard Blank}$$

Prepare a standard curve by plotting the $\Delta A_{436\text{nm}}$ of the Standard versus Ammonia (NH_3) concentration. Note that 1 mole of Ammonium Sulfate corresponds to 2 moles of Ammonia, therefore a 6 mM Ammonium Sulfate standard is equivalent to a 12 mM ammonium standard.

Sample Determination:

$$\Delta A_{436\text{nm}} \text{ Test} = A_{436\text{nm}} \text{ Test} - A_{436\text{nm}} \text{ Test Blank}$$

Determine the μmoles of NH_3 liberated using the standard curve.

$$\text{Units/ml enzyme} = \frac{(\mu\text{mole of NH}_3 \text{ liberated})(2.20)}{(0.2)(30)(0.1)}$$

2.20 = Volume of Step 1

0.2 = Volume of Step 1 used in Step 2

30 = Time of assay in minutes

0.1 = Volume (in milliliters) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will liberate 1.0 μmole of ammonia from L-asparagine per minute at pH 8.6 at 37°C.

FINAL ASSAY CONCENTRATION:

In a 2.20 ml reaction mix, the final concentrations are, 23 mM Tris, 8.6 mM L-asparagine and 0.2 - 0.4 units of asparaginase.

**Enzymatic Assay of ASPARAGINASE
(EC 3.5.1.1)**

REFERENCES:

Shirfrin, S., Parrott, C.L. and Luborsky, S.W. (1974) *Journal of Biological Chemistry* **249**, 1335-1340

NOTES:

1. This assay is based on the cited reference.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packaging slip.