



QUALITY CONTROL RANGES

MILLIPLEX® MAP Human Kidney Injury Magnetic Bead Panel 1 Catalog # HKI1MAG-99K

Control Catalog # HKI1-6099-1 Lot # HKI1-108 and HKI1-208

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Calbindin	Control 1	1.5 – 5.0	ng/mL	Renin	Control 1	0.16 – 0.33	ng/mL
	Control 2	14 – 30	ng/mL		Control 2	1.7 – 3.5	ng/mL
GSTα	Control 1	0.27 – 0.57	ng/mL	FABP-1	Control 1	50 – 105	ng/mL
	Control 2	2.8 – 5.7	ng/mL		Control 2	412 – 855	ng/mL
TIMP-1	Control 1	0.63 – 1.3	ng/mL	Collagen IV	Control 1	2.4 – 5.0	ng/mL
	Control 2	6.8 – 14	ng/mL		Control 2	24 – 50	ng/mL
KIM-1	Control 1	0.34 – 0.70	ng/mL	TFF-3	Control 1	3.3 – 7.1	ng/mL
	Control 2	3.4 – 7.0	ng/mL		Control 2	35 – 73	ng/mL
IP-10	Control 1	0.032 – 0.066	ng/mL				
	Control 2	0.35 – 0.72	ng/mL				

CAT# HKI1-6099-1, Lot# HKI1-108.208 (10 January 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA