

M6930 Mycoplasma enrichment supplement

A selective supplement used for the isolation of *Mycoplasma*.

Composition

(per vial sufficient for 100 ml medium)

Horse Serum 20.00 ml

Yeast Extract (25% w/v solution) 10.00 ml

Thallous Acetate 25.00 mg

Penicillin G 20,000 IU

Store at -20°C.

Directions

Aseptically add the contents of one vial to 70 ml of sterile, molten Mycoplasma Agar Base (Cat. No M2178). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

