



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-521-8956 • (314) 771-5765
Fax 800-325-5052 • (314) 771-5757
sigma-aldrich.com

Product Information

Tissue Section Slide, Single Organ, from rat brain

Product Number **R 9651**
Store at Room Temperature

Product Description

This slide contains a single paraffin embedded section prepared from rat brain. The tissue section is fixed in paraformaldehyde after excision and then a section approximately 5 μ m thick is embedded in paraffin. The tissue slide is suitable for most immunohistochemistry applications and as an immunohistochemistry control slide. The slide may also be used for *in situ* hybridization and *in situ* PCR studies, including tissue specific protein and mRNA expression comparing the location of novel genes among different tissues.

Component

Slides containing tissues of single organs or single systems are available in sets of five slides.

Storage/Stability

It is recommended to store the slide at room temperature.

Precautions and Disclaimer

Consult the MSDS for information regarding hazards and safe handling practices.

LCM/MAM 7/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.