



QUALITY CONTROL RANGES

HUMAN NEUROPEPTIDE Y (NPY) ELISA CATALOG # EZHNPY-25K

STANDARD, LYOPHILIZED CATALOG # E8025-K

STANDARD	HUMAN NPY VALUE
LOT NUMBER: E825-14K	1011 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2.0 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6025-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL CATALOG #E6025-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HNPY 114	25 – 51	pg/mL
QC 2	HNPY 214	198 – 411	pg/mL

E6025-K: HNPY-114.214 (20 September 2024)

EMD MILIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA