



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Anti-Sodium Channel ASIC3

(DRASIC)

Developed in Rabbit, Affinity Isolated Antibody

Product Number **S 5070**

Product Description

Anti-Sodium Channel ASIC3 is developed in rabbit using a synthetic peptide KPRSGLEEAQRRQAS(C) corresponding to residues 2-16 of rat ASIC3 as the immunogen. The sequence has 14/15 residues identical in mouse and 11/15 identical in human. The antibody was affinity isolated on immobilized immunogen.

Anti-Sodium Channel ASIC3 specifically recognizes ASIC3 in rat dorsal root ganglion lysate by immunoblotting.

The non-voltage gated sodium channel superfamily encodes a group of conserved proteins involved in diverse biological processes, including sodium homeostasis, salt taste, nociception, touch sensation, and response to acidic pH.¹ Members of this family include the Acid-Sensing Ion Channels (ASICs), the Epithelial Sodium Channel (ENaC) and the Degenerins; most of these channels are blocked by the diuretic amiloride.

Acid-sensing ion channels (ASICs) are non-voltage-gated Na⁺ channels that are transiently activated by a rapid drop in extracellular pH.² To date four ASIC genes have been identified in mammals, two of which exist in different splice variants.³ Functional ASICs are formed by homo- or heterotetrameric assembly of ASIC subunits 1a, 1b, 2a, 2b and 3. Each subunit has two transmembrane domains separated by a large extracellular loop; four subunits form the functional channel.³

ASIC1a, 2a and 2b are expressed in the peripheral and the central nervous system, whereas expression of ASIC1b and 3 is restricted to the peripheral nervous system.² In the periphery, they are found in sensory neurons involved in pain perception, suggesting a role in pain perception following tissue acidosis caused by tissue injury, ischemia, or inflammation.⁴⁻⁶ ASICs in central neurons might play a similar role and might also contribute to the neuronal death associated with brain ischemia or epilepsy, both of which are accompanied

by extracellular acidification. In addition, studies with knockout mice indicate a central role for ASICs in learning and memory.^{3,7} Some of the evidence is contradictory, and the exact functional significance of the ASICs in pain perception, mechanoreception, plasticity, learning and memory have yet to be determined.

Reagent

The antibody is supplied lyophilized from phosphate buffered saline (PBS), pH 7.4, 1% BSA, and 0.05% sodium azide as preservative.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling.

Preparation Instructions

Reconstitute the lyophilized vial with 0.05 ml or 0.2 ml deionized water, depending on the package size purchased. Antibody dilutions should be made in buffer containing 1% bovine serum albumin.

Storage/Stability

Lyophilized powder can be stored intact at room temperature for several weeks. For extended storage, it should be stored at -20 °C or below. The reconstituted solution can be stored at 4 °C for up to 2 weeks. For longer storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Centrifuge all antibody preparations before use (10000 x g 5 min). Working dilution samples should be discarded if not used within 12 hours.

Product Profile

The recommended working dilution is 1:200 for immunoblotting.

Note: In order to obtain best results in different techniques and preparations we recommend determining optimal working concentration by titration test.

References

1. Kellenberger, S. and Schild, L., Epithelial sodium channel/degenerin family of ion channels: a variety of functions for a shared structure., *Physiol. Rev.*, **82**, 735–767 (2002).
2. Poirrot, Olivier, et al., Selective regulation of acid-sensing ion channel 1 by serine proteases, *J. Biol. Chem.*, **279**, 38448–38457 (2004).
3. Krishtal, O., The ASICs: Signaling molecules? Modulators?, *Trends in Neurosciences*, **26**, 477–483 (2003).
4. Deval, E., et al., ASIC2b-dependent regulation of ASIC3, an essential acid-sensing ion channel subunit in sensory neurons via the partner protein PICK-1., *J. Biol. Chem.*, **279**, 19531–19539 (2004).
5. Mamet, J., et al., Proinflammatory mediators, stimulators of sensory neuron excitability via the expression of acid-sensing ion channels., *J. Neurosci.*, **22**, 10662–10670 (2002).
6. Hesselager, M., et al., pH dependency and desensitization kinetics of heterologously expressed combinations of acid-sensing ion channel subunits., *J. Biol. Chem.*, **279**, 11006–11015 (2004).
7. Wemmie, J.A., et al., Overexpression of acid-sensing ion channel 1a in transgenic mice increases acquired fear-related behavior, *Proc. Natl. Acad. Sci. USA*, **101**, 3621–3626 (2004).

MCT,PHC 02/05-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.