

**RABBIT ANTI-PRESENILIN-1
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB5308

LOT NUMBER:

QUANTITY: 100 µL

SPECIFICITY: Presenilin-1, Loop.

IMMUNOGEN: Human Presenilin-1 (amino acids 275-367) fused to C-terminus of maltose binding protein.

APPLICATIONS: Western blot: 1:500-1:3,000.
Immunocytochemistry: Cos7, SY5Y cell lines.
Immunoprecipitation.
Optimal working dilutions must be determined by end user.

SPECIES REACTIVITIES: Human, monkey, mouse and hamster.

FORMAT: Rabbit serum.

PRESENTATION: Liquid in PBS with 0.02% sodium azide.

STORAGE/HANDLING: Maintain frozen at -20°C in undiluted aliquots for up to 12 months after date of receipt.
Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*