Sigma-Aldrich<sub>®</sub>

**Data Sheet** 

## Oli-neu Mouse Oligodendrocyte Precursor Cell Line

Immortalized Cell Line Cat. # SCC261

FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURES NOT FOR HUMAN OR ANIMAL CONSUMPTION pack size: <u>></u>1x10<sup>6</sup> viable cells/vial

Store at: liquid nitrogen

## Background

Oligodendrocytes, comprising 5-10% of the total population of glial cells in the central nervous system, are responsible for producing the myelin sheath encasing axons which is essential for efficient electric conduction of neural impulses.<sup>1</sup> Oligodendrocytes arise from precursor cells that differentiate just before birth. Disorders of myelination of axons can lead to emergence of neurodegenerative diseases such as Alzheimer's and multiple sclerosis, underscoring the critical importance of understanding oligodendrocyte biology and interactions. The myelination process occurs within the complex milieu of neuronal cell types, thus is difficult to reconstruct without cellular models.

Oli-neu is an immortalized mouse oligodendroglial cell line that recapitulates morphology and in vivo physiological activity of oligodendroglial cells. Oli-neu cells differentiate into oligodendrocytes upon daily addition of 1 mM dibutyrl-cAMP.<sup>2</sup> Oli-neu cells have been utilized in a broad range of studies of oligodendrocyte biology and has greatly contributed to the understanding of the mechanisms of myelination and molecular players in myelination disorders. Differentiated Oli-neu oligodendrocyte express O4 antigen, an oligodendrocyte marker, and manifest cell processes characteristic of oligodendrocyte morphology.<sup>2</sup> The Oli-neu cell line is a versatile cellular model for oligodendrocyte biology and associated neurophysiological investigations.

### Source

Oli-neu was derived from oligodendrocytes and oligodendrocyte precursor cells isolated from the cerebellum of a 6day-old NMRI mouse and infected with retrovirus containing the *t-neu* tyrosine kinase oncogene.<sup>2</sup>

## Short Tandem Repeat (STR Profile)

M18-3: 17	M15-3: 20.3
M4-2: 19.3	M6-4: 16, 17, 18
M6-7: 11, 12	M11-2: 15
M19-2: 11, 13	M17-2: 15, 16
M1-2: 14	M12-1: 19, 20
M7-1: 25.2, 29	M5-5: 12
M1-1: 14, 15	MX-1: 25
M3-2: 10, 13	M13-1: 15
M8-1: 16	D8S1106: NA

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



# Submit your published journal article and earn credit toward future EMD Millipore purchases. Visit www.emdmillipore.com/publicationrewards to learn more! EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited. Chemicon®, Upstate®, and all other trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2021 Merck KGaA, Darmstadt, Germany. All rights reserved. 20306518; Version 2.0

M2-1: 9 D4S2408: NA Quality Control Testing

- Each vial contains  $\geq 1X10^6$  viable cells.
- Cells are tested negative for infectious diseases by a Mouse Essential CLEAR Panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of mouse origin and negative for inter-species contamination from rat, human, chinese
  hamster, Golden Syrian hamster, and Non-human Primate (NHP) as assessed by a Contamination Clear panel by
  Charles River Animal Diagnostic Services
- Cells are negative for mycoplasma contamination.

## Storage and Handling

Oli-neu mouse oligodendrocyte precursor cells should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting functionality.

## **Representative Data**



**Figure 1.** Bright-field images of oli-neu cells undifferentiated (A) and differentiated (B) for 4 days with 1 mM dibutyrl cAMP. Actin staining of undifferentiated cells (C). Cell confluency (66%) was assessed using MiliCell Digital Cell Imager (D, Sigma MDCI 10000). Differentiated cells express O4 oligodendrocyte marker (E, Sigma 07139).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



## Protocols

Thawing Cells

- Prepare poly-L-lysine coated flasks: Coat flasks overnight at room temperature with poly-l-lysine 0.01% solution (Sigma A-005-M). When cells are ready to be plated, aspirate the coating solution and rinse twice with sterile water before plating Oli-neu cells.
- Prepare Oli-neu Expansion Medium comprised of DMEM High Glucose Medium (Sigma D5796) containing 10 μg/mL Insulin (Sigma I9278-5ML), 10 μg/mL Apo-transferrin (Sigma T4382), 1X B27<sup>™</sup> Supplement (50X) (ThermoFisher 17504044), 100 μM Putrescine dihydrochloride (Sigma P5780), 200 nM Progesterone (Sigma P0130), 500 nM 3,3',5-Triiodo-L-thyronine (Sigma T6397), 520 nM L-Thyroxin (Sigma T1775), and 220 nM Sodium Selenite (Sigma S5261). Note: Expansion medium should not contain FCS as this may push the cells towards the astroglial lineage.
- 3. Remove the vial of frozen Oli-neu cells from liquid nitrogen and incubate in a 37°C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells. **IMPORTANT: Do not vortex the cells.**
- 4. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
- 5. In a laminar flow hood, use a 1 or 2 mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
- 6. Using a 10 mL pipette, slowly add dropwise 9 mL of cold Oli-neu Expansion Medium (Step 2 above) to the 15 mL conical tube. **IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.**
- 7. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles. **IMPORTANT**: Do not vortex the cells.
- 8. Centrifuge the tube at 2-8°C at 300 x g for 2-3 minutes to pellet the cells. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
- Resuspend the cells in 10 15 mL of Oli-neu Expansion Medium. Transfer the cell mixture to a poly-I-lysine coated T75 tissue culture flask. Incubate the cells at 37°C in a humidified incubator with 5% CO<sub>2</sub>. Exchange medium every other day.
- 10. Oli-neu should be passaged at ~80-85% confluency to a new poly-I-lysine coated flask. Typical split ratio is 1:6. Oli-neu cells usually do not grow as dense as for example HEK cells. Thus at full confluency, there will be some space between the cells. Monitor the cells carefully as they may change growing behavior from time to time. Adjust the plating density accordingly (i.e. slowly growing cells have to be seeded at a higher density and vice versa).

### Oli-neu Differentiation

Differentiation can be induced by daily addition of 1 mM dibutyrl cAMP. Maximal differentiation is estimated after one week of treatment, but normally 3-4 days are sufficient to achieve considerable differentiation. A stop in treatment or medium change may lead to de-differentiation.

- 1. Prepare a subconfluent 10-cm plate by plating proximately 4-5 million cells per plate.
- 2. Initiate differentiation by adding dibutyrl cAMP to 1 mM final concentration.

#### Cryopreservation of Cells

Oli-neu cells may be frozen in RPMI-1640 medium containing 20% FBS and 10% DMSO using a Nalgene slow freeze Mr. Frosty container.

### References

1. Curr Pharm Des. 2016, 22(6): 656-679.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



2. Eur J Neurosci. 1995, 7(6): 1245-1265.

#### ACADEMIC USE AGREEMENT (subject to local law)

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA ("SELLER") CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means Oli-neu Mouse Oligodendrocyte Precursor Cell Line (SCC261).

"Academic Research" means any internal *in vitro* research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- · Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means *any* product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

#### COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

### GMO

- This product contains genetically modified organisms.
- Este producto contiene organismos genéticamente modificados.
- Questo prodotto contiene degli organismi geneticamente modificati.

Dieses Produkt enthält genetisch modifizierte Organismen.

Ce produit contient des organismes génétiquement modifiés.

Dit product bevat genetisch gewijzigde organismen.

Tämä tuote sisältää geneettisesti muutettuja organismeja. Denna produkt innehåller genetiskt ändrade organismer.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.