

SIGMA QUALITY CONTROL TEST PROCEDURE

ProductInformation

Enzymatic Assay of GLUTATHIONE REDUCTASE (EC 1.6.4.2)

PRINCIPLE:

 β -NADPH + GSSG ^{Glutathione Reductase} > β -NADP + 2 GSH

Abbreviations used: β -NADPH = β -Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form β -NADP = β -Nicotinamide Adenine Dinucleotide Phosphate, Oxidized Form GSSG = Glutathione, Oxidized Form GSH = Glutathione, Reduced Form

CONDITIONS: $T = 25^{\circ}C$, pH = 7.6, A_{340nm} , Light Path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 100 mM Potassium Phosphate Buffer with 3.4 mM Ethylenediaminetetraacetic Acid (EDTA), pH 7.6 at 25°C
 (Prepare 200 mI in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379 and Ethylenediaminetetraacetic Acid, Dipotassium Salt, Sigma Stock No. ED2P. Adjust to pH 7.6 at 25°C with 1 M KOH.)
- B. 30 mM Glutathione Substrate Solution (GSSG) (Prepare 5 ml in deionized water using Glutathione, Oxidized Form, Disodium Salt, Sigma Prod. No. G-4626.)
- C. 0.8 mM β-Nicotinamide Adenine Dinucleotide Phosphate, Reduced Form Solution (β-NADPH) (Prepare 5 ml in cold Reagent A using β-Nicotinamide Adenine Dinucleotide Phosphate, Tetrasodium Salt, Sigma Prod. No. N-1630.)
- D. 1.0% (w/v) Bovine Serum Albumin (BSA) (Prepare 100 ml in Reagent A using Albumin, Bovine, Sigma Prod. No. A-4503. This solution should be kept cold.)

Enzymatic Assay of GLUTATHIONE REDUCTASE (EC 1.6.4.2)

REAGENTS: (continued)

E. Glutathione Reductase Enzyme Solution (Immediately before use, prepare a solution containing 0.30 - 0.60 unit/ml of Glutathione Reductase in cold Reagent D.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	Test	Blank
Deionized Water	0.65	0.65
Reagent A (Buffer)	1.50	1.50
Reagent B (GSSG)	0.10	0.10
Reagent C (β-NADPH)	0.35	0.35
Reagent D (BSA)	0.30	0.40

Mix by inversion and equilibrate to 25°C. Monitor the A_{340nm} until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent E (Enzyme Solution)	0.10	
Reagent E (Enzyme Colution)	0.10	

Immediately mix by inversion and record the decrease in the A_{340nm} for approximately 5 minutes. Obtain the ΔA_{340nm} /minute using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

 $(\Delta A_{340nm}/min \text{ Test} - \Delta A_{340nm}/min \text{ Blank})(3)(df)$

Units/ml enzyme =

(6.22) (0.1)

3 = Total volume (in milliliters) of assay df = Dilution factor 6.22 = Millimolar extinction coefficient of β -NADPH at 340 nm 0.1 = Volume (in milliliters) of enzyme used

units/ml enzyme

Units/mg solid = mg solid/ml enzyme

Enzymatic Assay of GLUTATHIONE REDUCTASE (EC 1.6.4.2)

CALCULATIONS: (continued)

units/ml enzyme

Units/mg protein = mg protein/ml enzyme

UNIT DEFINITION:

One unit will reduce 1.0 µmole of oxidized glutathione per minute at pH 7.6 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 75 mM potassium phosphate, 2.6 mM ethylenediaminetetraacetic acid, 1 mM glutathione, 0.09 mM β -nicotinamide adenine dinucleotide phosphate, reduced form, 0.13% (w/v) bovine serum albumin, and 0.03 - 0.06 unit of glutathione reductase.

REFERENCE:

Mavis, R.D. and Stellwagen, E. (1968) Journal of Biological Chemistry 243, 809-814

NOTES:

- 1. This assay is based on the cited reference.
- 2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.