



QUALITY CONTROL RANGES
MILLIPLEX[®] MAP Human Apolipoprotein Magnetic Bead Panel Kit
Catalog # AP0MAG-62K

Standard Catalog # HAP0-8062
Lot# 862M-11K

Apolipoproteins	Standard Value		Apolipoproteins	Standard Value
Apo AI	54033 ng/mL		Apo AII	4518 ng/mL
Apo CIII	3917 ng/mL		Apo B	34639 ng/mL
Apo E	2199 ng/mL		Apo CII	1913 ng/mL

Control Catalog # HAP0-6062
Lot# AP0M-111 and AP0M-211

<u>Apolipoproteins</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Apo-AI	Control I	AP0M-111	229 - 475	ng/mL
	Control II	AP0M-211	1115 - 2317	ng/mL
Apo-CIII	Control I	AP0M-111	16 - 33	ng/mL
	Control II	AP0M-211	82 - 170	ng/mL
Apo-E	Control I	AP0M-111	6.1 - 12.6	ng/mL
	Control II	AP0M-211	28 - 58	ng/mL
Apo-AII	Control I	AP0M-111	9.7 - 20	ng/mL
	Control II	AP0M-211	68 - 140	ng/mL
Apo-B	Control I	AP0M-111	96 - 200	ng/mL
	Control II	AP0M-211	518 - 1076	ng/mL
Apo-CII	Control I	AP0M-111	6.2 - 12.9	ng/mL
	Control II	AP0M-211	33 - 68	ng/mL

HAP0-6062: AP0M-111.211 (29 November 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA