



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

Anti-Avidin
Developed in Rabbit
Delipidized Whole Antiserum

Product No. **A5170**

Product Description

Antiserum is developed in rabbit using purified egg white avidin as the immunogen. The antiserum has been treated to remove lipoproteins. Rabbit anti-Avidin is supplied as a liquid containing 15 mM sodium azide as preservative.

Antiserum is determined to be immunospecific for egg white avidin by immunoelectrophoresis versus egg white - avidin. No reactivity is observed versus streptavidin.

Precautions

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult

the MSDS for information regarding hazards and safe handling practices.

Product Profile

Protein concentration is determined by Biuret.

Storage

For continuous use, store at 2-8 °C for up to one month.

For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Kaa/Pcs 9/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.