

don't use printout for evaluation!!!



1.14410
MQuant®
Silicate Test
Silicat-Test
Test Silicates
Test Silicatos


Si

Color card
Farbkarte
Carte colorimétrique
Tarjeta colorimétrica




© 2025

1.



2.



7.69144.1008-XXXXXX 09/2025

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.sigmaldrich.com/colorimetry

EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321

MilliporeSigma Canada Ltd., 2149 Winston Park Dr,
Oakville, Ontario, L6H 6J8, Canada,
Phone: +1 800-565-1400

0 0.01 0.02 0.04 0.06 0.08 0.10 0.15 0.20 0.25 mg/l Si



20 ml



3 x Si-1



pH 1.2 - 1.6



3 min




3 x Si-2



10 x Si-3



2 min



0 0.021 0.043 0.086 0.13 0.17 0.21 0.32 0.43 0.53 mg/l SiO₂

